

#### **MODULAR PERSONAL CARE SOLUTIONS**

SAFE | INDEPENDENT | ERGONOMIC

Distributed by:





# Customer insight-led design to enhance independence, safety and flexibility.

## Personal care is highly private activity and is an important part of a person's life and wellbeing.

The Aquatec Ocean Ergo Family of shower chair commodes have been designed with user comfort, independence and dignity in mind. Based on the well established Aquatec Ocean range, the Ocean Ergo family comes packed with a range of new features, designed to make personal care activities easier for both the user and carer.

#### Safety, independence and dignity

A newly contoured seat surface has been developed to encourage a more upright and stable seating position, promoting user independence and overall comfort. Enhanced stability at the front of the base allows for a smaller footprint and offers users added safety and reassurance when seated in the chair. Furthermore, the optimised wheelbase makes it even easier for a carer to manoeuvre the chair in small spaces as well as allowing for a better fit over standard toilets for independent toileting.

#### Modular and flexible design

The modular design of the Ocean Ergo family means that the chairs can be adapted to suit a wide range of user needs by making some simple adjustments, or by adding some of our easy to fit, off-the-shelf-accessories.



### **Tilt in space models**





## **GREATER ADAPTABILITY**

With an increased height adjustment range on the footrests and armrests.

#### **COMPACT AND EASY TO HANDLE**


Optimised wheel base makes it easy to manoeuvre in small spaces.




#### **INNOVATIVE ERGONOMIC SEAT**

Insight-led design encourages a more upright seating position for enhanced independence and comfort.

#### **OTHER** features include:







#### **EXTREMELY STABLE**

Enhanced stability at the front for additional safety and reassurance.

## **Highly configurable**



Positioning



Safety



toileting and seating accessories that can be used to support a wide range of conditions and postural needs.

Available with an extensive range of positioning, safety,

See pages 5 and 6 for full range of accessories.



Toileting



**Seating and Backrests** 




# Introducing the Ocean Ergo family

From a standard transit model, through to tilt and recline, meet the flexible range of shower chairs that offer personal care solutions for a wide range of physical conditions.



### **Tech Specs**



## **Accessories**

## **Seating and backrests**



Standard soft seat BTC087674

With key shaped cut out.



Ergonomic soft seat **BTC087642** 

Extra soft, textured surface for pelvic stability.



Solid plastic backrest **BTC087688** 

Offers a stable upright position and is easy to clean.



Soft seat insert BTC087676

Can be used with the soft seat when commode function is not required.



Variable soft seat BTC087680

Hygiene cut out can be positioned in all 4 directions.



Wider XL armrests BTC087628

Increases the distance between the armrests by 40mm on each side.



Small soft seat BTC087678

With small, oval opening for smaller and young adults.



Universal soft seat BTC087684

Comfortable, closed seat surface.



Soft backrest cushion **BTC087702** 

For additional stability and comfort. Will also benefit smaller users and young adults.



Small soft seat insert

Can be used with the small soft seat when the commode function is not required.



XL back BTC087690

Increases the width between the armrests by 80mm. Can be retrofitted to standard Ocean Ergo models.



Comfort backrest BTC087694 - Ergo BTC087692 - Ergo XL

For additional comfort and easy cleaning.

## **Toileting**



### Splash guard NEW DESIGN BTC087686

Clips on easily for use with the standard plastic seat, standard soft seat and small soft seat.



#### Sanitary pan w/lid - BTC087670 Toilet pan holder - BTC087672

- Ergonomic handle allows carer to remove pan from holder easily
- Scale inside the pan
- · Lid holder when lid is not required



## **Accessories**

### **Safety**



Safety bar - BTC087606 XL Safety bar - BTC087608

For additional safety and reassurance.



Hip belt - BTC087632 Chest belt - BTC087630

For safe positioning. Individually adjustable, high load capacity and easy to assemble. Machine washable.



Armrest lock BTC087624

Locks armrest into the lower or higher position.



Front and rear anti-tippers
For additional safety.

## **Positioning**



### Lateral supports BTC087700

For safe positioning. Width and height adjustable – also for small users and young adults.



Headrest holder

Safe and comfortable with various adjustment positions (pad not included).



Ergonomic headrest/neckrest

headrest/neckrest
Offers exceptional
management of the head
position. Neckrest shown.



Calf support with foot support BTC087704 - Left

BTC08//04 - Left BTC08/706 - Right

Padded, both height and angle adjustable with additional foot support.



Calf and amputee support Left - BTC087696 Right - BTC087698

Padded, both height and angle adjustable. For use with the standard plastic seat, standard soft seat and small soft seat.



# **Compatibility**

Most accessories are cross compatible with all Ocean Ergo models.	Product Code	Ocean Ergo	Ocean Ergo XL	Ocean VIP Ergo	Ocean Dual VIP Ergo
Standard soft seat	BTC087674	1	1	1	1
Soft seat insert	BTC087676	1	1	1	1
Small soft seat	BTC087678	1	1	1	1
Small soft seat insert	BTC087682	1	1	1	1
Splash guard	BTC087686	1	1	1	1
Universal soft seat	BTC087684	1	1	1	1
Ergonomic special soft seat	BTC087642	1	1	1	1
Variable soft seat	BTC087680	1	1	1	1
Safety bar	BTC087606	1	X	1	1
Safety bar XL	BTC087608	X	1	X	X
Chest belt	BTC087630	1	1	1	1
Pelvic belt	BTC087632	1	1	1	1
Armrest locking device	BTC087624	1	1	1	1
XL Back	BTC087690	1	_	X	X
Fixed backrest	BTC087688	1	X	1	1
Soft backrest cushion	BTC087702	1	1	1	1
24" wheels (set)	BTC083716	1	Х	Х	X
Lateral	BTC087700	1	1	1	1
Toilet pan	BTC087670	1	1	1	1
Toilet pan holder	BTC087672	1	1	1	1
Comfort soft backrest Ocean Ergo models	BTC087694	1	X	X	X
Comfort soft backrest Ocean Ergo XL models	BTC087692	X	1	X	X

✓ Compatible 

✓ Not Compatible 

— Included



# **Adjustments**





FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com Thailand: Phone +66 0 2108 8555 or visit www.invacare.com



www.aidacare.com.au

Distributed by: