## PRISM MEDICAL UK





# **Bariatric** C380 Ceiling Hoist

Bariatric Transfers with Respect & Dignity

#### Exclusive Distribution with High Quality Service & Support

Prism Medical UK are the global leader in innovation and safety focused medical products. This trusted brand has a proven history in Australia and is the preferred choice of ceiling hoist solutions for many Aged Care and Hospital groups. Key safety features and stringent product testing ensure Prism Ceiling Hoists are the top performers in their category.

Aidacare's partnership with Prism Medical UK ensures that extremely high quality and internationally tested products are readily available in Australia and backed with excellent sales, installation and preventative maintenance support nation wide.

#### C380 Bariatric Ceiling Hoist

Bariatric clients often exceed a hoist's safe working load and this can create an unsafe working environment. In cases where a client's weight, size or medical condition requires precise control when transferring or re-positioning, the Prism C380 Bariatric Ceiling Track Hoist provides a safe and reliable solution.

The C380 is specifically designed to lift clients that exceed a standard hoist's maximum weight capacity. As the strongest Ceiling Track Hoist in its class, the C380 safely transfers bariatric clients with respect and dignity.



#### Emergency stop/emergency lower 3

- 4. Manual or powered traverse available
- 5. Smooth stop start to ensure user comfort
- Bariatric four point carry bar 6.
- 7. Carry bar fitted with Quick Release System for simple interchange of carry bars
- 8. Easy to read LCD display

1.

2

Low profile hoist unit provides a great lift height 9.

### C380 Ceiling Hoist Range

Code	Description	S.W.L
LIC350250	C380 Bariatric - Manual Traverse	380kg
LIC350260	C380 Bariatric - Powered Traverse	380kg
LIC350270	C380 Bariatric - Powered Traverse and Return to Charge	380kg
LIC350280	C380 Bariatric - XY Powered Traverse	380kg

#### **Quality & Compliance**

- All systems are tested and installed in compliance with Australian Standard AS/NZS ISO 10535-2011: Hoists for the transfer of disabled person
- Aidacare is an ISO9001 compliant company
- Aidacare uses an accredited network of trained installers, engineers and load testers.

