

Ceiling Hoist Systems




Aidacare is an experienced supplier and installer of a range of ceiling hoist systems. We understand the risks involved in the installation and use of ceiling hoists and can help during the early-stage specification and design process to reduce risk and improve user-acceptance.


Our systems are accompanied by accredited installation guidelines which allow us to work with a range of installers. Aidacare is an ISO 9001 accredited company and consistently provides high-quality outcomes for Australia's largest Aged Care and Hospital providers.

We can supply and install a combination of motor, track and accessory options to suit client specification.


Examples of the options are provided include:




H / XY / Traverse



Linear Track with Turn Table




Straight Monorail



J Track

Ceiling Hoist Considerations

- A** Structural Attachment Point
- B** Threaded Rod/Steel Pendant
- C** Ceiling Bracket
- D** Ceiling Track
- E** Traverse (Track)
- F** Charging Station
(Various charging options)
- G** X-Y/Traverse Rail Carrier
- H** Ceiling Lift Motor
- I** Spreader/Sling Bar
- J** Lift Strap
- K** Sling
- L** Gait Assembly
- M** Linear/Straight Track
- N** Modified Door Header



Installation Solutions

No two installations are alike. Aidacare has extensive experience dealing with many different installation substrates and scenarios. Aidacare utilises manufacturer-certified installation guidelines. Some of the more complex examples are shown below. This leads to a consistent high-quality outcome for our clients.

Maintenance Safety

Aidacare fully-supports its installations with ongoing maintenance programs. Aidacare also assists our clients to develop their own internal maintenance regime.

Training & User Acceptance

Our Equipment Specialists offer user training and ongoing support. We also provide WHS signage with all equipment.

Quality & Compliance

- All systems are tested and installed in compliance with Australian Standard AS/NZS ISO 10535 Hoists for the transfer of disabled persons.
- All systems are maintained in accordance with manufacturer's guidelines and in accordance with AS/NZS 3551 Management Programs for Medical Equipment.
- Aidacare is an ISO 9001 compliant company.
- Aidacare uses an accredited network of trained installers, engineers, and load testers.



For more information about our Ceiling Hoists, Engineering solutions and wide range of manual handling products scan the QR and visit our website

