

DIGITAL CHAIR SCALE BW-3136AK



The BW-3136AK Digital Chair Scale is a durable, comfortable, and cost-effective solution for accurate patient weighing. With a 300kg capacity and smooth portability, it's designed to make patient transfers easier and everyday weighing simpler. A clear backlit display and sturdy aluminium frame ensure reliability in any clinical setting.



KEY FEATURES

- | 300kg x 50grams Max Capacity
- | Extra Wide Comfortable ABS Seat
- | Wheels & Handles for Portability of Scale
- | Four Wheel Brakes for Safety
- | Hinged Armrests for Easy Patient Transfer
- | Large LCD Backlight Display
- | Powered by AC Adapter, or AAA batteries
- | Robust Aluminium Frame
- | Cost Effective Weighing Solution
- | Hinged Footrest for Comfort and Easier Patient Transfer
- | **TGA Approved: ARTG ID 520168**

SPECIFICATIONS

Model / Capacity
Display Reading
Finish

BW-3136AK - 300kg x 0.1g
KG

Operating Temperature
Operating Humidity
Storage Temperature
Weighing Surface Area
Physical Foot Print
Carry Bag Supplied
Power Supply

Robust Aluminium Frame,
ABS Chair and Display

0° to 40°C

85% R.H.

0° to 50°C

530mm x 350mm

880mm x 650mm

N/A

100-240V AC Adapter
(included), or 4 x AAA Batteries
(non-rechargeable)

Shipping Weight

26kg

Distributed by

aidacare™
HEALTHCARE EQUIPMENT

1300 133 120 | aidacare.com.au