



Redefining Infection Control

Disposable Medical Pulp Containers



Disposable medical pulp products for everyday patient care



TO PLACE AN ORDER

TEL: 1300 133 120

EMAIL: sales@aidacare.com.au WEB: www.aidacare.com.au

Disposable Medical Pulp Range

DDC Dolphin medical pulp products are ergonomically designed for ease of use with a smooth high quality finish, helping to maintain hygiene standards and reduce the risk of cross contamination and include:-

- · Good product rigidity for ease of use
- · Careful and considerate design for controlled carrying and pouring
- Suitability for use with warm fluids
- · Easy to dispose of and manufactured from fully maceratable materials

The products are designed for general purpose use in surgeries, hospitals and care facilities. We understand the importance of products that are not only cost effective, but which also give healthcare providers the confidence to perform tasks secure in the knowledge that safety and quality are key elements of the design and manufacturing processes allowing infection control strategies to be met.

DDC Dolphin medical pulp products are made from recycled and maceratable pulp material in compliance with PAS29:1999 for healthcare disposable products and Medical Devices Directives 93/42/EEC - Class 1 as amended by Directive (2007/47EC) where required. Consequently, the products specifically designed to contain human waste will do so for a minimum of 4 hours. It is recommended that these products are disposed of to avoid the risk of leakage after this period.

Bed Pan Support and Liner

Suitable for general toileting of patients and can be used in bed or commode chair. The liner is placed inside the support before use.

Code	Description	Box Qty	Pallet Qty
UPP100210	Bed Pan 2 Litre Liner	100	18
UPP100110	Bedpan Support (figure 8)	1	N/A



Male Urinal Bottles

Enables independent patient toileting. The traditional round design is ergonomically shaped for use in bed, whilst the square design stands upright if required.

Code	Description	Box Qty	Pallet Qty
UPP100250	Male Urinal Round 875 ml	100	5



Fluid Weighing Scales

Electronic scales, quick and easy to use with easy clean surface. Maximum 5000 ml (for indication purposes only).

	Description	Box Qty	Pallet Qty
F	Fluid Weighing Scales 23 x 16 x 2 cms	1	N/A

^{*} Not a stocked item POR.









Midi & Maxi Slipper Pan Supports and Liners

Recommended for general toileting purposes of patients with limited mobility in bed.

The liner is placed inside the support and then can be slipped underneath the patient without rolling or lifting.

Code	Description	Box Qty	Pallet Qty
UPP100215	Midi Slipper Pan 1.3 Litre Liner	100	17
UPP100115	Midi Slipper Pan Support	1	N/A
UPP100220	Maxi Slipper Pan 2 Litre Liner	100	16



Detergent Proof Washbowl

Holds warm water and detergent for up to 2 hours with a maximum 4 litre capacity.

Code	Description	Box Qty	Pallet Qty
UPP100242	Detergent Proof Washbowl 4 Litre	100	12



General Purpose Bowls

General purpose bowls designed for ease of use.

Code	Description	Box Qty	Pallet Qty
UPP100235	General Purpose Bowl 1 Litre	200	24
UPP100240	General Purpose Bowl 3 Litre	100	20



Kidney Dish

General use receptacle for holding dressings, wipes, surgical instruments and other medical items.

Code	Description	Box Qty	Pallet Qty
UPP100230	Kidney Dish 700 ml	300	24



Commode Bowl

Designed for general toileting purposes in a commode chair.

Code	Description	Box Qty	Pallet Qty
UPP100225	Commode Bowl 1.7 Litre	200	10



Graduated Jug

General purpose jug with ergonomic handle for enhanced grip.

Code	Description	Box Qty	Pallet Qty
-	Graduated Jug 1 Litre	100	32



Urine Sampling Cup

Designed for a multitude of functions, but primarily used for urine samples.

Code	Description	Box Qty	Pallet Qty
-	Urine Sampling Cup 300 ml	200	36



We reserve the right to modify products without prior notice



Medical pulp and the modern sluice/dirty utility room

In the modern environment with an increased focus on healthcare acquired infections (HCAI's) and the seemingly endless demands on healthcare budgets, paper pulp is considered to be a clean and highly cost effective solution.

Disposable pulp products are considered a sustainable product as they are produced from recycled materials such as phone books, newspapers, magazines - in fact any product made of paper fibres.

Using single use pulp products to dispose of human waste certainly reduces the risk of cross-contamination and re-contamination from products moved from ward to ward.

Although it differs from hospital to hospital, or care facility to care facility, the most prominent intended use of the DDC Dolphin range of pulp products, are concerning the management of various fluids and liquids, including, but not limited to body fluids and human waste. Our product line includes items that are purpose built for managing a wide number of different, but entirely typical nursing situations that occur on an everyday basis, such as holding / collecting vomit, urine or faeces. After the designated pulp product has been used, it can safely be destroyed together with the human waste material.

DDC Dolphin medical pulp is made up of 100% biodegradable and compostable raw material, and by using along side one of our Pulpmatic macerators the waste material can be disposed of effectively through the existing drainage system.

This reduces the handling of contaminated products to the minimum, for the clear benefit of all staff involved in

By switching to pulp utensils, instead of stainless steel or plastic; staff and healthcare providers can reduce costs by avoiding costly equipment / utensil purchases, eliminating washing and storage of re-usable utensils.

DDC Dolphin Pulp Macerator Solutions



All our medical pulp products are made from 100% recycled paper and are fully compostable, recyclable and biodegradable.





Human waste including the receptacle is disposed of directly into the macerator, the entire waste and bedpan is broken down and flushed into the drainage system.

Pulp Macerators

The growing popularity of single use pulp items has been driven by the need for higher infection control standards across all healthcare facilities.

One of the key advantages of using single-use pulp items and pulp macerators is that all of the human waste, including the receptacle, are disposed of directly into the machine and there is no further clinician interaction with the pulp receptacle. The entire waste and pulp is macerated and flushed into the drainage system in one simple procedure.

DDC Dolphin's Pulpmatic macerators offer greater efficiency made possible by their unique blade systems which allow for faster and quieter operation as well as breaking items down into consistently smaller particles allowing for complete waste removal with no remaining residue left in the chamber. A narrow gap between the blade and the upper chamber prevents large particles entering the lower chamber as they cannot physically pass through the gap, helping to prevent blockages within the unit or the main drainage system.

The Pulpmatic macerator cycle uses both a high pressure water spray and carefully designed cutting blades, which macerate the waste into a fine slurry, which is then flushed into the drainage system. The chamber is then rinsed thoroughly with water and DDC Dolphin's EcoCleanse+ macerator disinfectant leaving the chamber completely clean and smelling fresh, ready for the next use.



Our macerators offer completely hands free operation by means of a sensor foot cup at the base of the machine and an optical sensor to close and start the machine, there is also an auto close option.

They are available to suit all requirements ranging from large capacity machines able to cope with 4-6 pulp items, suitable for busy sluice / dirty utility rooms, and smaller, compact machines designed for low-demand areas or single bed wards and en-suite facilities.

Using pulp macerators with single use pulp products eliminates any potential contamination risk arising from reusable products. Our macerators ensure complete disposal with faster cycle times and the minimum use of water and electricity which in turn translates into reduced running costs and greater staff efficiencies.

All of the macerators we design, manufacture, and supply are WRAS certified and use the highest quality materials and design principles, anti-microbial surfaces and hands-free operation.

DDC Dolphin are the only manufacturer in the UK able to offer the full range of sluice/dirty room products including stainless steel furniture and consumables. Nearly 30 years of experience allows us to provide comprehensive and meaningful advice and support throughout planning, design and specification processes for new or existing facilities enabling you to meet modern infection control standards as efficiently as possible.







Certificate No. 1161

TO PLACE AN ORDER

TEL: 1300 133 120

EMAIL: sales@aidacare.com.au WEB: www.aidacare.com.au



Redefining Infection Control



Accreditations





















1300 133 120

www.aidacare.com.au

