



KANJO ECO BATHLIFT

The Kanjo Eco provides perfect hygiene and safety when bathing.

With it's electric lift and lowering functions, reduced mobility during bathing is no longer an issue.

An ergonomically designed seat and electric backrest allow the user to bathe in comfort all with the press of a hand control.

- · Independent lowering and raising of the user in a bathtub
- Floating hand control (lithium battery)
- Flexible side panels maintain contact with the bath wall at all times ensuring gap free transfers
- · Up to 40 degrees of backrest recline
- Suitable with tub and corner shape baths
- · Comfort seat padding included as standard



Comfort Seat Padding



Floating Hand Control



KANJO bathlift
BTS112240

(mm)	
1105	69

Overall	Overall
Height	Depth
(mm)	(mm)
1105	690 - 900

verall epth	Seat Height
mm)	(mm)

Seat
Depth
(mm)

500

Seat Width (mm)

375 - 710

Backrest Height (mm)

680



350



GSP872728 Height Adapter - 2cm 65 - 455