









AVAILABLE TRIM COLORS







ocean





power ' pink



red



sunburst orange

SPECIFICATIONS

galactic

Maximum User Weight	136 kg
Maximum Speed	10 km/h 8 km/h without anti-tippers
Range	21 km*
Base Width	24"
Base Length	40" with anti-tippers 36" without anti-tippers
Minimum Turning Radius	26.5"
Minimum Ground Clearance	3"
Recommended Battery Type	Group 34 & Group 24 Gel
Weight, including batteries	175 kg with Group 34 batteries 183 kg with Group 24 batteries
Drive Electronics	R-Net 120A

Seat-to-Floor Height (STFH)	17.5", 18.5", 19.5"
Optional Seat Elevator	12"
Tilt Options	0° - 50° Posterior 5°, 10°, 20° & 30° Anterior
Recline Options	85° - 120° (manual) 85° - 180° (power)
Legrest Angle	85° - 170°
Backrest Height	20", 23" - 28" (by 1" incr.)
Seat Depth	14" - 22" (by 1" incr.)
Seat Width	17", 19", 21", 23"
Armrest Pad Lengths	10", 13", 16", 18"

^{*}Actual driving range will vary based on driving and battery conditions.

Permobil Australia © 2015 Permobil. Ph: 1300 845 483 Email: orders.au@permobil.com www.permobil.com.au F3 Corpus Flyer AU RevA







THE F3 BASE offers the smallest footprint of Permobil's redesigned series of front-wheel drive power wheelchairs. Its size and agility let users move quickly in and out of tight spaces like busy city streets and crowded areas like malls... or even parties. The platform has a lower, repositioned center of gravity, which greatly improves stability while traversing inclines or declines. With its fully independent suspension, the F3 provides a comfortable ride and allows users to navigate the obstacles of life.

THE CORPUS® SEATING SYSTEM. The F3 Corpus comes standard with all of the ergonomic and technological advances of the Corpus seating system from its sleek and stylish design that fits the contours of the wheelchair user's body to its revolutionary tilt and recline features. Unique to the F3 Corpus are the 12" seat elevator and the A/P tilt, which allows for optional 30° of anterior tilt and 50° of posterior tilt. There's no need to comprise comfort with a smaller chair, and the F3 is proof of that.



All things Corpus power features like elevating legrests 85° - 170°

180° of power recline 50° of posterior tilt along with 5°, 10°, 20° and 30° of anterior tilt



Compact maneuverability with an option to remove anti-tippers

fully independent for a more comfortable ride



that maintains three-points of contact for increased stability



