



GoTo

VINYL

Introducing
the new

Our lightweight,
portable seat for
kids with special
needs just **got**
even easier to
use!

Promoting inclusion
and participation in
everyday life, GoTo's
new wipe-clean vinyl
cover set make it
perfect for nursery,
school, therapy and
busy homes.



NEW FEATURES

- Wipe-clean
- Hard-wearing
- Soft-touch
- Water-resistant
- Vibrant new colour

www.fireflyfriends.com.au



Headrest Choice
(Standard or Advanced)



GREAT FUNCTIONS

5-point harness

Adjustable
lateral
support

Vinyl pressure
relief cushion

Floorsitter
(can also be used with
Splashy bath seat)



GoTo Size 1
and 2 suit children
aged between 1-8
years (approx) who
would benefit from
extra postural
support.

	Seat base to top of headrest Max	Seat base to top of headrest Min	Lateral support width Max	Lateral support width Min	Seat Depth	Max Weight	Age Range (approx)
S1	540mm	380mm	260mm	190mm	157mm	15 kg	1-3
S2	680mm	510mm	300mm	210mm	244mm	30 kg	3-8

Distributed by:



1300 133 120

www.aidacare.com.au