




# FormFit Lift Recliner

## DESIGN

The Sleek form of the backrest built with layering of foam and fibre offers outstanding lateral support to help maintain an optimal seating position, and has concealed behind it an impressive range of adjustments to give the option of finetuning your comfort. It also helps accommodate adjustments in the case that your seating needs change over time.



## ACTIONS



The chair comes standard with the Oscar **Dual Motor** Platform, allowing you to recline the backrest independently of the footrest, allowing the greatest range of flexibility and comfort.

There is the option to enhance your functionality with the **Comfort Pack** that adds a power adjustable headrest functionality that is designed to follow the natural arc of the neck for those who need additional neck support and variable positioning options, This option also includes a pneumatic lumbar adjustment.

For those needing a bit of extra lift there is the option of the **Heavy Duty** two function system that has the same functionality with an increased user capacity of 180Kg.



Exclusively distributed by  
**aidacare**<sup>TM</sup>  
HEALTHCARE EQUIPMENT

1300 133 120 | [aidacare.com.au](http://aidacare.com.au)



# FormFit Lift Recliner

The FormFit is an outstanding lift recliner chair proudly made in Australia by Oscar Furniture. It is made to be scripted and highly customisable so that anyone can get the perfect Form and Fit for their unique needs.

When we are seated in a chair for a significant amount of time we all know fit is critical to comfort. That's where the FormFit has you covered.

Take your seat on Australia's finest seat foams with confidence that it will continue to feel this great for years to come.

All this comes upholstered in your choice of fabric from the vibrant Versotela range with vinyl and Leather options also available to choose from. There is a colour and finish to suit every room, style and individuals need.



## FEATURES

- 130kg weight capacity
- Controller available on either side of the chair
- Emergency battery backup system 2X 9v batteries
- 50mm castors with Large foot brake on rear
- Onsite adjustable Seat depth
- Made to measure, Scriptable sizing options
- Available in wide range of Versotela Fabrics
- 2 Year Warranty on mechanism, electrics and fabric
- 5 Year Warranty on timber frame
- Manufactured in Australia by Oscar Furniture

## SPECIFICATIONS



	Seat Height	Seat Width	Seat Depth	Arm Height	Backrest Height	Footrest Length
<b>SMALL SIZE</b> CHP231800						
<b>MIN ADJUSTMENTS</b>	-25mm	-50mm	-75mm	-60mm	-25mm	-25mm
<b>MEASUREMENT</b>	475mm	550mm	480mm	150mm	640mm	375
<b>MAX ADJUSTMENTS</b>	+50mm	+50mm	+50mm	+60mm	+120mm	+25mm
<b>MEDIUM SIZE</b> CHP231810						
<b>MIN ADJUSTMENTS</b>	-25mm	-50mm	-75mm	-60mm	-40mm	-50mm
<b>MEASUREMENT</b>	475mm	560mm	505mm	170mm	690mm	390
<b>MAX ADJUSTMENTS</b>	+50mm	+50mm	+25mm	+60mm	+90mm	+50mm
<b>LARGE SIZE</b> CHP231820						
<b>MIN ADJUSTMENTS</b>	-25mm	-50mm	-75mm	-60mm	-60mm	-50mm
<b>MEASUREMENT</b>	500mm	560mm	530mm	180mm	740mm	390
<b>MAX ADJUSTMENTS</b>	+50mm	+50mm		+60mm	+40mm	+50mm

## MODIFICATIONS

- HYGROFLEX FOAM
- LEATHER
- SEAT WIDTH
- BACK HEIGHT
- SEAT HEIGHT
- FOOTREST LENGTH
- SEAT DEPTH
- NO CASTORS
- ARM HEIGHT
- COMFORT PACK (ADJUSTABLE HEADREST AND LUMBAR)
- 180KG WEIGHT CAPACITY
- CONTROLLER LEFT SIDE
- 24 VOLT RECHARGEABLE BATTERY
- 75MM WHEELS
- PRESSURE CUSHION CUTOUT
- HEAD AND ARM REST COVERS



## FABRIC SELECTION

All VersoTela Fabrics come standard with a suite of high-end specifications, such as, antimicrobial coating, mould resistance, scratch resistance and fire resistance. Our wide range of collections each boast unique properties that are carefully designed to suit every customer's needs. Whether it be built for use in a care facility, or as a functional but stylish addition to the living room, our VersoTela Fabric options are sure to tick all the boxes when it comes to crafting your perfect chair.

Please get in touch to find out more about which fabric is right for you.

