

# ICON 30 FAF

FOR OUTDOOR ACTIVITIES



REHASENSE



**2** YEAR  
WARRANTY

Exclusively distributed by  
**aidacare™**  
HEALTHCARE EQUIPMENT

# ICON 30 FAF

## THE LIGHTWEIGHT WHEELCHAIR

ICON 30 FAF is a stable, foldable wheelchair specifically designed for outdoor activities. It's fixed anterior frame (FAF) and multiple adjustment options allow it to be tailored specifically for the users needs to maximise performance and comfort. Constructed from high-strength aluminium tubing with reinforced side frames, it provides a stable and safe connection for PAWS Power Add-ons and the Track Wheel.



**Foldable Rear Handles**  
Enables the chair to be propelled by a carer or family member



**Camber**  
Up to 6° of rear wheel camber



**Drop-in Armrests**  
Facilitate side transfers



**Footplates**  
Height adjustable to suit multiple user heights



### ACCESSORIES

**FATBOY 20"**  
Rear Wheels S  
Smooth Tread  
MWT481800

**FATBOY 20"**  
Rear Wheels  
Offroad Tread  
MWT481810

**EASYRIDER 24"**  
Rear Wheels  
Smooth Tread  
MWT481820

**FATMAN 24"**  
Rear Wheels  
Smooth Tread  
MWT481830

ICON 30 FAF	Seat Width	Seat Depth	Anterior Seat Height	Posterior Seat Height	Overall Width	Backrest Height	Rear Wheel Size	Front Castor Size	Min Chair Weight	SWL
MWC451400	400mm	400mm	435-485mm	390-455mm	630mm	310-340mm	24"	5"	13.6kg	125kg
MWC451410	400mm	450mm	435-485mm	390-455mm	630mm	310-340mm	24"	5"	13.6kg	125kg
MWC451420	450mm	400mm	435-485mm	390-455mm	680mm	310-340mm	24"	5"	13.6kg	125kg
MWC451430	450mm	450mm	435-485mm	390-455mm	680mm	310-340mm	24"	5"	13.6kg	125kg

Rehasense's quality management complies with both the ISO 9001 and ISO 13485 standards, as confirmed by TÜV NORD. Rehasense's medical devices are CE-marked and meet the requirements of EU Regulation 2017/745 (MDR).