

ICON 60

THE ACTIVE WHEELCHAIR



REHASENSE



2 YEAR
WARRANTY

Exclusively distributed by
aidacare™
HEALTHCARE EQUIPMENT

ICON 60

THE ACTIVE WHEELCHAIR

ICON 60 is an adjustable, lightweight wheelchair made from aluminium. To make the chair more sporty or comfortable, users can make many adjustments including the length of the from, seat angle, and seat position to suit their needs and preferences. ICON 60 also features a backrest that can tilt forward up to 20°, or lean back up to 10°, which can be folded down completely to transport easily.



ICON 60 with Eco Assist

Active wheelchair in ‘open frame’ design
Multi-adjustable: Seat depth, front and rear seat heights, and centre of gravity



ICON 60 with Track Wheel

Angle adjustable front forks
Angle, height, and tension adjustable backrest



Backrest options, straight closed with or without foldable handles



Basic or sports seat made of 3D breathable mesh fabric



Side panel in carbon/aluminium, with/without fenders, and optional padded armrests



Open footrest, also available in reinforced version for Track Wheel



Horizontal scissor brake, with vertical sports brake option available

ICON 60	Seat Width	Seat Depth	Anterior Seat Height	Posterior Seat Height	Backrest Height	Rear Wheel Size	Front Castor Size	Min Chair Weight	SWL
MWC451500	360mm	420-465mm	490-510mm	490-510mm	300mm	24"	5"	From 6.9kg	125kg
MWC451510	380mm	420-465mm	490-510mm	490-510mm	300mm	24"	5"	From 6.9kg	125kg
MWC451520	400mm	420-465mm	490-510mm	490-510mm	300mm	24"	5"	From 6.9kg	125kg
MWC451530	420mm	420-465mm	490-510mm	490-510mm	300mm	24"	5"	From 6.9kg	125kg
MWC451540	440mm	420-465mm	490-510mm	490-510mm	300mm	24"	5"	From 6.9kg	125kg
MWC451550	460mm	420-465mm	490-510mm	490-510mm	300mm	24"	5"	From 6.9kg	125kg
MWC451560	480mm	420-465mm	490-510mm	490-510mm	300mm	24"	5"	From 6.9kg	125kg

Rehasense’s quality management complies with both the ISO 9001 and ISO 13485 standards, as confirmed by TÜV NORD. Rehasense’s medical devices are CE-marked and meet the requirements of EU Regulation 2017/745 (MDR).