Invacare® | Action® 1NG





Simple, robust and efficient!

The *Invacare®* **Action®1NG**, although made of steel, remains a lightweight wheelchair of exceptional quality. The rigidizer blocks, commonly used with the **Action®** Family range, make the product extremely stable and, consequently, very easy to propel.

The **Action® 1NG** is available in standard configurations for immediate delivery, but Invacare also offers a wide range of accessories to meet the differing clinical needs of users and carers for optimal, everyday use.

So if you're looking for a high quality wheelchair that's practical, economical and convenient, then look no further than the **Action®1NG**.

Code: CN154XXXX

(Refer to customer service for the code size you require)

Benefits

Cross compatible to suit every need

The Invacare® Action®1NG shares multiple components with other wheelchairs from the Action®NG Family, including brakes, wheels, legrest assortment, cane holder and drip stand. That means you can adapt your **Action®1** to fully meet the user needs. Enabling the sharing of common components reduces the complexity of handling, prescribing and adjusting the product.

Manoeuvrable, and easy to transport and store

The Action®1 has been designed for optimal manoeuvrability, whether being pushed or selfpropelled, indoor or out. It even remains mobile in tight spaces or crowded areas because of its modest dimensions and compact frame. Plus, thanks to it's single crossbar system, compactness and quick release rear wheels axles, the Action® 1NG is easy to handle, transport and store. Simply take off the rear wheels and fold it away.

Simple to customize

The Action®1NG is available in both transit and self propulsion versions to accommodate individual needs. The wheelchair can easily be adjusted to the user morphology. A wide range of accessories is available to enhance patient's comfort and customize the chair to their specific requirements.

Features and Options



Transit version A transit version is available with 12" solid rear wheels for dependent users.



Legrest Fixed 80° Legrest.



Front castors 8" x 1.3/8" solid.



Detachable swing-away armrest Fixed height desk



Break lever The folding lever of push to lock brakes makes transfer easier.

Technical Data

Warranty:

2 years frame & crossbraces. Limited 1 year components.



40.5 45.5 51 cm



42.5 cm





(240mm).





41 cm



80 / Nº to 80



33 - 45.5 cm



24 cm





Transit 26.5 cm



Transit 69.5cm (inc footrest)



Self Propel 81.5/85.5 Transit 69.5



86.5/91.5cm

All product dimensions are in millimetres to apply EN 12183 standards. For more pre-sales information about this wheelchair, please consult the www.invacare.com website It can be found under the section Low Active Wheelchairs / Action1 / Sales & Images in the product catalogue. For information about typical user description, intended use. adjustments and transportation, please refer to the user



127 kg



Self Pronel 16 5/18 5/ 18.8kg Transit 16.3kg



Self Pronel 10 5/11 5/ 12.9kg Transit 14.6kg

Frame Colour



Night Blue





Black Nylon

Invacare® | Action® 2NG





The next generation has arrived

The all-new **Action 2NG** has been developed from the tried-and-tested **Action 3NG** concept.

With new cutomisation features such as depth-adjustable seat upholstery, the specifications meet the basic needs of the users while offering greater adjustability. The next generation *Action 2NG* shares the same essential platform as the *Action 3NG*, offering a low total weight for easy handling combined with good driving abilities. Due to the shared use of common components, the *Action 2NG* fits perfectly into the Action family.

Code: BQ152XXXX

(Refer to customer service for the code size you require)

Benefits

Light and comfortable

With a weight of just 15kg, the *Invacare®* **Action 2NG** fulfills all essential functional requirements. The stable folding frame means the wheelchair is easy to handle, and the light weight makes transportation trouble-free. Cross-compatibility with other members of the Action family ensures that the **Action 2NG** can be easily adapted to meet individual user needs.

Quality of life

The tried-and-tested technology of the Action family products allow the user to go wherever they desire. Cross-compatibility with *Action 3NG* make this an ideal, basic, fleet product.

Features



NEW Depth-adjustable seat upholstery

Decrease seat depth by 50mm for customised seat size.



Swing in and out footrest The 80° footrest swings in and out for easy transfers in or out of



Adjustable center of gravity The Axle position can be adjusted in height by 5cm.

Technical Data

Warranty:

2 years frame & crossbraces. Limited 1 year components.



41, 45, 48 cm



40 & 43 cm



the chair.

80°/0°-80°



48.5 cm



7° / 0 - 30°*



43 cm



35 - 48 cm



26 cm



58 5 - 68 5 cm



41 cm



95 - 109 cm



82 / 84.5 cm



93.5 cm



max. 125 kg

^{*} Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.



16.2kg Self Propel 14.9kg Transit

Frame Colour



Black

Upholstery



Black Nylon

^{*} Requires reclining back canes.

Invacare® Action® 3NG





Greater driving performance

The *Invacare®* **Action®3NG** has been re-designed to ensure greater driving performance; a lightweight frame and high quality components make for easy rolling and manoeuvrability. An appealing 'next generation' design together with improved functionalities in terms of stability and comfort means that the **Action®3NG** marks the rebirth of our flagship product in style.

Codes: CS155XXXX or A3LD-xxXxx L or R (Refer to customer service for the code size you require)

Benefits

Functionality

Invacare® always seeks to innovate. The **Action®3NG** introduces a new cross bar with a folding mechanism under the seat which improves stability and makes folding and unfolding even easier. Together with a new-style footplate, the Action®3 is a great choice for all-round maximum comfort.

Durability

Adaptable and modular, the **Action®3NG** accommodates progressive changes in user needs. It is also simple to refurbish, with a wide range of options. This model is a great example of a second-generation wheelchair, with a proven track record behind it.

Transportable

With a low transfer weight and foldable backrest, the **Action®3NG** is easy to fold and to dismantle, making it easily transportable in a small car.

Features



14kg weight

The total weight of 14kg addresses the need for improved handling, storage and transportation. The lowest transportation weight (without legrest, armrest and rear wheels) is 7.6 kg.



Brake and Brake extension

Fold the brake handle for easier side transfers.



Length compensating leg support

New option elevating legrest offers optimal length compensation with the elevation stepped in approx. 3 – 5 degrees.



Footrest

The proven design of the swing in / swing out footrest has an easy-to-operate quick release mechanism.



20 rear wheel positions

Driving performance is enhanced with a choice of 20 rear wheel positions as standard. Users can set the rear wheel to Standard, Medium, Active or Amputee position.



Magnesium castor attachment

A castor attachment provides ease of angle adjustment and enhanced support. The use of the magnesium attachment reduces the weight of the chair by almost 400 gram.

Technical Data

Warranty:

2 years frame & crossbraces. Limited 1 year components.



40.5. 45.5. 51 cm



40.5, 45.5, 51 cm



80° / 0° - 80°



38.5, 41, 43.5, 46, 48.5, 51 cm



0° / 7° / 0° - 30° /- 6° - +24°



40, 43, 51 cm



35 - 48 cm 34 - 47 cm



20 - 23 - 26 cm, 24 cm, 19 - 29.5 cm



58.5 - 68.5 cm



sw+ 20.5 cm



100 cm



72 - 85 cm



81.5 - 102 cm



max. 125 kg





14.2kg Self Propel 13.5kg Transit



7.6 kg (without armrest, footrest & rearwheels)





Anthracite Black

Upholstery



Black Nylon

Invacare Action Junior



A wheelchair that grows with the child

Invacare Action³ Junior is a lightweight, foldable paediatric wheelchair designed for children aged between 3 and 15 years. Action³ Junior has been developed to match the individual needs of the child and can grow as they grow. Because a child's needs will change the Action³ Junior offers a large range of options to accompany the child in their development and to provide the necessary clinical support. The Action³ Junior comes in attractive, modern styles, designed to appeal to kids.

Codes: A3JRSP-EVO (Self-Propel) A3JRTR-EVO (Tranist)







Grow with Invacare Action³ Junior

As children grow and develop, a key feature of a children's wheelchair is that it can be adjusted to accommodate the child's individual growth and needs. With the *Invacare Action*³ *Junior*, the width of the seat and backrest can be adjusted, as can the height of the arm and footrests.

An individual fit with extensive size combinations

When choosing a paediatric wheelchair, a key factor is ensuring a perfect and comfortable fit for the child. The **Action³ Junior** comes in one of the most extensive ranges of size combinations with seat width from 205 to 380 mm, offering the right fit for all children, from toddlers to teenagers.

Customise your chair!

A wide range of options is available on the **Action³ Junior** to meet every child need in term of comfort and positioning. In addition the child can customize his wheelchair and appropriate it thanks to the wide range of frame colours, spoke guards and other options.

Features and Options



Seat-width adjustment

The seat width of the chassis can be adjusted with an innovative telescopic crossbar over 4 seat positions from 205 to 380 mm.



Upholstery adjustment

As the chassis is made wider, so the upholstery on the back and on the seat of the chair can be made wider with Velcro flaps.



Table tray

In addition to the standard slide-on table tray, we offer a swing-away tray which can easily be folded away. Code: 1492956 - Acrylic



Adjustable seat height and balance

For improved driving performance and different height options, the rear wheel can be mounted in 20 different positions.



Front stabilizer kit

A front stabilizer kit is available to increase forward stability for active children who are likely to climb all over the wheelchair and be especially energetic when propelling the wheelchair.

Code: 1531160



Spoke guard

Four different types of spoke guards are offered to match the different wheelchair colours and tastes. (#4 = transparent)
Code: 1556990

Features and Options



Hangers

① 90° hanger, ② 80° hanger and ③ elevating legrest adjustable from 0° to 80° are available.

Codes: Refer to Script Form



Armrests

Removable and height adjustable desk ① and junior ② armrests as well as skirt- ③ and mudguards ④ are available.

Codes: Refer to Script Form



Abduction legrest position

A 30° abduction position is possible with 80° hangers.



Transit Version

We also provide a transit version, with a seat height of 460 mm and 12" plastic rear wheels.



Front wheel assortment

3 different types of front wheels are available: ①5" x 1 ¾" solid tyre and ② 5" x 1" solid tyre.
Codes: Refer to Script Form



Cushion

A 50 mm bi-density foam cushion is available and follow depth adjustment of the seat upholstery.





Ergonomic push handles

The push handles can be adapted to all kinds of backrests and are adjustable in a variety of ways. A width-adjustable back brace can be added (or used on its own) to make the backrest more rigid. It can also be utilized as a headrest support, ensuring greater user comfort.

Codes: Refer to Script Form

Technical data



20.5 / 23/25.5 / 28 / 30.5 /33 / 35.5 / 38 cm



20 - 25 / 20 - 30 / 30 - 35 / 27.5 - 32.5 / 32.5 - 37.5 cm



1.5° / 5° / 8.5°



40 / 43.5 / 46 cm



0° / 7° / 0° - 31



35 - 43 cm



80°/90° 0°-80°



12 5 - 45 cm



15 - 25 / 18 - 28 cm



+ 19 / + 17* cm



31.5 / 27* cm



74 - 96.5 mm



51 - 74 cm



87 - 113 cm



80ka



15 / 14.5* kg



8.5 / 12* kg



^{*} Transit Version

Frame colours

Standard Colours







Sapphire Blue

Please note colors may vary slightly from those displayed here.



Upholstery



Black Nylon

