



**MANUAL
CHAIRS**

ADULT								PAEDIATRIC		
Folding				Rigid			Tilt in Space	Folding		
Kuschall Compact	Kuschall Ultralight	Kuschall Champion	Kuschall Advance	Top End Crossfire T6	Kuschall KSL	Crossfire T7A	Solara 3G	Action 3 Junior	Freedom NXT	
Seat Width (mm)	280 - 500	280 - 500	340 - 460	340 - 460	300 - 460	340 - 420	300 - 500	305 - 610	300 - 360	255 - 510
Seat Depth (mm)	320 - 500	320 - 500	340 - 460	360 - 520	360 - 460	375 - 475	360 - 500	305 - 560	200 - 360	305 - 510
Floor to Seat Ht (mm)	380 - 530	430 - 530	430 - 540	450 - 550	410 - 530	450 - 510	400 - 530	320 - 480	410 - 460	320 - 495
Max User Wt (kg)	130	130	120	125	113	100	113	113 - 118	80	91
Weight (kg)	11	9.9	10.7	7.5	8.8	6.7	8	15	7.2	10



**POWER
CHAIRS**

ADULT					
Mid Wheel Drive			Rear Wheel Drive		
Rovi X3	TDX SP	TDX SR	3G Torque	3G ARROW	
Seat Width (mm)	280 - 550	300 - 600	300 - 600	300 - 550	300 - 610
Seat Depth (mm)	250 - 580	300 - 550	300 - 550	300 - 550	300 - 560
Floor to Seat Ht (mm)	410 - 430	400 - 450	420 - 520	420 - 530	440 - 530
Max User Weight (kg)	136	136 (Std), 180 (HD)	136 (Std), 180 (HD)	136 - 180	159 (Std), 181 (HD)
Base W x L (mm)	600 x 1150	610 x 900	650 x 900	630 x 800	630 x 800
Max Speed (kph)	8	9	11.2	9	10
Features	Patented ARC (Active Ride Control) suspension. Full power seat functions available	Stability Assist, Sure Step, G-Trac	Stability Assist, Sure Step, True-Trac	Rear-wheel suspension. Full power seat functions available	Rear wheel suspension, gearless & brushless motors. Full power seat functions available

CLICK ON THE PRODUCT IMAGES TO BE DIRECTED TO MORE INFORMATION ON THE AIDACARE WEBSITE