



Junior+ by Jenx

A comfortable, contoured and stable modular seating system, designed to suit the needs of children between 18 months and 12 years



Always Adaptable...

Jenx Limited, Wardsend Road, Sheffield S6 1RQ, United Kingdom tel: +44 (0) 114 285 3376 email: export@jenx.com www.jenx.com



Junior+ by Jenx

The Junior+ offers highly flexible supportive seating for children from 18 months - 12 years.

An outstanding level of versatility is achieved through a wide variety of head, trunk, leg and foot support options as well as a choice of bases, meaning each Junior⁺ can be exactly tailored to the needs of any individual child. The contoured seat unit constructed from super-soft polyurethane and memory foam provides unparalleled comfort and the perfect base for true pelvic stability.

Created in conjunction with children, parents, therapists and care staff, the Junior⁺ has been thoughtfully and carefully designed to meet a wide range of needs - meaning every child can be comfortable and perfectly supported, allowing them to get on with playing and learning.



Outstanding Pelvic Control System

Contoured seating with memory foam seat pad, independently adjustable lumbar/sacral support and a 4 point lap strap all provided as standard make the Jenx Junior+ outstanding when it comes to pelvic control. Plus if you find this isn't enough, the optional **Ischial Bar** will help prevent low tone or very active children from slipping into sacral sitting.

INNOVATIVE FEATURE

PU Support Pads

Included as standard for hip support with easy, flip-away action for leg support. All pads are made from super-soft polyurethane – fully flexible, easy to clean and impervious to fluids, making them hygienic and robust.

INNOVATIVE FEATURE

Adjustable Seat Depth

To ensure the fit is exactly right for each child, the seat depth is fully adjustable. This is a tool free, quick and easy to access adjustment.

INNOVATIVE FEATURE

Leg Control Pads

Use 2 or 4; these unique **Leg Control Pads** allow you to control leg positioning independently, creating abduction or adduction or controlling windsweeping. Simple to adjust, they flip completely out of the way for smooth, safe, dignified and unobstructed transfers. INNOVATIVE FEATURE

Flip-up Angle Adjustable Footplate

To accommodate or encourage dorsi or plantar flexion. Adjusts quickly and simply with the turn of a button, flips up for safe and easy standing transfers with the footplate held out of the way with a concealed magnet. Now with easy fit slots making attaching sandals and other foot support options a breeze.

Improved Shoulder Protractors

Comfort style **Shoulder Protractors** to encourage a functional mid-line hand position and prevent excessive retraction.

Improved Flexible Thoracic Supports

. .

Designed with comfort in mind these super-flexible thoracics, available in different sizes, hug to provide just the right amount of support without restricting function.

NEW Super-soft PU

New, super-soft PU - still the same great hygienic, easy-to-clean seat cushion but now with added comfort.

NEW Ankle Snug

An optional foot positioning accessory to the standard sandal. The **Ankle Snug** offers a unique magnetic buckle fastening mechanism, which can be operated single handed.

Always Supportive

Every position a child adopts has an effect on their development and functional movement now and in the future.



Always Configurable Accessories

The Junior+ by Jenx has a wide range of accessories and optional extras meaning it can be highly customised for every child's needs, making it a highly flexible supported seating system for children from 18 months to 12 years old.

Junior+ Seat Unit

Suitable for children from approx 18 months – 12 years across 3 sizes, each Junior⁺ seat unit provides all the basics for a solid seating foundation, including:

- Contoured PU seat with removable abduction pad and memory foam seat pad
- Adjustable, flip-away PU hip pads
- Angle adjustable & flip-up footplate
- 4 point lap strap
- Height adjustable PU back pad
- Independently adjustable PU lumbar/sacral pad
- Leg control brackets

X-Base Black

A compact, robust and cost effective base with high-low and tilt-in-space functionality. With lockable, swivel castors and a height range of 345-515mm. Tilt-In-Space 30° recline, 10° prone.



X-Base White

The white based version of the X-Base offers a higher height range of 395-595mm, meaning it can now allow children to engage with family members and friends at dining table height with ease. Tilt-in-space 30° - 10° prone.



Y-Base

The most flexible base in the Jenx range, the Y-base offers a huge range of height adjustment from 335 - 755mm alongside tilt-in-space functionality and lockable, swivel castors. Tilt-in-Space 30° recline, 10° prone.



Multigrip[™] Head Support

The original and innovative Jenx **Multigrip™ Head Support** is available in 3 sizes for all Jenx seating systems. Provides adjustment for height, depth and angle.



Oval Head Support

Simple headrest available in a range of colours for children who require minimal head support. Comes mounted on a flexible fixing bracket giving adjustment for height, depth and angle.



Adjustable Push Handle

Makes moving equipment safer and easier for carers, flips away in seconds when not in use so ideal for busy classroom environments or where space is at a premium. Has a comfortable padded foam handle and is height adjustable.



Medium Back Pad

This accommodates larger children who fit the length of the traditional Junior+ but need extra width for comfort and support. Standard on size 3 Junior+.



5-port Back Extension Post

Gives extra height to the **Back Pad**, for those who require it.

Shoulder Protractors

Comfort style shoulder protractors to encourage a functional mid-line hand position and prevent excessive retraction. Completely flip-away when not required. Independently adjustable for width and angle. Available in 2 sizes.



Tray

A useful surface for work, feeding or play. Made from food grade plastic providing a hygienic and easy to clean surface. Contrasting tray edging gives a definitive boundary particularly useful for users with visual impairments. Plus, enables toys or other items to be securely fastened in place. Available in two sizes. Requires armrests for use.

NB: Suction toys, anchors or grab bars will not work with this tray – to use these items with the Junior⁺, you will also require a **Reduced Cut-Out Tray Infill**.

Tray Elbow Blocks

Used to encourage arms into a mid-line position and/or to block excessive retraction.



Reduced Cut-Out Tray Infill

Designed to work with the standard Junior⁺ **Trays** by reducing the aperture for smaller children. Perfect surface for use with suction toys, anchors or grab bars. Prevents even the slimmest of arms from falling into gaps.



Clear Tray Overlay

Ideal for use with communication aids, makaton symbols and other teaching/learning aids. Simply sandwich the teaching aid between the tray and clear overlay to provide a protective, waterproof cover.



Tray Edge Padding

Extra cushioning for the tray aperture. Attaches with hook and loop fastening and is covered in wipe-clean PU fabric. Ideal for protecting bony elbows and arms from banging on the tray edge or support surface.



Full and Half Tray Padding

Extra cushioning for the tray. Attaches with hook and loop fastening and is covered in wipe-clean PU fabric.



Activity Frame

Great for encouraging head control and hand/eye co-ordination. Easy to attach and adjust. Comes with two clips for attaching toys or other items.



Grab Bar

Perfect for children who need to fix one hand to increase functionality.



Armrests

Provided as a pair the **Armrests** are independently adjustable for height and angle. Clear markings make it easy to set the **Armrests** at the correct height. Padded armrest covers in a range of colours are also available.



Seat Extension

Provides extra length on one side of the seat for children with leg length discrepancy.



NB. Cannot be used in conjunction with **Leg Control Pads**.

Ischial Bar

Simple yet highly effective positioning aid for children with a strong tendency towards sacral sitting. Position beneath the memory foam pad in front of the ischial tuberosities.

Waistcoat Harness

Provides additional trunk support to children with low tone. Attach, detach and flip-away in seconds with a single hand – no tools required. Available in 4 sizes.



Flexible Thoracic Supports Pads

Available in a range of sizes, these slimline and super flexible **Support Pads** (with straps) are designed to provide outstanding positioning with comfort in mind. All sizes feature single handed flip-away action and are independently adjustable for height and width to ensure maximum trunk control.



Also available: High-rise Flexible Thoracic Supports, offering a greater range of height adjustment.

Thoracic Support Pads

Available in a range of short, long and slimline sizes (all without straps) with regular, or high/low rise brackets, these trunk support pads are made from flexible and hygienic PU. There is a single handed push button flip away action on all



supports and they are independently adjustable for height and width. Can be offset to brace a scoliosis or accommodate other skeletal deformities.

Hip Width Reducing Pads

Designed to effortlessly slot in between the hip pad and cover, the infill reduces the hip width by 50mm, accommodating a narrower hip-width.



Abduction Blocks

Available in three sizes providing different levels of abducted seating.

Leg Control Pads

Unique Leg Control Pads

allow you to control leg position independently of pelvic position, creating abduction or adduction or controlling windsweeping. Flip completely out of the way to facilitate smooth, safe and dignified transfers.



Box Style Footplate

Perfect for children whose condition precludes them from having their feet fixed. Provides a level of control while still allowing freedom of movement.



Additional Footplates & Brackets

Provides support for feet. 3 different sizes of **Bracket** allow for different leg lengths.



Footplate Pad

Attach to the **Footplate** to give a softer surface for children who don't need or can't tolerate more structured foot positioning.



Sandals

Available in 4 different sizes to accommodate any size of foot, even when wearing splints, Jenx **Sandals** are interchangeable across the entire range of Jenx seating and standing products.



Ankle Snugs

An optional foot positioning accessory to the standard sandal. The **Ankle Snug** offers a unique magnetic buckle fastening mechanism, which can be operated single handed.



Booties

Available for the smallest Jenx seating systems, **Booties** provide gentle and comfortable support for tiny feet even without shoes. They are soft and extra padded with no fiddly straps, press studs or buckles.



Sandal Raising Blocks

Available for every size of Jenx **Sandals**, the raising block kits give the **Sandals** that extra bit of height where needed! Use to accommodate for leg length discrepancy or use to raise the height of the **Footplate** on any product to maximise the available growth.



NB: One kit needed per individual sandal.

Washable Padded Covers

A full range of colour co-ordinated covers are available for the seat and accessories. Easily removed and machine washable at 40°. Available in Red, Blue, Pink, Green, Yellow or Grey. Optional memory foam padding is also available for additional pressure relief where needed. A wipe-clean PU seat hygiene cover is also available in black.



visit jenx.com for more information



Junior⁺ by Jenx Technical Specifications

IMPORTANT INFORMATION: Please be aware, some minimum and maximum measurements are achieved by the use of additional accessory components.	Measurement	Units	Junior+ 1 JUC21	Junior+ 2 JUC22	Junior+ 3 JUC23
Age Range	approx.	mths - yrs	18mths - 5yrs	3yrs - 8yrs	5yrs - 12yrs
User Weight	max	kg	50	50	50
Back Angle	recline - forward	degrees	15° - 15°	15° - 15°	15° - 15°
Back Pad Height from Seat to Top	min - max	mm	345 - 435	345 - 435	410 - 500
Chest Width	min - max	mm	185 - 280	185 - 280	210 - 300
Hip Width	min - max	mm	150 - 375	150 - 375	150 - 375
Seat Depth	min - max	mm	210 - 305	240 - 355	260 - 440
Lower Leg Length	min - max	mm	195 - 270	235 - 350	250 - 405
Footplate Angle	dorsi - plantar	degrees	5° - 35°	5° - 35°	5° - 35°
X and Y Base Tilt-In-Space Angle	supine - prone	degrees	30° - 10°	30° - 10°	30° - 10°
X Base Footprint	width x length	mm	625 x 625	625 x 625	625 x 625
Y Base Footprint	width x length	mm	660 x 700	660 x 700	660 x 700
X Base Black Seat Height From Floor	min - max	mm	345 - 515	345 - 515	345 - 515
X Base White Seat Height From Floor	min - max	mm	395 - 595	395 - 595	395 - 595
Y Base Seat Height From Floor	min - max	mm	335 - 755	335 - 755	335 - 755

Find Your Distributor

To locate one of our carefully selected distributors near you, visit: jenx.com/where-to-buy





JENX Positioning for life

Jenx Limited, Wardsend Road, Sheffield S6 1RQ, United Kingdom tel: +44 (0) 114 285 3376 email: export@jenx.com www.jenx.com