









Monkey by Jenx

A simple, yet highly adjustable standing system, that provides prone to upright positioning for children from 9 months to 4 years old



Always Supportive...



Monkey Prone Stander by Jenx

Specifically designed to appeal to young children and their families, the Jenx Monkey is an ideal fun and gentle introduction to standing. Offering smooth and easy angle adjustment, even with the child in the product, the Monkey can go from upright at 90° right down to 20° prone – perfect for gradually building tolerance to standing as part of a therapy programme. The capacity for outstanding thoracic and pelvic support coupled with flexible leg positioning make the Monkey prone stander a great option for almost every child.



De-rotational Thoracic and Pelvic Straps, with Lateral Support Pads.

Straps fasten in the centre and adjust independently from either side, enabling a mid-line position to be achieved and maintained. Buckles and Velcro mean that, once adjusted correctly, the straps can be fixed making it difficult for siblings or classmates to tamper with and keeping everything looking neat and tidy!

Integral Angle Adjustable Tray

Provided as standard, with both bowl and smooth tray infill for a wide variety of different activities.

INNOVATIVE FEATURE

Symmetrical Adjusting Hip Straps

Fastened in the centre and then adjusted from either side, which keep a child's hips in mid line. This is the key feature that keeps your child close to the standing support.

Removable Abduction Block

Provided as standard to help create and support an abducted standing posture.

INNOVATIVE FEATURE

Sandal Straps

Sandals fully open at the back, making self transfers easier for the child and carer.



Independently Adjustable Thoracic and Pelvic Boards

Allowing the exact proportions of each individual child to be catered for, whilst also giving scope to avoid pressure on delicate areas such as feeding buttons.

Always Encouraging...

Every position a child adopts has an effect on their development and functional movement now and in the future.

Improves blood pressure & circulation





Aids respiration & speech

Helps digestion, bowel & bladder emptying



0

Helps relieve pressure created during sitting

Encourages bone & muscle development





Supports pelvis, thighs & feet so children can join in everyday activity

Assists freedom to explore & participate





Provides new independence & increases social skills

Children gain huge psychological benefits from standing, as well as important physical and development changes.

Smooth & Simple Angle Adjustment

Angle can be adjusted with the child in the standing frame.



Always Thoughtful Outstanding features and benefits

Bright and Friendly

The first diagnosis of a child's special needs is very difficult for parents who have to learn a whole new vocabulary and deal with equipment and treatments that are new and bewildering. With this in mind, the Jenx Monkey has been carefully designed to offer the best in early intervention support in a gentle, bright and friendly way. Parents often comment on the huge difference this thoughtful design has had in helping them accept and use standing equipment for their very young children.

Angle Adjustment

Each child will have a slightly different posture in standing and need slightly different support. The angle of the support surface can be adjusted from 20° prone through to upright, 90°. And, because the adjustment is infinite, Monkey can be set to exactly the right angle for each individual child. This means a child's progress can be monitored and the angle of the board altered in line with changes in their ability.

Comfortable Upholstery

Fully upholstered chest, hip and leg boards help make standing more comfortable for children, which helps to increase their tolerance of standing. A happier and more comfortable child will spend longer standing and more time reaping the benefits of all that standing offers.



Flexible Growth

The chest and hip boards can be independently altered and fixed allowing you to accommodate the exact proportions of each child.

Adjustable and Safe Hip & Chest Strap

The innovative straps close in the centre and are adjusted at each end for length. By adjusting the straps in this way you can ensure the child is positioned symmetrically – helping to keep them upright and comfortable and maximising the benefits of standing.



Quick, Simple and Easy to Use

There are many reasons why the Jenx Monkey is probably the most popular infant prone stander available! Monkey comes as a complete unit meaning that what you receive out of the box is ready to be quickly and easily set up (with very few adjustments required) so your child can start using it straight away!

Visit **jenx.com** for additional product details.

Always Configurable Accessories

The Jenx Monkey standing frame comes complete with:

- · Chest, hip and leg support boards
- Chest and hip straps
- Footboard
- 4 lockable, swivel castors

Sandals

Available in 4 different sizes to accommodate any size of foot, even when wearing splints, Jenx sandals are interchangeable across the entire range of Jenx seating and standing products.



Straps open from the back facilitating easy stand transfers into prone standing.

NB: Monkey comes with size 1 sandals as standard – larger children or children wearing splints may find size 2 sandals more appropriate. Sizes 3 & 4 are not appropriate for use with Monkey.

Grab bar

Perfect for children who need to fix one hand to increase functionality.



- · Adjustable chest, hip and leg positioning pads
- Intergrated angle adjustable tray with bowl
- Size 1 sandals

Sandal Raising blocks

Available for every size of Jenx sandals, the raising block kits give the sandals that extra bit of height where needed! Use to accommodate for leg length discrepancy or use to raise the height of the footplate on any product to maximise the available growth.



NB: One kit needed per individual sandal.

Activity Frame

Great for encouraging head control and hand/eye co-ordination. Easy to attach and adjust. Comes with two clips for attaching toys or other items.



Always Adaptable





Jenx Monkey Technical Specifications

	Units	NR-50 Prone
Angle Range	degrees	20°-90°
Age Range (approx)	mths/yrs	9mths-4yrs
Base Size	mm/inch	575 - 880 / 22½ - 34½
Chest Height (footplate to top of chest pad) (min/max)	mm/inch	510-750 / 20 - 29½
Chest Pad Width	mm/inch	150-235 / 6 - 91/4
Hip Pad Width	mm/inch	150-235 / 6 - 91/4
Max User Weight	kg/lbs	25 / 55
Tray Depth	mm/inch	360 / 14
Tray Width	mm/inch	530 / 21

Always here

To arrange a demonstration or assessment with the Jenx Monkey standing system, just get in touch!

My name is:

Please ring me on:

Or email me at:

