KERA SIT2SIT HOIST A revolution in Manual Handling





Solve your transfer issues with the Kera

The Kera's unique design uses an easy to fit backstrap and chestpad to securely hold a non weight bearing person. Now it's possible for a single caregiver to quickly and safely transfer and perform hygiene processes for people that have no standing strength or balance.

Give time, energy and dignity back to people living with a disability by eliminating the need for multiple transfers during hygiene processes. The Kera SIT2SIT provides a comfortable and safe option thereby eliminating the need for traditional hoists or manual transfers.

The Kera SIT2SIT has improved the lives of hundreds of people living with disability for more than five years.

- No holding on, no standing/sitting strength or balance required
- Only one caregiver needed
- No weight bearing ability needed
- Users are held firmly by the upper body making the SIT2SIT ideal for amputees and tetraplegics
- Replacement for a standing hoist or sling hoist

Contact your Aidacare representative 1300 133 120 | aidacare.com.au



Exclusively Distributed in Australia by



FEATURES AND BENEFITS





•

Up to 120kg Fits through doors >610mm

h Fits around m chairs up to 800mm wide

Proven by people living with MD, MND, MS, stroke, amputations, spinal injuries and many other conditions







Patented and innovative transfer position

- Access to clothing eliminates multiple transfers and enables hygiene processes
- Manual operation means no waiting to charge

Unique moving knee supports

- Takes pressure off the knees during the lift
- Helps to reposition deep into a chair

Unique Chestpad and Backstrap combination

- Securely lifts by the upper body, no holding on, no standing/sitting strength or balance required
- Can replace a standing hoist or hold off progression to a sling hoist
- No sling to fit under the legs for quick transfers between chairs
- 4 x configurable options for customisation to the user and environment
- Accessory range to adapt to user ability and size, including children
- Independently laboratory tested to international hoist standards (AU/NZS 10535)
- Entered in the ARTG (388973)

ACCESSORIES RANGE

LSS398350	Backstrap extenders - (pair)
LSS398360	Hook Lowering Bracket
LSS398370	Chest Pad Assembly - Standard - Configurable (includes cradle & slider bar)
LSS398380	Chest Pad Assembly - Short - Configurable (includes cradle & slider bar)



Steve's Story

Steve has MND. He has limited mobility and is dependent on those who care for him to go about his daily life. He sees the Kera as a way to stave off the progression to the inconvenience of a full hoist.

Steve has a full life, filled with family and adventure. Having the Kera has allowed him to travel far more easily. Most importantly, the Kera gets him on to his scooter where he can enjoy the freedom and independence of a ride around the neighbourhood.

Find out more with videos of Steve and other similar stories at htsystems.com.au

Contact your Aidacare representative 1300 133 120 | aidacare.com.au



Exclusively Distributed in Australia by

