

BEDSIDE CABINET & OVERBED TABLE



OVERBED TABLE

- Lightweight height adjustable overbed table with two locking castors
- Extendable tray surface ensures meals and drinks can easily be enjoyed in bed
- Smooth height adjustment is easily accessible to users and carers
- Tray removes for easy cleaning and management of dishes
- Extendable tray offers compatibility with integrated drive systems, maintaining full patient access to the table-top





BEDSIDE CABINET

- Dual-fronted bedside cabinet with ample storage for personal belongings
- Raised edges on the top surface for fluid retention and to prevent spills
- Pull out towel rail
- Provision to add draw lock (lock sold separately - BEA828200)

| Product | Code | Width | Depth | Height |
|-----------------|-----------|-------|-------|--------------|
| Bedside Cabinet | BEA827165 | 500mm | 480mm | 770mm |
| Overbed Table | BEA827170 | 630mm | 410mm | 880 – 1150mm |

NESTING BASE DESIGN



Powder-coated steel with plastic finishings, wipe clean and ideal for infection control management

Nesting base design, reduces overall floor space when not in use