

BRERA RX PATIENT TROLLEY



ARTG - 295859 *Pictured with Optional Mattress

FEATURES AND BENEFITS:

- Foot operated hydraulic lift, tilt, back raise, Trendelenburg and reverse Trendelenburg functions
- Full length sliding X-ray deck ensures complete spinal coverage
- Quick Release CPR with access from both sides for backrest and knee break ensures the trolley is emergency ready at all times
- Collapsible vertical side rails fold away below the patient surface
- Fold-away carer push handles to maximise workspace
- 200mm castors minimise rolling resistance with central and directional locking functions
- Molded base cover facilitates cleaning and supports infection control processes
- Large corner bumpers minimise wear and tear on equipment and facilities
- Frame is finished with Biocote, a technologically advanced surface coating that promotes infection control via its anti-bacterial properties





OPTIONAL ACCESSORIES



IV Pole (TSE674130)

Oxygen Bottle

Holder (TSE674110)



Folding Monitor Tray (TSE674120)



Mattress (TSE674315)



BIOCOTE Antibacterial frame coating for improved infection control





Code	TSE674300
Overall Length	2100mm
Overall Width	830mm
Overall Height	560 - 880mm
Patient Surface Length	2000mm
Patient Surface Width	650mm
Backrest Angle	0° - 90°
Trendelenburg	-17° - 17°
Castor Size	200mm
SWL	270kg