

### THE BED I WAS WAITING FOR







### THE EVOLUTION FROM HOSPITAL BED TO WORKING TOOL

**DELTA4** is the new generation of beds developed according to the needs of healthcare professionals, with an idea in mind: if the bed allows the operators to work at their best, the patient's well-being follows naturally. With **DELTA4** the bed becomes an active platform that simplifies and maximises the performances of healthcare professionals.

# A SOLUTION FOR EACH DEPARTMENT

Working alongside people who work in critical situations, we realised how crucial the bed is.

We asked ourselves, how much can care improve if the bed increases the operator's skill? Designing a bed around a caregiver is no longer a limitation, but a powerful and versatile tool. From this fieldwork we have made innovations that make the **DELTA4** bed the perfect ally in all care conditions.

"DESIGNED TO MAKE ME WORK BETTER"

# **AN EVOLUTIONARY** LEAP THAT STARTS WITH PROFESSIONALS

# **DESIGNED FOR THE MOST** DEMANDING SITUATIONS,





ongoing problem.

optimal for ventilation.



### How much time and effort can you save if the patient doesn't need continual attention in bed?

- Those who work in this field know this well: the patient sliding towards the foot of the bed is an
- The continuous slipping and repositioning cause discomfort to the patient and produce a friction that in the long run can lead to pressure injuries. Futhermore, repositioning involves not one, but two people at a time, in strenuous and dangerous
- operations for both caregiver and patients.
- We solved this critical point with the INCENTRO
- SYSTEM: the combined back and pelvic translation nullifies horizontal forces, prevents slipping and keeps the body at a 30° angle,
- A patient who does not slip reduces fatigue and the risk of injury, and frees the staff to focus on more useful activities for the patients.

### "MY PATIENTS ARE SAFE EVEN WHEN UNATTENDED"

### SAFE EXIT



Safe exit position button brings the bed to the minimum height to help the patient stand up. The patient can autonomously exit from the bed or be attended by the nursing staff.









### **NIGHT LIGHT LED**

The patient can safely leave the bed during the night without the help of the nursing staff, thanks to the night light.

### **HIGH CONTAINMENT SAFETY SIDES**

High containment safety sides (42 cm) with dampened descent.

(M)MALVESTIO<sup>®</sup>

### **SAFETY HEIGHT ALWAYS EVIDENT EVEN IN THE DARK**

Overnight I can immediately verify if the bed height guarantees a safe and autonomous bed exit for the patient.

The orange light indicates the unsafe mattress platforn height, the green light indicates safe mattress platform height.





### **"BED EXIT DETECTION:** MY EYES EVEN WHEN I'M NOT THERE"



The sensor system monitors the patient's presence and their possible exit by alerting the nursing staff with an acoustic and luminous alarm. On request, the warning can be sent remotely. The desired parameters can be set up through the touch screen on the side rail. (System available with active pressure sensors).

### **FALLS FROM BED REPRESENT MORE** THAN 60% OF **TOTAL FALLS**

### **AUTOMATIC SELF-BRAKE**

Once connected to the power net, a sound alarm alerts the nursing staff about the "un-braked" bed and after 20 seconds the brake is automatically activated.



20'

# PATIENT WEIGHT UNDER **CONTROL AT ALL TIMES**



The exclusive weighing system was designed to guarantee a high level of precision independently from the mattress platform's position.

Thanks to the 4 sensors located at the corners of the bed, the load is better distributed and the reading level more accurate.

The weighing data can be recorded manually whenever necessary or automatically every 24 hours.

While the weight of the patient is monitored it is possible to set up different types of alarm: exceeded threshold weight alarm or "bed exit" alarm.









### " INCREASED MANEUVERABILITY: I CAN MOVE THE PATIENT SAFELY AND EFFORTLESSLY "





With the user-friendlyTraction Control Switch you can control both the forward/ backwardmovement andwalking speed.

### MOTORISED 5<sup>TH</sup> WHEEL



Thanks to the support, the ascent and descent brake assist of the **DELTADRIVE** moving system, I can easily face ramps with up to a 10% incline.







The system is equipped with an anti-crushing sensor that ensures the operator's safety during backwards maneuvers.

### A MULTI-FUNCTIONAL WORK PLATFORM



#### 13



### "VERSATILE & SAFE: EVEN THE MATERIAL HELPS ME"

### X-RAY EXAMINATION WITH NO NEED TO MOVE THE PATIENT

Wide area of radiotransparency. The HPL radiotransparent backrest has an X-Ray cassette holder for chest x-ray examinations, without moving the patient and at any angle of the backrest.



SOLID YET LIGHT MATTRESS PLATFORM



The mattress platform can be easily removed; all sections are designed for easy cleaning.

(MMALVESTIO<sup>®</sup>



SMOOTH AND ROUND SURFACES WITHOUT SHARP EDGES

14

The X-Ray tray can also be pulled out with backrest in horizontal position, avoiding patient movement.

### "THE BED THAT SIMPLIFIES MY JOB"







Minimum height.

Maximum height.



The antishock position is reached quickly, with the simultaneous lowering of the backrest, zeroing of the mattress platform sections and Trendelenburg positioning.



Reverse Trendelenburg position.



The examination position height adjustment facilities medical staff in the examination, medication and monitoring operations, by standing

The CPR adjustment allows to bring the mattress platform immediately down to the horizontal.

## FUNCTIONAL FEATURES





**4 SECTOR SAFETY SIDES WITH ELECTRIC CONTROLS** The system of 4 sector safety sides with integrated controls allows full autonomy of movement to the patient.



#### **INCENTRO SYSTEM**

The combined back and pelvic movement nullifies horizontal forces, prevents slipping and keeps the body at a 30° angle, optimal for ventilation.



**BED LENGTHENER** The bed lengthener can be controlled by just one person.











#### SAFE EXIT

The safe exit control moves the mattress platform to the minimum height and lifts the backrest up to facilitate the patient to get up from the bed.



#### CPR

The CPR control on both sides allows the bed platform to be brought to the horizontal position immediately; accessible even with the side rail lowered



#### REMOVABLE HEADBOARD AND FOOTBOARD

The bed ends can be easily removed for emergency operations.



WATCH THE VIDEO!







**BRAKING SOUND ALARM + ORANGE / GREEN LIGHT** The orange light indicates the unsafe mattress platform height, the green light indicates the safe one. Sound alarm for wheels not breaked.



AUTOMATIC SELF-BRAKE Once connected to the power net, a sound alarm alerts the nursing staff about the "un-braked" bed and after 20 seconds the brake is automatically activated.



**ELECTRIC FOOTREST** Electric adjustable footrest.



**BLANKET AND "SUPERVISOR" HOLDER** Extractable tray for blankets and supervisor.





X-RAY TRAY The X-ray tray together with HPL backrest allows extremely easy x-ray operations to be carried out

The INCENTRO SYSTEM is suitable for any type of hospital mattress. By using a dedicated mattress, specific functions are enhanced





All LIFECOMFORT® Mattresses are fitted with superior, hospital-grade polyurethane coated covers designed to withstand the rigors of usage and cleaning requirements in Home Care, Aged Care and up to Acute Care environments. Multiple foam layering combinations provide maximum comfort, support and pressure redistribution across a wide range of user shapes and sizes.

### 🛷 aspire **ACTIVE AIR<sup>®</sup>**



Acute 8 Active Air Mattress

Bari 10 Mattress

Pressure care alternating air mattresses are designed to offer comfort and support for patients at high risk of pressure sores. The alternating air mattress provides increasing levels of support for medium to very-high-risk acute care patients.

For more information on our static and dynamic mattress ranges please see our website http://www.aidacare.com.au or Call 1300 133 120 to be connected with your local branch.

**TURNING CYLINDER HOLDER** 

### Aspire LIFECOMFORT® Static Mattress Range

### Aspire ACTIVE AIR® Dynamic Acute Mattresses



Upon delivery of the **DELTA4** beds, Aidacare professionals fully train the staff in taking full advantage of the features of the new work tool. The support of specialized personnel covers all areas of use of the bed in its various phases of use: from on-call assistance for any need to scheduled and preventive maintenance to ensure that your **DELTA4** bed is always in full working order.

### "ALL THE SUPPORT AND TRAINING I NEED'



### **PRODUCT INFORMATION**

BEB050140 - DELTA 4 WARD BED BEB050150 - DELTA 4 WARD BED W/DELTADRIVE & WEIGH SCALE

DELTA 4 beds come in two variations and have a range of accessories available to further improve their performance.

	B
Feature	Delt
Incentro System	Ì
150mm Castors	
DeltaDrive Electric 5th Wheel	
Weighing Scale	
Bed Lengthener	
Blanket Holder	
XRay Tray	
4 Section Side rails with electric controls	
LED Night Light	
Automatic Self Brake	
Safe Exit	
Removable Bed Ends	
Removable Mattress Platform	
CPR	
Brake Alarm (Yellow/Green lights)	
Backlit Nurse Control Panel	
Ocygen Bottle Holder	
Radio Translucent Mattres Platform	
Urine Bottle Holder	
Accessory Bar Holder	
Accessories (available seperately)	
Self Help Pole (BEA825090)	
IV Pole (BEA825101)	



EB050140	BEB050150
a 4 Ward Bed	Delta 4 Ward Bed W/Delta Drive & Weigh Scale
$\checkmark$	√
$\checkmark$	$\checkmark$
-	$\checkmark$
-	$\checkmark$
$\checkmark$	√
$\checkmark$	$\checkmark$
$\checkmark$	√
$\checkmark$	$\checkmark$
$\checkmark$	



Exclusively distributed by:







The Company reserves the right to introduce all the technical and/or aesthetic modifications deemed necessary to their products without prior notice.