



ELEKTRA PATIENT TROLLEY

WITH DELTA DRIVE



*Pictured with Optional Mattress

FEATURES AND BENEFITS:

- Foot operated hydraulic lift, tilt, back raise, Trendelenburg, and reverse Trendelenburg functions
- Full length sliding X-ray deck ensures complete spinal coverage
- DeltaDrive electric centre 5th wheel with integrated controls in headboard provides powered maneuverability for hospital staff
- Quick Release CPR with access from both sides for backrest and knee break ensures the trolley is emergency ready at all times
- Collapsible vertical side rails fold away below the patient surface
- Fold-away carer push handles to maximise workspace
- 200mm castors minimise rolling resistance with central and directional locking functions
- Molded base cover facilitates cleaning and supports infection control processes
- Large corner bumpers minimise wear and tear on equipment and facilities
- Frame is finished with Biocote, a technologically advanced surface coating that promotes infection control via its anti-bacterial properties



OPTIONAL ACCESSORIES

ĺ		ľ
		j

V Pole (TSE674130)

「「

Oxygen Bottle Holder (TSE674110)



Folding Monitor Tray (TSE674120)



Mattress (TSE674315)



BIOCOTE Antibacterial frame coating for improved infection control









Code	TSE674310	
Overall Length	2270mm	
Overall Width	830mm	
Overall Height	580 - 860mm	
Patient Surface Length	2000mm	
Patient Surface Width	650mm	
Backrest Angle	0° to 90°	
Trendelenburg	-16° to 16°	
Castor Size	200mm	