

IDEA2 TREATMENT CHAIR




STANDARD INCLUSION



Battery Back-up

OPTIONAL ACCESSORIES



IV Pole CHT100720

COLOUR OPTIONS



Blue Vinyl



BIOCOTE
Antibacterial frame coating for improved infection control



FEATURES AND BENEFITS

- Premium therapy chair that effortlessly meets the needs of oncology departments and the demands of infusion and transfusion medicine
- Gas-assisted, height and laterally adjustable swing away armrests to suit individual user needs
- Ideal for procedures such as chemotherapy, dialysis, blood donation, and blood transfusion
- Extra-long leg rest for patient comfort during extended periods of treatment
- Intuitive handset allows for quick and easy chair repositioning
- Ergonomic push handle with inbuilt paper roll holder
- Electric adjustment allows for infinite positioning of the backrest, legrest, seat height, seat tilt and Trendelenburg (CPR shock position)
- Seamless patient surface for maximised infection control compliance
- 125mm twin-wheel castors with central locking that provides both total and directional locking

2 YEARS WARRANTY

SWL 200kg

Code	CHT100710
Overall Length	1600mm
Overall Width	840mm
Seat Height	570 - 890mm
Patient Surface Length	2130mm
Patient Surface Width	580mm
Backrest Angle	80°
Trendelenburg	12°
Castor Size	125mm