



ROOMING-IN CRIB



Exclusively Distributed by:





Rooming-in Crib

Rooming-in crib with vertically lowering padded safety sides which can be removed for cleaning

Full +/- trendelenburg functionality

Name plate for child identification

Ventilation/breathing holes located at both ends

Epoxy powder coating with antibacterial finish for improved infection control

Base frame is equipped with 100mm twin wheel braking castors, which are activated by the push handle

Gas assisted height adjustment

Padded sides ensure continuity between the crib and the mattress for a safe transfer of the baby to and from the crib

Can be placed either side of the bed



PED490200 Matris Paediatric Crib

PED490210 Clothes Basket

✓ Matris Cot is registered in the ARTG, AUST R ARTG 294878

✓ IEC60601-2-52 Compliant

✓ AS/NZS 2130:1998 Compliant

Code: PED490210

Clothes Basket





Colours



Blue



For a mother, being close to their new born is important to offer them reassurance, protection and love.

For a new born, it is important to be able to recognise their mother, to watch her, listen to her voice and feed from her.

The Matris Rooming-in Crib meets these needs, making it possible to be placed in complete safety close to the mother's bed.

Technical Data

Overall dimensions	89 cm x 52 cm x 91-99 cm
Mattress platform height	from 49 cm to 71 cm
Trendelenburg & Reverse Trendelenburg	12° & 6°
Safe Working Load & Max Patient Weight	20 kg & 12kg
Mattress in flexible expanded polyurethane. Flame retardant cover in unshrinkable polyester fibre (Trevira cs) HOMOLOGATION 11M D.M. 26/6/84.	





Malvestio2018

The Company reserves the right to introduce all the technical and/or aesthetic modifications deemed necessary to their products without prior notice.



malvestio.it



Exclusively Distributed by:



Building 3A/1 Moorebank Ave,
Moorebank NSW 2170
T – 1300 133 120
W – www.aidacare.com.au