



ELECTRIC PAEDIATRIC BED




Exclusively Distributed by:













Code: PED490120 Oxygen Bottle Holder













Electric Paediatric Bed

3 x section electrically operated paediatric ward bed

Epoxy powder coating with antibacterial finish for improved infection control

Base frame is equipped with 150mm twin wheel castors featuring a central locking system. The full width foot pedal is accessible from either side of the bed end

Adjustable mattress platform with electric length adjustment from 158-198cm. This function allows the bed to grow with the patient and is suitable from young children through to young adults

Blow moulded ABS folding side rails

Bed is equipped with battery backup as standard, allowing for power supply when transporting around wards or in the event of a power failure/outage

Emergency release levers on both sides of the bed

PED490100	Smile Paediatric Ward Bed
PED490110	Foam Mattress w/Extension Bolster
PED490120	Oxygen Bottle Holder

Overall dimensions	183-218cm x 89cm x 55-95cm
Trendelenburg	15°
Reverse Trendelenburg	15°






Trendelenburg and Reverse-trendelenburg



adjustable height











Malvestio2018

The Company reserves the right to introduce all the technical and/or aesthetic modifications deemed necessary to their products without prior notice.









Exclusively Distributed by:



Aidacare P/L

4 Noonan Rd, Ingleburn NSW 2565 T – 1300 133 120 W – www.aidacare.com.au