

TOPIC Sanitisers

Top of the range for hygiene and economy



At last more time for patients. TOPIC takes the strain





The new MEIKO TOPIC is the perfect and cost-effective MEIKO cleaning and disinfection appliance to meet the hygiene requirements of patient utensils in care facilities.

Depending on the situation in the building, MEIKO TOPIC appliances can be supplied for wall-mounting so that they are clear of the floor or for installation directly on the floor. As care units complete with everything that is necessary and beneficial, TOPIC stand-alone appliances are specially designed for the rapid and uncomplicated replacement of out-dated cleaning and disinfection technology. They are modified to fit into the position of the old appliance whatever the make of the original – and to use the same connections.

We will work with you to develop an ergonomic sluice room design which meets the most stringent demands of cost-effectiveness, optimisation of work processes, safety in the workplace and hygiene. In this way you will receive your own solution customised to meet the hygiene standards you require.

The innovative TOPIC range of appliances works with rigorous efficiency to create painstaking cleanliness and disinfection, releasing nursing teams from this work so they can devote themselves to their real job – interaction with their patients.

TOPIC is sophisticated technology. Robust technology made from stainless steel which is the starting point for a long working life, easy cleaning and meticulous hygiene. Top design: the deep-drawn cleaning chamber is designed to be self-cleaning with no corners and sharp edges.

The water jet system is arranged to provide intensive cleaning, guaranteeing optimal hygiene for care utensils.

Bed pans, urine bottles, commode buckets etc. are automatically and reliably emptied, cleaned and disinfected. And always safely behind a sealed door, of course

This means modern hygiene, convenience and a reduced work load for the entire nursing staff. The accountants can also breathe a sigh of relief – TOPIC's reliable, low maintenance technology makes only miserly demands on resources and budgets.

Your TOP TEN: the *TOPIC* hygiene advantages

- Ideal operation as a result of the "One Touch"operating system with large-sized membrane key-pads and unambiguous icons
- Hygiene security thanks to the A₀ value control system complying with DIN EN ISO 15883
- Very economical due to reduced use of water, detergent and rinse agents
- Short, individually adjustable and effective program cycles
- Maximum flexibility provided by multi-functional brackets for care utensils
- Patented, steam-proof cleaning and disinfection system
- All-round cleaning by the powerjet system
- Electronic controls with versatile cleaning programs
- Hygiene security created by thermal disinfection
- Multiple uses from functionally efficient options and extensive accessories





2 MEIKO TOPIC Sanitiser

Distributed by aidacare 3

TOPIC cleaning and disinfection appliances provide maximum ease of operation. Starting the program, appliance control and operational analysis are all managed via the simple TOPIC membrane key-pad.

Simple "One Touch" technology and the innovative M-Commander analysis software help you to help your patients – only with TOPIC is patient care so simple and easy.





Everything plainly in view:

TOPIC with M-Commander and M-Commander W, the analysis software for MEIKO cleaning and disinfection appliances.



MEIKO Topic 10 & Topic 20





Topic 10

- For the clear-of-the-floor wall mounting
- The easy way of installing a cleaning and disinfection appliance.
- Space-saving and ready for connection.
- Available with automatic doors

Topic 20

- Floor mounted cleaning and disinfection appliance with plinth
- Just place in position and connect.
- The ideal way of utilising existing connections
- Available with automatic door

Code	Description	Width	Depth	Height
UPM200500	Sanitiser - Meiko - Topic 10	500 mm	450 - 600 mm	1630 mm
UPM200510	Sanitiser - Meiko - Topic 10 AT - with Auto Door	500 mm	450 - 600 mm	1630 mm
UPM200520	Sanitiser - Meiko - Topic 20	500 mm	450 - 600 mm	1730 mm
UPM200530	Sanitiser - Meiko - Topic 20 AT - with Auto Door & Foot	500 mm	450 - 600 mm	1730 mm

4 MEIKO TOPIC Sanitiser

Distributed by aidacare 5

Accessories and range of fixtures



MEIKO Topic 20 XL



Topic 20 XL

- Double floor mounted cleaning and disinfection appliances with plinth
- Just place in position and connect.
- The ideal way of utilising existing connections
- Available with automatic door

Code	Description	Width	Depth	Height
UPM200540	Sanitiser - Meiko - Topic 20 XL (side-by-side)	1000 mm	450 - 600 mm	1730 mm
UPM200550	Sanitizer - Meiko - Topic 20 XL AT (side-by-side) - with Auto Door & Foot	1000 mm	450 - 600 mm	1730 mm

In addition to the Meiko cleaning and disinfection technology, we can also supply you with a comprehensive range of accessories for safe and ergonomic working and the hygienic storage of care utensils. We illustrate below 6 examples of practical solutions. We would be happy to develop a range of suggested solutions specially for you. Just give us a call. We would be pleased to help you at a moment's notice.































ACCESSORIES

	Code	Description	Width	Depth	Height
0	UPA400100	Basket insert with hinged lid - Meiko	-	-	-
2	UPA400110	Basket with hinged lid for small items - Meiko	-	-	-
3	UPA400200	Bottle Rack - Meiko TU50 - 8 Urine Bottles (drip tray sep)	500 mm	350 mm	650 mm
4	UPA400210	Bottle Rack - Meiko TU70 - 11 Urine Bottles (drip tray sep)	700 mm	350 mm	650 mm
5	UPA400250	Pan & Bottle Rack - Meiko TUS 140 - 9 Pan & 9 Bottle (drip tray sep)	1400 mm	350 mm	650 mm
6	UPA400260	Pan & Bottle Rack - Meiko TUS 50 - 6 Pan & 6 Bottle (drip tray sep)	500 mm	350 mm	650 mm
7	UPA400300	Pan Rack - Meiko TS100 - 12 Pans (drip tray sep)	1000 mm	350 mm	650 mm
8	UPA400310	Pan Rack - Meiko TS50 - 3 Pans (drip tray sep)	500 mm	350 mm	650 mm
9	UPA400320	Pan Rack - Meiko TS70 - 6 Pans (drip tray sep)	700 mm	350 mm	650 mm
10	UPA400330	Pan Rack - Meiko TS90 - 9 Pans (drip tray sep)	900 mm	350 mm	650 mm
•	UPA400340	Pan Rack (Not Mounted) - Meiko - For vertical storage of bedpans - 4 Pan	-	-	-
12	UPA400350	Pan Rack (Not Mounted) - Meiko - For vertical storage of bedpans - 8 Pan	-	-	-
13	UPA400400	Rack (Internal) - Meiko - SUH 1 Internal rack (Non-standard in Aust)	-	-	-
14	UPA400410	Rack (Internal) - Meiko - SUH 2 Quick-change rack (Standard in Aust)	-	-	-
15	UPA400420	Rack (Internal) - Meiko - TopCut B - Quick Change Rack with Bag Empty	-	-	-
16	UPA400430	Rack (Internal) - Meiko - TopCut H - Bag Emptying Rack (internal)	-	-	-
17	UPA400500	Drip tray (For all wall-mounted racks) - Meiko - Stainless Steel (order with rack)	-	-	-





DISTRIBUTED BY



1300 133 120 www.aidacare.com.au