


















## MANUAL CHAIRS

|                       | ADULT   |   |   |   |   | ADULT / PAED  | PAEDIATRIC  |
|-----------------------|---|---|---|---|---|---|---|
|                       | Folding   |   |   |   | Rigid   | Rigid   | Folding   |
|                       | Start Range (M1, M4, M5)  | Motus CV, CS  | Avantgarde (CV, CS, CLT)  | Avantgarde XXL  | Voyager Evo   | Ventus  | Start JNR   |
|                       |  |  |  |  |  |  |  |
| Seat Width (mm)       | 380 - 580   | 355 - 555   | 320 - 500   | 460 - 620   | 280 - 460   | 280 - 500   | 280 - 380   |
| Seat Depth (mm)       | 400 - 490   | 400 - 460   | 360 - 540   | 400 - 560   | 360 - 520   | 300 - 500   | 240 - 400   |
| Floor to Seat Ht (mm) | 380 - 510   | 420 - 520   | 390 - 550   | 460 - 530   | 440 - 530   | 400 - 540   | 395 - 520   |
| Max User Wt (kg)      | 125 - 160   | 140   | 140,140,100   | 180   | 100 - 140   | 140   | 90  |
| Weight (kg)           | 14.5 - 17.5   | 13  | 12,11,10  | 17  | 7.73  | 12  | 13.5  |

## POWER CHAIRS

|                       | ADULT  |  |   |   |  |   |  |   |   | PAEDIATRIC  |
|-----------------------|--|--|---|---|--|---|--|---|---|---|
|                       | Rear Wheel Drive   |  | Rear / Front Wheel Drive  |   |  | Front Wheel Drive   | Rear / Front Wheel Drive   | Rear Wheel Drive - Sit to Stand   |   | Rear Wheel Drive  |
|                       | A200   | B400   | B500  | B500S   | B600   | C1000DS   | C2000  | Xeno  | Paragolfer  | Skippi  |
|                       |       |             |  |  |  |                  |       |  |  |  |
| Seat Width (mm)       | 340 - 480  | 340 - 480  | 340 - 560   | 340 - 560   | 340 - 560  | 380 - 560   | 380 - 560  | 430 - 480   | 460   | 300 - 380   |
| Seat Depth (mm)       | 340 - 500  | 340 - 500  | 340 - 530   | 380 - 460   | 360 - 560  | 360 - 560   | 380 - 560  | 370 - 530   | 356 - 510   | 280 - 380   |
| Floor to Seat Ht (mm) | 430 - 480  | 400 - 500  | 450 - 550   | 450   | 450 - 550  | 450 - 550   | 470 - 530  | 530 - 630   | 510 - 660   | 430   |
| Max User Wt (kg)      | 100  | 140  | 140   | 140   | 140 - 180  | 140 - 200   | 160 - 200  | 136   | 140   | 50  |
| Base W x L (mm)       | 570 x 1000   | 580 x 1080   | 640 x 1080  | 640 x 1080  | 640 x 1080   | 650 x 1200  | 680 x 1390   | 650 x 1010  | 900 x 1490  | 570 x 850   |
| Max Speed (kph)       | 6  | 10   | 10  | 6   | 10   | 10  | 10   | 10  | 10  | 6   |
| Features              | LED single function control. Manual back recline. 20km range. Detachable control panel | LED single function control. Manual elevating footrest. Curb climb assist option. 30km range | Expandable controls. 5 power seat options. 35km range                             | Expandable controls. 5 power seat options. 35km range                             | Expandable controls. Telescoping frame. Lights and rear suspension. 35km range     | Power drive castors. Expandable controls. Telescoping frame. Lights and rear suspension. 35km range | Rugged terrain. Power drive castors. Seat raise 40cm. Centre mount leg rests. 60km range | Power sit to stand 35km range. Tilt & recline seating optional                      | Powerful sit-to-stand used in sports. Wide base for stability                       | Expandable controls. Detachable seat for travel. Attendant control & lighting       |