

PAEDIATRIC CATALOGUE



About Aidacare

Aidacare is one of Australia’s leading healthcare equipment and service providers. Operating since 1987, our longevity is the result of our commitment to excellence and continuous improvement.


We are a trusted product partner and service provider across the Acute Care, Residential Aged Care, Home and Community Care, and Complex Rehabilitation sectors. Our team deliver on our core values of providing unparalleled customer service and equipment expertise. Aidacare staff are highly trained and skilled professionals who respond sensitively and expertly to our clients’ needs.

We are dedicated to supplying the highest quality products. We partner with industry leading brands and design and manufacture a suite of our own equipment solutions that serve as a benchmark for innovation, quality and reliability. Our national service teams provide maintenance and asset management across all product categories and to all our customer segments. We also provide a range of high quality training services to customers and allied health partners

Aidacare offers a national retail, wholesale and distribution network providing truly local support from our many locations and partnerships.

NDIS

Aidacare is proud to be an official contracted provider under the NDIS program. We supply equipment and support, working closely with families, coordinators, therapists, and Veterans to achieve exceptional outcomes. Codes are listed in product specifications.



Visit Aidacare Funding & Contracts



Complex Rehab


Our specialised seating and customisation service is provided by an expert team of Physiotherapists, Occupational Therapists, and skilled technicians. The customised seating can be built in-house to individual specifications.




Aidacare Complex Rehab Website

Online Learning Centre

Visit aidacare.com.au for a user-friendly guide to our extensive product range, clinical assistance guides, and educational tools. Explore videos, PDFs, manuals, and brochures to make informed decisions with ease. Our platform empowers healthcare professionals and users alike with concise insights and detailed specifications for confident decision-making.




Visit our Aidacare Learning Centre



Retail Showrooms

Visit any of our 90 showrooms, where our friendly staff are ready to help you identify the perfect solution for your needs.

Our showrooms also serve as convenient training centres and venues for equipment trials with your clients.



Aidacare Stores & Network Partners

Product Catalogue

Explore our full range of high-quality assistive technology healthcare equipment via our product catalogue.



Aidacare Product Catalogue



Service & Maintenance

Implement effective Risk Management by scheduling regular servicing and maintenance for your equipment, reducing breakdown risks and aiding compliance with accreditation requirements. Benefit from asset tracking and equipment audits as part of our comprehensive services.




Aidacare's Preventative Maintenance & Asset Management Programs

Complex Rehab Equipment Catalogue

Explore our full range of highly specialised complex rehab products via our Complex Rehab Equipment Catalogue.




Aidacare Complex Rehab Equipment Catalogue





Scan to browse our range of Paediatric products online or visit aidacare.com.au to see our full range of products and services



Contents

Standing Aids	6
Seating Aids	10
Strollers	19
Walking Aids	23
Manual Wheelchairs	27
Power-Assist Systems	33
Power Wheelchairs	35
Head Supports	40
Night Time Positioning	41
Bathroom	44

Standing Aids



Leckey Pronestander



- The Pronestander is a stander with a choice of 5 prone positions ranging from 45° to 85°
- Assists in standing and weight bearing
- Highly adjustable to suit individual users
- Includes chest and knee supports, hip belt, angle adjustable tray, activity bowl and a mobility base with 4 locking castors

LECKEY°

3 YEARS WARRANTY

Code	Description	User Standing Height	Chest Pad Height	SWL
PED479250	Size 1, 1.5 - 5 years (as per script)	800 - 1150	470 - 680	60kg
PED479260	Size 2, 4 -10 years (as per script)	1100 - 1450	680 - 975	80kg
PED479270	Size 3, 9 - 18 years (as per script)	1400 - 1850m	975 - 1350	90kg

NDIS 05_044839811_0103_1_2

Leckey Freestander



- The Freestander is a robust upright stander ideally suited to the classroom, therapy room or home
- Designed to assist in standing and weight bearing development
- Highly adjustable to suit individual users
- Includes chest and knee supports, hip belt, angle adjustable tray, activity bowl and a mobility base with 4 x locking castors for stability

LECKEY°

3 YEARS WARRANTY

Code	Description	User Standing Height	Chest Pad Height	SWL
PED479170	Size 1,1.5 - 5 years (as per script)	800 - 1151	470 - 680	60kg
PED479180	Size 2, 4 - 10 years (as per script)	1100 - 1448	680 - 975	80kg
PED479190	Size 3, 9 - 18 years (as per script)	1400 - 1849	975 - 1350	90kg

NDIS 05_044839811_0103_1_2

Medifab Jenx Monkey Standing Frame

- Jenx Monkey is a gentle and fun introduction to standing, providing prone to upright positioning
- De-rotational thoracic and pelvis straps/pads fasten in the centre and adjust independently from either side
- Polyurethane support surfaces with antimicrobial agent
- Integral angle-adjustable tray with both bowl and smooth tray infill for a wide variety of different activities
- Quick, simple and easy to use
- Comes as a complete unit and very few adjustments are required

SWL 25KG 2 YEARS WARRANTY

Code	Description	Angle Range	Base Size	Chest Pad Height
PED479375	Jenx Monkey Standing Frame (as per script)	20° - 90°	880 x 575	510 - 750

NDIS 05_120603831_0103_1_2



medifab

Medifab Jenx Multistander Abduction

- Versatile standing frame, suitable for children from around 9 months right up to 13 years
- PU boards and support pads are flexible, easy clean and hygienic
- De-rotational thoracic and pelvis straps/pads fasten in the centre and adjust independently from either side
- Independently adjustable thoracic and pelvic boards cater for the child's individual proportions
- Deliberately open design ensures you can always see precisely how the user is positioned each and every time

2 YEARS WARRANTY

Code	Description	User Standing Height	Chest Pad Height	SWL
PED479630	Abduction Stander (as per script)	1160 - 1750	490 - 1170	35 - 60kg

NDIS 05_120603831_0103_1_2



medifab

Medifab Jenx Multistander

- Versatile standing frame suitable for children from around 9 months up to 13 years
- Flexible support pads which are easy to clean
- Thoracic and pelvic boards can be adjusted independently for each child's individual proportions

2 YEARS WARRANTY

Code	Description	Angle Range	Max. User Height	Chest Pad Width	SWL
PED491530	Size 1, 9 months - 6 years (as per script)	10° - 90°	1160	90	35kg
PED491540	Size 2, 3 - 13 years (as per script)	0° - 80°	1750	100	60kg



medifab

Leckey Squiggles Stander

LECKEY[°]



- Squiggles Standing System supports the child in a comfortable, secure, prone, supine or upright standing posture
- Chassis can clip on and off in seconds for folding and storage
- Adjustable to accommodate children with differing needs
- Accessories available include tray, knee blocks, supports, headrests which all attach without tools

3 YEARS WARRANTY

SWL 22KG

Code	Description	User Standing Height	Chest Pad Height
PED479379	1 - 5 years (as per script)	750 - 1100	550 - 820

NDIS 05_044839811_0103_1_2

Leckey Squiggles + Stander

LECKEY[°]



- Hip abduction of up to 60° in prone, upright and supine standing.
- Fully horizontal loading in supine
- New wrap around knee supports to better accommodate all sizes, with added comfort and support
- New, improved, pivot chassis
- Height-adjustable, foldable push handle for ergonomic use
- Individually adjustable footplates to accommodate leg length discrepancy

3 YEARS WARRANTY

SWL 22KG

Code	Description	User Standing Height	Chest Pad Height
PED479380	1 - 5 years (as per script)	750 - 1110	550 - 820

NDIS 05_120603831_0103_1_2

Leckey Mygo Stander

LECKEY[°]



- The Mygo Stander provides a 3-in-1 prone, supine and upright standing positions
- Comes with height adjustable chest pad that is ergonomically designed to imitate the shape of the rib cage
- Includes hip support, chest support, 4 castors with brakes and height adjustable footplates and knee supports

3 YEARS WARRANTY

Code	Description	Width	Length	Chest Width	Hip Width	Height of User Range	SWL
PED479335	Size 1, 4 - 10 years (as per script)	630	978	160 - 280	220 - 330	1080 - 1450	50kg
PED479336*	Size 2, 8 - 14 years (as per script)	650	1130	200 - 320	280 - 390	1350 - 1750	60kg

NDIS 05_044839811_0103_1_2 *05_120603831_0103_1_2

Leckey Horizon Stander

LECKEY[°]



- Durable 3-in-1 stander which can be positioned in prone, supine or upright
- Fully adjustable to suit the individual user and is ideal in multiuser environments
- Available with headrest, adjustable sandals, lateral supports and pelvic positioning belt

3 YEARS WARRANTY

Code	Description	User Standing Height	Chest Pad Height	SWL
PED479105	Size 1, 4 - 10 years (as per script)	1000 - 1400	820 - 1150	50kg
PED479130	Size 2, 9 - 18 years (as per script)	1400 - 1800	1000 - 1450	80kg
PED479135	Size 3, 14 years - adult (as per script)	1550 - 1800	1000 - 1450	100kg

NDIS 05_044839811_0103_1_2

Apex Mobility Rifton Stander

APEX MOBILITY



- Versatile configuration as supine, prone, or multi-position
- Allows hip abduction up to 30° to promote hip joint development
- Independent adjustments for hip and knee flexion and leg length
- Broad range of positioning components for comfort
- Frame can be raised and lowered for eye-level interaction
- Elegant and approachable design centred on the child

1 YEAR WARRANTY

Code	Description	Knee to Hip Height	Sandal to Knee Height	User Standing Height	SWL
PED491780	Size 1 (as per script)	159 - 267	172 - 279	711 - 1066	34kg
PED491790	Size 2 (as per script)	216 - 343	254 - 387	965 - 1346	100kg

NDIS 05_044839811_0103_1_2



Seating Aids



Leckey Early Activity System



- Strong, durable, practical and adjustable seating system, providing moderate support for kids and young adults aged 1-18 years
- System supplied as a “kit” to support a young child to sit, lie on their tummy, back or side
- Can be used by therapists as well as parents
- Designed to suit children up to 3 years old

LECKEY°

3 YEARS WARRANTY

SWL 30KG

Code	Description	Seat Depth	Hip Width	Backrest Height
PED480100	Early Activity System	190 - 260	150 - 200	100

NDIS 05_120603831_0103_1_2

Leckey PAL Classroom Seat



- Attractive, robust and easy to use classroom and nursery seating system for kids with mild to moderate postural needs
- Matching tables available
- Improves stability, reduces fatigue and longer periods of concentration and fine motor activity

LECKEY°

3 YEARS WARRANTY

Code	Description	Seat Depth	Hip Width	Backrest Height	SWL
PED478800*	Size 1, 1 - 3 years (as per script)	180 - 240	150 - 250	265	30kg
PED478806	Size 4, 7 - 12 years (as per script)	300 - 400	280 - 380	450	50kg

NDIS 05_044839811_0103_1_2 *05_122503111_0105_1_2

Leckey Everyday Seat

- Designed to help physical, cognitive and sensory development for babies with special needs
- 4-point harness gives a secure, stable base of support
- Height and angle adjustable flip-up footplate to facilitate transfer
- Highly configurable with 3 sizes, a range of chassis options and 5 colour options

3 YEARS WARRANTY

Code	Description	Seat Depth	Hip Width	Backrest Height	Handle Height	SWL
PED490010	Size 1 (as per script)	150 - 280	160 - 280	300 - 400	200 - 300	30kg
PED490015	Size 2 (as per script)	230 - 380	260 - 380	400 - 520	200 - 300	50kg

NDIS 05_122500821_0105_1_2



LECKEY°

Leckey Squiggles Seating System

- Squiggles Seating System is designed for younger children requiring postural support
- Adjustable seating system with tilt-in-space mechanism
- Designed to suit children from 1 to 5 years and allows for growth
- Seating accessories include thoracic supports, sandals, tray, harness plus others
- Can be ordered with a choice of base including wheelchair base, stroller base or high low indoor base
- Seating system fits a number of transporting systems

SWL 22KG

3 YEARS WARRANTY

Code	Description	Seat Depth	Hip Width	Backrest Height
PED479369	1 - 5 years (as per script)	160 - 260	190 - 290	320 - 380

NDIS 05_120603831_0103_1_2



LECKEY°

Leckey Mygo Seating System

- Multifunctional seating system designed for children and adolescents requiring moderate to complex postural support
- Includes leg guides, four-point pelvic positioning harness and adjustable sacral support
- Seat width and depth are continuously adjustable allowing for a large degree of growth
- Available with unique Leckey pelvic cradle for improved postural control
- Range of additional postural supports available
- Seating system fits a number of transporting systems

3 YEARS WARRANTY

Code	Description	Seat Width	Seat Depth	Backrest Height	SWL
PED479339	Size 1, 3 - 10 years (as per script)	200 - 325	270 - 420	360 - 470	50kg
PED479338	Size 2, 8 - 14 years (as per script)	220 - 345	350 - 470	460 - 570	60kg

NDIS 05_122503111_0105_1_2



LECKEY°

Leckey BeMe

LECKEY^o



- Hi-low mobility base enhances participation and inclusion with easy manoeuvrability
- Modular seating that is customisable to fit the child's changing needs as they grow
- Contoured and adjustable for ongoing support and comfort

3 YEARS WARRANTY

Code	Description	Seat Width	Seat Depth	Backrest Height	SWL
PED479640	Small (as per script)	275	170 - 325	360 - 470	35
PED479642	Medium (as per script)	325	275 - 390	390 - 520	50
PED479644	Large (as per script)	415	350 - 475	500 - 650	70

NDIS 05_120603831_0103_1_2

Spex Wonderseat Space Seating System

spex^o



- Solves common challenges with a small wheelchair system
- Works as a cohesive unit
- Suitable for children needing moderate to high postural support
- Ideal for children using a tilt-in-space wheelchair
- Manages the influence of gravity on posture
- Provides postural support for children with little to no sitting control
- Supports can be easily moved away to aid with transfers

2 YEARS WARRANTY

SWL 47KG

Code	Description	Seat Width	Seat Depth	Back Support Width	Back Support Height	Product Weight
PED491665	Size 0 or 1, 6 months - 10 years	150 - 300	150 - 380	200 - 250	250 - 430	7.9 - 9.5kg

NDIS 05_180909811_0103_1_2

Spex Wonderseat Bounce Seating System

spex^o



- Solves common challenges with a small wheelchair system
- Designed to work as one cohesive unit
- Suitable for children needing moderate to low postural support
- Ideal for children transferring in and out of wheelchairs with support
- Provides postural support for extended periods in the chair
- Features fewer moving parts and swing-away options for durability
- Offers lower support compared to the Space model
- Not suitable for active users or self-propelling children

2 YEARS WARRANTY

SWL 47KG

Code	Description	Seat Width	Seat Depth	Back Support Width	Back Support Height	Product Weight
PTC626792	Size 0, 6 months - 5 years	150 - 250	150 - 280	200	250 - 380	7.6kg

NDIS 05_043303511_0103_1_2

Medifab Jenx Junior Seating System

medifab^o



- Flexible supportive seating for children aged 18 months to 12 years
- Wide variety of head, trunk, leg, and foot support options
- Contoured seat unit made of super-soft polyurethane and memory foam
- Designed for true pelvic stability and unparalleled comfort
- Developed with input from children, parents, therapists, and care staff to meet diverse needs

SWL 50KG

2 YEARS WARRANTY

Code	Description	Back Angle	Seat Width	Seat Depth	Backrest Height
PED479346	Size 1, 18 months - 5 years (as per script)	15° recline	150 - 375	210 - 305	345 - 510
PED479348	Size 2, 3 - 8 years (as per script)	15° recline	150 - 375	240 - 355	345 - 510
PED479353	Size 2, 5 - 12 years (as per script)	15° recline	150 - 375	260 - 440	410 - 575

NDIS 05_120603831_0103_1_2

Rifton Activity Chair

APEX MOBILITY



- Versatile and adaptable for pediatric to adult clients
- Provides ultimate clinical positioning, revolutionizing active indoor seating
- Available in standard and hi-lo base options for maximum convenience
- Features tool-free adjustments for home or classroom activities
- Accommodates growth with three sizes suitable
- Offers a wide range of accessories to meet various seating needs

1 YEAR WARRANTY

Code	Description	Frame Width	Frame Depth	Seat Height	Seat Angle Forward / Backward	Seat Width With Hip Guides	SWL
PED491750	Small	533 - 597	546 - 724	241 - 546	15° / 15°	178 - 229	34kg
PED491760	Medium	591 - 648	572 - 730	318 - 584	15° / 15°	216 - 292	68kg
PED491770	Large	667 - 705	775 - 895	419 - 610	15° / 15°	279 - 356	113kg

NDIS 03_090000911_0103_1_1

Medifab Jenx Corner Seat

medifab^o



- Adjustable height back, with optional upper back support and abduction block for proper positioning
- Corner shape for midline maintenance; optional triangular cushions for flat back with side wings
- Ideal first seat for children, compatible with Jenx Nursery Table
- Secured to the table with touch and close tape straps
- Promotes long sitting with legs stretched out to maintain muscle length behind the knee

2 YEARS WARRANTY

Code	Description	Base Size Foot Print	Seat Depth	Back Pad Height	SWL
PTA535105	6 months - 10 years (as per script)	490 - 630	350 - 465	315 - 680	25 - 55kg

NDIS 05_122506111_0105_1_2

Seating Matters Phoenix Paediatric Chair



- Gives maximum support to the upper body through lateral and head supports, also adapt to the changing needs of a person being rehabilitated or meet the long term needs of a progressive condition
- Adjustable lateral, shoulder and head supports can conform to posture & position easily
- Tilt-in-space and adjustable components
- Meets range of individual support needs, sizes, and levels of mobility

7 YEARS WARRANTY

Frame

1 YEAR WARRANTY

Upholstery

SWL

160KG

Code	Description	Hip Width	Seat Depth	Backrest Height	Floor-to-seat Height
CHP197800	Phoenix Paediatric Chair (as per script)	290 - 380	250 - 400	680 - 830	540
NDIS 03_090000911_0103_1_1					

Medifab Gravity Chair




- The chair is supportive and comfortable, protecting posture while providing relaxation indoors, provides a snug and 'safe' feeling, catering to users who need increased sensory feedback
- It follows the principle of 'Tension and Relaxation,' offering a stable alternative to floor or sofa seating
- Optional Accessories available to accommodate user needs



Optional Side Filler Pads Optional Depth Adj. Liner Optional Abductor Pommel



Optional Mobile Base



Optional Therapy Tray



Optional Head Support



Optional Foot Support

1 YEAR WARRANTY

Code	Description	Seat Depth	Backrest Height	Height	Width	Depth	User Height
PED470082	Small (as per script)	240	600	550	460	600	900 - 1100
PED470086	Medium (as per script)	300	690	610	490	710	1100 - 1300
PED470088	Large (as per script)	380	740	670	570	820	1300 - 1500
PED470090	Extra Large (as per script)	430	810	710	600	910	1300 - 1700
PED470092	Extra Large Wide (as per script)	430	810	720	710	910	1300 - 1700

Seating Matters Sorrento Paediatric Chair



- Designed to improve and maximise comfort as well as postural support and pressure management, the Sorrento is simple to operate
- Adaptability and effortless adjustability
- Tilt-in-space, forward tilt and adjustable back angle recline can be used to improve posture and redistribute pressure
- Optional lateral supports and pillows

SWL

160KG

7 YEARS WARRANTY

Frame

1 YEAR WARRANTY

Upholstery

Code	Description	Hip Width	Seat Depth	Backrest Height	Floor-to-seat Height
CHP197819	Sorrento Paediatric Chair	250 - 400	290 - 380	650	540
NDIS 05_180921121_0103_1_2					Available in multiple manual & electric combinations

Rolapal Snooza Chair



- All-foam Snooza Chair uses gravity to keep the occupant safe with good postural stability while they are out of their wheelchair
- Soft snug close support helps the occupant to relax & high armrests promote better breathing
- Carry handles make transportation easy and large storage pockets located on each side
- Wipe-clean medical grade vinyl cover
- The centre comfort wrap supplied is removable and washable
- Standard Snooza Chair kit includes:
 - 1 x Wedge for angle adjustment
 - 1 x Adjustable headrest
 - 1 x Padded hip belt

18 MONTHS WARRANTY

SWL 50KG

Code	Description	Floor to Seat	Seat Depth	Back Height	Width at Knees	Width at Hip	Width	Length
PED472000	Size 1, Apple	210	260	610	260, 320*	200, 250*	460	610
PED472010	Size 1, Royal Blue	210	260	610	260, 320*	200, 250*	460	610
PED472020	Size 1, Orange	210	260	610	260, 320*	200, 250*	460	610
PED472030	Size 1, Hot Pink	210	260	610	260, 320*	200, 250*	460	610
PED472040	Size 2, Apple	240	320	700	300, 350*	230, 280*	480	720
PED472050	Size 2, Royal Blue	240	320	700	300, 350*	230, 280*	480	720
PED472060	Size 2, Orange	240	320	700	300, 350*	230, 280*	480	720
PED472070	Size 2, Hot Pink	240	320	700	300, 350*	230, 280*	480	720
PED472080	Size 3, Apple	270	400	750	370, 420*	270, 320*	570	870
PED472090	Size 3, Royal Blue	270	400	750	370, 420*	270, 320*	570	870
PED472100	Size 3, Orange	270	400	750	370, 420*	270, 320*	570	870
PED472110	Size 3, Hot Pink	270	400	750	370, 420*	270, 320*	570	870
PED472120	Size 4, Apple	300	450	810	380, 430*	300, 340*	600	910
PED472130	Size 4, Royal Blue	300	450	810	380, 430*	300, 340*	600	910
PED472140	Size 4, Orange	300	450	810	380, 430*	300, 340*	600	910
PED472150	Size 4, Hot Pink	300	450	810	380, 430*	300, 340*	600	910

*with infill pads

Firefly Goto Seat



Optional Firefly Backpack



Optional Pressure Relief Cushion

PINK/PURPLE

LIGHT/DARK BLUE

2 YEARS WARRANTY

- Makes daily tasks easier for kids and ensures kids don't miss out on activities
- Used worldwide in shopping carts, swings, dining tables, classrooms, trikes, floor sitting, and more
- Suitable for ages 1 to 8 years
- Height-adjustable headrest
- Wrap-around laterals
- Optional floor sitter with 26 recline settings
- Available in wipe-clean vinyl version



Code	Description	Seat Base to Top of Headrest	Lateral Support Width	Seat Depth	SWL
PED490540	Size 1 (1 - 3 years) Standard Headrest - Pink / Purple	380 min. 540 max.	190 min. 260 max.	175	15kg
PED490545	Size 1 (1 - 3 years) Standard Headrest - Light/Dark Blue	380 min. 540 max.	190 min. 260 max.	175	15kg
PED490555	Size 1 (1 - 3 years) Advanced Headrest - Pink / Purple	380 min. 540 max.	190 min. 260 max.	175	15kg
PED490560	Size 1 (1 - 3 years) Advanced Headrest - Light/Dark Blue	380 min. 540 max.	190 min. 260 max.	175	15kg
PED490590	Size 2 (3 - 8 years) Standard Headrest - Pink Purple	510 min. 680 max.	210 min. 300 max.	244	30kg
PED490595	Size 2 (3 - 8 years) Standard Headrest - Light/Dark Blue	510 min. 680 max.	210 min. 300 max.	244	30kg
PED490600	Size 2 (3 - 8 years) Advanced Headrest - Light/Dark Blue	510 min. 680 max.	210 min. 300 max.	244	30kg
PED490610	Size 2 (3 - 8 years) Advanced Headrest - Pink/Purple	510 min. 680 max.	210 min. 300 max.	244	30kg
PED490650	Size 1 (1 - 3 years) Pressure Relief Cushion - Grey				
PED490655	Size 2 (3 - 8 years) Pressure Relief Cushion - Grey				
PED491250	Firefly Backpack for use with GoTo & Splashy				

NDIS 05_122500821_0105_1_2

Firefly Goto Vinyl Seat



Optional Firefly Backpack



Optional Pressure Relief Cushion

RED VINYL

2 YEARS WARRANTY

- Lightweight and portable
- Ideal for kids with special needs
- Promotes inclusion and participation in everyday life
- New wipe-clean vinyl cover for effortless maintenance
- Suitable for nursery, school, therapy, and busy homes



Code	Description	Seat Base to Top of Headrest	Lateral Support Width	Seat Depth	SWL
PED490535	Size 1 (1 - 3 years) Standard Headrest - Red Vinyl	380 min. 540 max.	190 min. 260 max.	175	15kg
PED490550	Size 1 (1 - 3 years) Advanced Headrest - Red Vinyl	380 min. 540 max.	190 min. 260 max.	175	15kg
PED490585	Size 2 (3 - 8 years) Standard Headrest - Red Vinyl	510 min. 680 max.	210 min. 300 max.	244	30kg
PED490605	Size 2 (3 - 8 years) Advanced Headrest - Red Vinyl	510 min. 680 max.	210 min. 300 max.	244	30kg
PED490660	Size 1 (1 - 3 years) Pressure Relief Cushion - Red Vinyl				
PED490665	Size 2 (3 - 8 years) Pressure Relief Cushion - Red Vinyl				
PED491250	Firefly Backpack for use with GoTo & Splashy				

NDIS 05_122500821_0105_1_2

Leckey Carrot 3000 Car Seat

- Carrot 3000 caters for safety and the importance of keeping the child in a supported posture for travel
- Various extension accessories and multi-adjustable features allow longer use during the child's life and suit the wide spectrum of different needs
- Ergonomically designed head support holds the head and protects from side impacts and centrifugal forces
- 2 positions of seat depth adjustment for a better fit
- Shoulder protector wings: Extendable wings hold and protect the body from side impacts and centrifugal forces
- Large chest support holds the body securely
- Recline function with angle adjustable seat and backrest

1 YEAR WARRANTY

Code	Description	Age Range	Seat Depth	Seat Width	Shoulder Height	SWL
PED478199	Carrot 3000 (as per script)	3 - 15	250 - 400	310	320 - 520	15 - 36kg

NDIS 05_121212811_0103_1_2



medifab

Medifab Ipai-LGT Car Seat

- The Ipai LGT is lightweight, easy to use, and provides mild body support for children with additional needs.
- It comes standard with ISOFIX, thoracic supports, and a padded 5-point harness
- A large, stable tee pad in the harness ensures firm support across the child's pelvis

SWL 36KG

1 YEAR WARRANTY

Code	Description	Seat Depth	Seat Width	Head Support Height	Backrest Height
PED478478	Medifab Ipai-LGT Car Seat	320	Front: 430 Rear: 310	590 - 790	380 - 580
PED478473	Padded Tray Table				

NDIS 05_121212811_0103_1_2



medifab

Medifab Starlight-NXT Car Seat

- Fully adjustable car restraint with excellent posture support and an inbuilt tilt function
- Suitable for young children from 9kg and grows with the child up to 36kg
- Offers comprehensive support and adaptability as the child grows

1 YEAR WARRANTY

Code	Description	Seat Depth	Hip Width	Backrest Height	SWL
PED478477	Medifab Starlight-NXT Car Seat	260	300	270 - 390	9 - 36kg
PED478471	Padded Tray Table				

NDIS 05_121212811_0103_1_2



medifab

Medifab Hero NXT Car Seat

medifab



- Modified to meet the postural needs of children with additional requirements
- Fits most vehicles with ISOFIX connections
- Secure attachment without car seat belts, ensuring easy handling

1 YEAR WARRANTY

SWL 18KG

Code	Description	Seat Depth	Hip Width	Backrest Height	SWL
PED478479	Rear facing (as per script)	250	210	210	18kg
PED470080	Front facing (as per script)	250	210	210	18kg
PED478472	Padded Tray Table				

NDIS 05_121212811_0103_1_2

Firefly Scooot 4 in 1



- Scooter Backrest - Our improved higher backrest provides kids with more support and reduces the amount of trunk control needed to use Scooot
- Push pole aids the child to propel themselves out of the chair
- Footplates and Wheels - Scooot's durable and easily gripped rubber tyres make is suitable for outdoor use
- The footplate can be easily adjusted to your child's exact leg length
- Advanced Backrest Cushion - The laterals move up and down to give support where it's needed
- Backrest Cushion Lateral Length: 790mm min. -1010mm max.
- Push Pole Height: 800mm min. - 890mm max.



2 YEARS WARRANTY

SWL 30KG

Code	Description	Backrest Width	Backrest Height	Seat Base Width	Wheel Arch Width	Hip to End of Footplate when Seated	Backrest Cushion
PED490890	Scoot 4 in 1, 2 - 6 years	260	300	230	230	Not extended: 520 Fully extended: 630	Lateral length: 420
PED490900*	Scoot 4 in 1 with Advanced Headrest, 2 - 6 years	260	300	230	230	Not extended: 520 Fully extended: 630	Lateral length: 420

PED490910 Advanced Headrest only, 2 - 6 years

NDIS 05_122500821_0105_1_2 *05_120603831_0103_1_2

Firefly Playpak Portable Activity Kit



- The various supports allow you to achieve an even wider variety of developmental postures and positions such as; tummy lying, floor sitting, back lying and 4-point kneeling
- Therapeutic and early intervention strategies allow your babies and children to develop healthy habits
- Playpak is light, handy and can go anywhere



2 YEARS WARRANTY

SWL 30KG

Code	Description	Length	Width	Depth
PED490800	Folded Bag	440	490	150
	Floor Mat	1340	750	
	Long Roll	520	90	70
	Short Roll	280	90	70
	Positioning Strap	660	80	
	Large Horseshoe	440	490	150
	Small Horseshoe	70	90	200

NDIS 05_122500821_0105_1_2

Strollers



Convoid EZ Rider Stroller

Convoid

- Ultra lightweight compact folding stroller
- 10° seat tilt allows correct positioning for children needing postural support, but also keeps them sitting upright
- Seating accessories include headrests, torso supports and tray
- Available as crash tested Bus Transport model
- Bus Transport model comes standard with tie-down system, h-harness, headrest extension, and foot straps
- Available in 4 different seat widths



LIFETIME WARRANTY

Code	Description	Seat Width	Seat Depth	Back Height	Seat to Footrest	Product Weight	SWL
PED477700	EZ Rider (as per script)	305 - 455	229 - 495	451 - 508	76 - 559	11.4 - 13.2kg	30 - 77kg

NDIS 05_122707821_0105_1_2

Convoid Cruiser Stroller

Convoid

- Lightweight compact folding stroller
- Set at 30° seat tilt to allow for more postural support
- Seating accessories include headrests, torso supports and tray
- Available as crash tested Bus Transport model
- Bus Transport model comes standard with tie-down system, h-harness, headrest extension, and foot straps
- Available in 5 different seat widths



LIFETIME WARRANTY

Code	Description	Seat Width	Seat Depth	Back Height	Seat to Footrest	Product Weight	SWL
PED477260	Cruiser (as per script)	255 - 455	152 - 533	533 - 630	130 - 580	8.6 - 13.6kg	30 - 77kg

NDIS 05_122707821_0105_1_2

Convoid Rodeo Tilt Stroller



- Versatile tilt-In-Space stroller
- Compact folding and lightweight design for easy transporting
- Available as crash tested Bus Transport model
- Available in 4 different seat widths
- Bus Transport model comes standard with tie-down system, h-harness, headrest extension, and foot straps

LIFETIME WARRANTY

Code	Description	Seat Width	Seat Depth	Back Height	Seat to Footrest	Product Weight	SWL
PED478195	Rodeo Tilt (as per script)	255 - 405	178 - 457	533 - 635	230 - 457	13.2 - 17.7 kg	30 - 63kg
NDIS 05_122707821_0105_1_2							

Medifab Bingo Evo Cobra



- Cobra indoor base with BINGO Evo seat offers an excellent indoor seating solution
- 45° of tilt and continuous adjustable gas spring for height adjustment
- Tilt and height adjustments can be made with the child in the seat
- Multiple accessories for proper positioning and support
- Customisable and easy-to-clean upholstery design
- Back support recline of 90° - 160° (90° - 135° for size 3)
- Growth adjustable in seat depth, width, back support height, and lower leg length
- Sleek and sporty appearance
- High-quality German-made product

1 YEAR WARRANTY

SWL 50KG

Code	Description	Seat Height	Tilt	Length	Width	Product Weight
WCA859680	Size 1, 1 - 5 years - indoor frame	279 - 610	45°	710	610	11.3kg
NDIS 05_122442871_0105_1_2						

Medifab Bingo Evo Stroller



- Modern, compact folding size and sporty design
- Multifunctional seat unit with quick front to rear facing change
- Centrally positioned seat positioning controls for easy operation
- Various adjustment options for perfect individual adjustability
- Numerous accessories for perfect individual care

1 YEAR WARRANTY

SWL 50KG

Code	Description	Width	Length	Seat Width	Seat Depth	Back Support Recline	Product Weight
PED491520	Size 2 (as per script)	680	920	230 - 350	220 - 400	90° - 160°	16.1kg
NDIS 05_122707821_0105_1_2							

Spex Bingo Stroller



- Back support shaping accommodates or corrects trunk spinal deformities, supported by multi-adjustable lateral supports and anterior chest harnesses
- Adjustable pre-ischial area can be conformed on-the-spot for varying anti-thrust needs, and strategic positioning pads can be added or removed from underneath the cushion
- Comfortable, padded forearm supports on quick-release hardware and includes feet positioning sandals

SWL 50KG

2 YEARS WARRANTY

Code	Description	Seat Width	Seat Depth	Back Height	Lower Leg Length
PED491300	Size 1 - TIS - 45° (as per script)	170 - 280	180 - 300	360 - 440	250 - 400
PED491354	Size 2 - TIS - 45° (as per script)	220 - 330	230 - 380	400 - 480	250 - 450

NDIS 05_122707821_0105_1_2

Firefly WeGo Stroller



- It reclines to 140°, providing support for your child to lie back for a snooze or to be changed, while easily fitting a GoTo Seat
- Comes with a full accessory pack, including a sun canopy, rain cover, cup holder, mosquito net, and foot muff, and features latex-free, breathable fabrics for outdoor use

SWL 30KG

2 YEARS WARRANTY

Code	Description	Height	Weight	Seat Depth	Seat Height	Seat Width	Footrest*
PED491200	Push Chair with Full Accessories Pack	1120	110	260	500	350	250 330
NDIS 05_122707821_0105_1_2							*Shortest Setting Longest Setting

Medifab Pixi Stroller



- Compact frame, can be quickly folded down to place in the boot of a vehicle
- Backrest recline, allows the child to relax or sleep
- Manoeuvrable, ball bearing castors for added convenience
- Robust design, excellent reliability

1 YEAR WARRANTY

Code	Description	Width	Seat Width	Seat Depth	Back Angle	SWL	Product Weight
PED491350	Size 1 (as per script)	585	280	250	90° - 125°	35kg	12.5kg
PED491352	Size 2 (as per script)	635	330	400	90° - 123°	50kg	15kg

NDIS 05_122707821_0105_1_2

Dejay Special Tomato EIO Pushchair



- Designed for mild to moderate positioning needs, suitable for early intervention to preadolescence
- Padded fabric seat, back laterals, and head cushion for head-to-toe comfort, along with a standard angle adjustment foot well
- 5-point comfort wrap harness ensures safety, and the seat has full recline capacity for added convenience

SWL 40KG

2 YEARS WARRANTY

Code	Description	Seat Width	Seat Depth	Back Height	Handle Height to Floor	Product Weight
PED491640	EIO Pushchair (as per script)	356	305	609	762 - 1092	10kg
NDIS 05_180909811_0103_1_2						

Medifab Kool Stride Running Stroller

medifab



- Serious sports buggy for outdoor use
- Sealed bearings on all wheels
- Stainless steel spokes with aluminium hubs and rims
- Adjustable handle and patented fold system
- Quick-release wheels for easy transport

1 YEAR WARRANTY

SWL 60KG

Code	Description	Seat Width	Seat Depth	Back Height	Width	Shoulder Height	Product Weight
PED491305	Kool Stride Running Stroller	430	300	600	680	600	15kg

NDIS 05_122707821_0105_1_2

Apex Mobility XRouter Multi Functional Outdoor Stroller



- Certified medical stroller, registered with TGA
- Complies with AS EN 12182:2025 for users up to 120 kg
- Lightweight and easy to operate, suitable for adults up to 120 kg
- Ideal for outdoor terrains including city streets, beach, and snow
- Can be attached to an adult bicycle for added versatility
- Available in three sizes for children to adults

2 YEARS WARRANTY

Code	Description	Max. User Height	Height	Width	Hip Width	SWL
PED491610	Small (as per script)	1219	870	630	410	68kg
PED491620	Medium (as per script)	1524	950	630	410	82kg
PED491630	Large (as per script)	1879	1030	630	410	120kg

NDIS 05_044839811_0103_1_2

Ottobock Kimba Outdoor Stroller

ottobock.



Without outdoor mobility base option

- Certified medical stroller, registered with TGA
- Complies with AS EN 12182:2025 for users up to 120 kg
- Lightweight and easy to operate, suitable for adults up to 120 kg
- Ideal for outdoor terrains including city streets, beach, and snow
- Can be attached to an adult bicycle for added versatility
- Available in three sizes for children to adults

2 YEARS WARRANTY

SWL 50KG

Code	Description	Seat Depth	Seat Width	Backrest Height	SWL
PED478433	Kimba Seating Unit - seat only (as per script)	190 - 400	200 - 350	410 - 700	40kg
PED478432	Kimba Outdoor Mobility Base - frame only (as per script)				

NDIS 05_122707821_0105_1_2

Hippocampe All Terrain Beach Wheelchair



- Double width wheels for use on sand
- Optional balloon tyres enable the chair to float
- Corrosion resistant lightweight frame - 14 kg
- Can be fitted with skis for use in the snow
- Removable push bar for easy steering
- Can be folded for transport

2 YEARS WARRANTY
Frame & Upholstery

SWL 130KG

Code	Description	Floor-to-seat Height	Frame Depth	Seat Width	Width
PED491740	Hippocampe All Terrain Beach Wheelchair (as per script)	Rear: 340 Front: 350	720 - 1070	420	670

NDIS 05_044839811_0103_1_2

Walking Aids



Nurmi Neo Walking Aid

ottobock.

- Innovative walking aid for children
- Positioned behind the child for open directional walking
- Lightweight, folding height-adjustable frame
- Features reverse roll locks and swivelling front casters
- Customisable options for individual needs



Code	Description	Width	Length	Product Weight	SWL	Warranty
PED479349	Size 1 (as per script)	630	650	5.5kg	25kg	2 years
PED479350	Size 2 (as per script)	680	760	6kg	40kg	2 years
PED479351	Size 3 (as per script)	720	830	6.5kg	55kg	2 years

NDIS 05_120603831_0103_1_2

Nimbo Lightweight Walking Aid

- Lightweight aluminium frame
- Tool-free height adjustment
- Soft rubber wheels for all surfaces, preventing backward sliding
- One-directional rear wheels with aluminium ratchet and pin mechanism
- Latex-free handles with flanged ends
- Folds easily for storage and transport
- Accessories available on request



Code	Description	Width	Length	Product Weight	SWL	Warranty
PED480780	Tyke, Size 1, 390 - 500mm, Goldenrod Yellow	530	365 - 470	4kg	45.5kg	2 years
PED480790	Junior, Size 2, 480 - 640mm, Lime Green	555	480 - 635	5kg	45.5kg	2 years
PED480800	Youth, Size 3, 500 - 770mm, Cornflower Blue	635	580 - 770	5.5kg	91kg	2 years
PED480805	Young Adult, Size 4, 710 - 910mm, Midnight Blue	645	720 - 910	6kg	91kg	2 years
PED480807	Adult, Size 5, 920 - 1040mm, Flame Red	740	910 - 1040	7.72kg	125kg	2 years
PED480810	Nimbo Seat for Lightweight					2 years

NDIS 05_120603831_0103_1_2

Leckey MyWay+ Gait Trainer

LECKEY^o



- Open frame design enables children to touch, feel, learn and play at their own pace
- MyWay castors are multifunctional to assist children with their individual needs
- Can be fine tuned to achieve the optimum position for mobility for each child
- Open frame also allows therapists to work closely with the child on their lower limbs for therapeutic benefits
- Frame is foldable, facilitating convenient storage and transportation

3 YEARS WARRANTY

Code	Description	User Height	Seat Height	SWL
PED479600	Size 1, Ages 1 - 5 years (as per script)	635 - 920	220 - 440	25kg
PED479610	Size 2, Ages 4 - 11 years (as per script)	870 - 1190	400 - 610	50kg
PED479615	Size 3, Ages 9 - 16 years (as per script)	1070 - 1490	550 - 810	80kg

NDIS 05_120603831_0103_1_2

Ottobock Kidwalk

ottobock.



- Partially weight bearing, posturally supportive mobility system which encourages mobility and moves with the child
- Easy, tool-free set-up and adjustment
- Child walks unaided while also keeping hands free

2 YEARS WARRANTY

Code	Description	User Height	Width	Length	SWL
PED479355	Size 1, Ages 1 - 7 Years (as per script)	610 - 1220	640	825 - 875	30kg
PED479356	Size 2, Ages 5 - 12 Years (as per script)	1010 - 1520	700	900 - 950	40kg

NDIS 05_120603831_0103_1_2

Medifab Flux Walker

medifab



- Promotes natural walking movements
- Optimal design for upright posture
- Folding mechanism with safety lock for easy transport and storage
- Utilizes and promotes residual functions efficiently
- Sturdy yet lightweight frame construction
- Cross member in rear frame for safety
- Reliable reverse-roll locks
- Adjustable to grow with your child

1 YEAR WARRANTY

SWL 75KG

Code	Description	Width	Length	Grip Height	Product Weight
PED470076	Size 2 (as per script)	690	810	510 - 880	6.8kg

NDIS 05_044839811_0103_1_2

Rifton Pacer Gait Trainer

rifton.



- Pacer accommodates all abilities and environments
- Configurable to meet every client's needs
- Options: dynamic or standard upper frames
- Bases: standard, utility, or treadmill/stability

1 YEAR WARRANTY

Code	Description	Length	Width	Height	Product Weight	SWL
PED470074	Rifton Pacer Gait Trainer (as per script)	572 - 1029	520 - 813	318 - 978	5.2 - 6.8kg	34 - 113.4kg

NDIS 05_044839811_0103_1_2

Apex Mobility Malte Outdoor Posterior Walker

APEX MOBILITY



- For extra safety: parking brake for stability; EVA rear wheels with drum brake; the optional tilt guard protects against tilting back
- For easy transport: the aluminium frame can easily be folded up
- For stability: the optional backwards stop prevents unwanted rolling back
- For outdoor adventure: special tyre dimensions for outdoor use; gait trainer with outdoor approval
- For extra comfort: the padded bench can be easily folded back

2 YEARS WARRANTY

Code	Description	Length	Width	Handle Height	Seat Height	Product Weight	SWL
PED491710	Size 1	660	560	340 - 580	330	7.1kg	30kg
PED491720	Size 2	750	660	510 - 750	450	8.3kg	50kg
PED491730	Size 3	820	730	680 - 940	600	9.8kg	70kg

NDIS 05_044839811_0103_1_2

Rebotec Walk-On Kids Walking Frame

REBOTEC



- Smooth and stylish powder coated aluminium finish
- Non-slip hand-grips for safety and confidence
- Height adjustable walking frame
- 5x Height Adjustments

SWL 150KG

2 YEARS WARRANTY

Code	Description	Seat Depth	Hip Width	Handle Height	Product Weight
PED479618	Small	430	520	700 - 810	3.1kg
PED479619	Extra Small	430	520	620 - 710	3kg

NDIS 05_120603831_0103_1_2

Kids Forearm Crutches


- Just for kids, these light weight yellow Forearm Crutches are height adjustable

SWL 100KG

1 YEAR WARRANTY

Code	Description	Handle Height	Ferrule
WAC691700*	Child, 4 - 7 years	406 - 559	18
WAC691600	Child, 6 - 10 years	457 - 656	18

*NDIS 05_120603831_0103_1_2



Underarm Crutches Red Dot



- Patented clip-lock adjustable handle and push button foot pieces for fast, precise measuring and fitting, 22mm ferrules
- High density rubber cushions and handgrips providing long wear and greater comfort
- Lightweight aluminium frame, child weighing 1.1kg each, and youth-adult tall weighing 1.25kg each


1 YEAR WARRANTY

SWL 136KG

Code	Description	User Height	Crutch Height
WAC692300*	Child	1220 - 1370	787 - 940
WAC692500**	Youth	1320 - 1570	940 - 1140

**DVA AP03 *NDIS 03_120000911_0105_1_1

Aspire Underarm Crutches Wing Nut



- Easy to use Wing-Nut height adjustment
- High density rubber cushions and handgrips providing long wear and greater comfort
- Lightweight aluminium frame
- 22mm ferrules

1 YEAR WARRANTY

SWL 80KG

Code	Description	User Height	Crutch Height
WAC692503	Youth - Small	1370 - 1575	940 - 1145

DVA AP03 NDIS 03_120000911_0105_1_1

Firefly Upsee



- Creates the possibility of social and emotional development by letting children interact with their peers at eye-level
- Upsee can be used to work on a variety of mobility and ability skills, from weight-bearing and walking to improving head control
- When your child is motivated to explore and move and weight-bear, it creates a real possibility for long-term development

1 YEAR WARRANTY

Code	Description	Shoulder to Hip	Chest Circumference	Lower Hip Circumference	SWL
PED490700*	Therapist Kit				
PED490703	Extra Small - Purple	280 - 320	500 - 570	500 - 570	15kg
PED490704	Extra Small - Green	280 - 320	500 - 570	500 - 570	15kg
PED490705*	Small - Blue	310 - 360	520 - 600	520 - 600	15kg
PED490706	Small - Pink	310 - 360	520 - 600	520 - 600	15kg
PED490710*	Medium - Blue	350 - 390	570 - 650	580 - 670	20kg
PED490715	Medium - Pink	350 - 390	570 - 650	580 - 670	20kg
PED490720	Large - Pink	370 - 420	630 - 710	530 - 730	25kg
PED490725	Large - Blue	370 - 420	630 - 710	530 - 730	25kg

NDIS 05_122500821_0105_1_2 *05_120603831_0103_1_2

The Upsee kit is made up of three parts: an adult hip-belt, a child harness and shared sandals

Manual Wheelchairs



Aspire Rehab RX2 JNR Wheelchair

TIS



- Provides increased postural support and pressure care for children
- 26° tilt and 30° backrest recline for weight redistribution and skin integrity
- Highly adjustable seat depth, width, armrests, and headrest
- Multi-layered foam support with breathable, multi-stretch covers
- Contoured backrest and seat cushion with lumbar and lateral supports

SWL 180KG

2 YEARS WARRANTY Frame

1 YEAR WARRANTY Components

CRASH TESTED APPROVED



Code	Description	Overall Width	Seat Width	Seat Depth	Seat Height	Backrest Height	Armrest Height	Min. Chair Weight
MWS449750	Rehab RX2 JNR - 320mm	600	320	440 - 490	400 - 475	460 - 710	310 - 460	29.5kg

NDIS 05_122203811_0105_1_2

Aspire Evoke 2 Wheelchair JNR & Youth

FOLDING



- Sway-away, height adjustable and removable armrests
- Depth adjustable seating with rear angle adjustability
- Easy access push-to-lock brakes
- Quick release 20" solid rear wheels for easy transport and storage
- Anti-tippers included as standard

SWL 150KG

2 YEARS WARRANTY Frame

1 YEAR WARRANTY Components

CRASH TESTED APPROVED

OVERLAY COMPATIBLE



Code	Description	Seat Width	Overall Width	Seat Depth	Seat Height	Backrest Height	Armrest Height	Min. Chair Weight
MWS449560	Evoke 2 JNR - 300mm	300	510	300 - 375	360 - 410	430 - 490	205 - 225	14kg
MWS449570*	Evoke 2 Youth - 350mm	350	560	350 - 425	440 - 480	430 - 490	205 - 225	14.8g

NDIS 05_122203107_0105_1_2 *05_122203811_0105_1_2

Permobil TiLite Pilot Wheelchair

RIGID



- Ultra-lightweight design
- Carbon fiber frame panels make it 50% lighter and a breeze for kids to navigate
- Built-in forward growth adjustment
- To grow with your child and maintain optimal rear wheel access
- Clean frame design
- Features a seat pan and side panels that provide easy access to the adjustment bolts, making changes simple and intuitive

Range of frame colours available

ACID GREEN METALLIC

NEON PINK

SAPPHIRE

RED PEPPER

OCEAN BLUE METALLIC

ELECTRIC PLUM

BLACK

ULTRAVIOLET

SWL 54KG

Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Product Weight
MWC427483	TiLite Pilot Wheelchair (as per script)	203 - 355	228 - 381	76 - 368	4.9kg

NDIS 05_122203821_0105_1_2

Permobil Progeo Tekna Tilt

TILT



- Lightweight, folding design
- Shortened folding frame and integrated front end for increased manoeuvrability
- Performance-oriented options for seating and positioning needs
- Open frame design for compact transportation
- Provides a performance-oriented ride
- Aluminium Frame

SWL 125KG

Code	Description	Seat Width	Seat Depth	Seat Height	Backrest Height	Backrest Angle
MWC451130	Progeo Tekna Tilt (as per script)	270 - 480	350 - 480	420 - 485	245 - 470	95° - 120°

NDIS 05_122203821_0105_1_2

Permobil Progeo Joker

RIGID



- Combines benefits of rigid frame with adjustable seating and drive performance
- Provides balanced and responsive ride
- Preferred choice for users with changing needs due to its adjustability
- Aluminium Frame

SWL 125KG

Code	Description	Seat Width	Seat Depth	Seat Height (Front / Rear)	Backrest Height	Backrest Angle
MWC451180	Progeo Joker (as per script)	330 - 480	350 - 500	400 - 590 / 380 - 460	245 - 470	74° - 94°

NDIS 05_122203821_0105_1_2

permobil

Permobil TiLite Twist

RIGID

- Tru-Fit System for infinite adjustability
- Adjustable Camber Mount with up to 3.5" rear seat height adjustment
- Speed loader adjustable front caster housing for quick fork angle adjustments
- Wide range of colours available

SWL 54KG

Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Product Weight
MWC427472	TiLite Twist Wheelchair (as per script)	203 - 381	203 - 406	76 - 368	4.9kg

NDIS 05_122203821_0105_1_2

permobil



Ottobock Start Junior Wheelchair

FOLDING

- Lightweight, durable pediatric chair with resilient upholstery
- Full range of child-specific accessories
- Uniform caster forks for all wheel sizes
- Compact folding with lock for easy storage
- Adjustable footrest and available in various colours

SWL 90KG

Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Product Weight
MWC427410	Start Junior 30cm (as per script)	280 - 380	240 - 400	395 - 520	13.5kg
MWC427420	Start Junior 35cm (as per script)	280 - 380	240 - 400	395 - 520	13.5kg

NDIS 05_122203821_0105_1_2

ottobock.



Ottobock Avantgarde 4 Wheelchair

FOLDING

- Lightweight and self-propelled for children, adolescents, and small adults
- Folding aluminium alloy frame
- Sturdy aluminium side panels
- Abduction frame for correct pelvic position
- Adjustable curved push handles and knee lever brake

SWL 180KG

Code	Description	Seat Width	Front Seat Height	Rear Seat Height
MWC427365	Avantgarde 4 DV (as per script)	320 - 520	380 - 540	360 - 520
MWC451570*	Avantgarde 4 DS (as per script)	320 - 520	380 - 540	360 - 510
MWC427113	Avantgarde 4 XXL2 (as per script)	460 - 620	460 - 530	410 - 530
MWC451580*	Avantgarde 4 DF-S (as per script)	320 - 520	380 - 540	360 - 510
MWC451590*	Avantgarde 4 DF Version (as per script)	320 - 520	380 - 540	360 - 510

*NDIS 05_122203821_0105_1_2

ottobock.



Quickie 2 & Quickie 2 Lite Wheelchairs

FOLDING



- Folding ultra lightweight manual wheelchair featuring a modular design
- Versatile, modular and lightweight
- Over 64,000 possible rehab configurations with TruFit programme

Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Product Weight	SWL
MWC427235	Quickie 2 Lite (as per script)	356 - 508	356 - 508	360 - 560	8kg	120kg
MWC427230	Quickie 2 (as per script)	305 - 508	305 - 457	360 - 560	13kg	136kg

DVA AP19 NDIS 05_122203821_0105_1_2

Ki Mobility Little Wave XP Klik Wheelchair

RIGID



- Index system for easy growth positioning and adjustment options
- Built in adjustment
- 5th wheel option for extra balance
- Swing away arms



Code	Description	Seat Width	Seat Depth	Seat Height	Backrest Height
MWC427554	Little Wave XP Klik Wheelchair (as per script)	200 - 460	200 - 460	290 - 510	229 - 457

NDIS 05_122203821_0105_1_2

Ki Mobility Little Wave XP Wheelchair

RIGID



- Ultra-light rigid aluminium wheelchair
- Crescent CG system keeps the side frame free of excess hardware while providing infinite adjustability
- Tubular component system provides individual seat height settings



Code	Description	Seat Width	Seat Depth	Front Seat Height	Rear Seat Height
MWC427039	Little Wave XP Wheelchair (as per script)	250 - 410	250 - 410	290 - 520	290 - 510

NDIS 05_122203821_0105_1_2

Ki Mobility Rogue 2

RIGID



- Maximum adjustment and growth capability for optimal fit, function, and ride
- Integrated swing-away armrests adjust with the backrest
- Durable, ultra-lightweight 7000 series aluminium frame
- Adjustable dynamic 5th wheel for range and spring rate
- Spinergy backrest release cable
- Max Performance spoke wheels for style and performance at an affordable price



Code	Description	Seat Width	Seat Depth	Front Seat Height	Rear Seat Height
MWC427048	Rogue XP Wheelchair (as per script)	250 - 460	250 - 510	340 - 520	290 - 510

NDIS 05_122203821_0105_1_2

Zippie X'Cape Wheelchair

FOLDING



- Rigid wheelchair performance and folding wheelchair convenience using the New Xlock Crossbrace
- Infinite footrest adjustment using the new z-finity system
- Versatile for limitless seating and accessory options



Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Product Weight
MWC427501	Zippie X'Cape	203 - 406	203 - 457	330 - 482	8.6kg

NDIS 05_122203821_0105_1_2

Zippie Zone Wheelchair

RIGID



- Ultra Lightweight
- Built in Growth Kit and Truefit growth Programme for growth capacity
- Inset backrest to facilitate reach to wheels
- Angle adjustable footplates
- Ultra lightweight aerospace aluminium
- Colourful options



Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Product Weight
MWC427487	Zippie Zone Wheelchair (as per script)	200 - 460	200 - 480	360 - 480	6.5kg

NDIS 05_122203821_0105_1_2

Zippie Iris Wheelchair

TIS



- Rotation-in-space wheelchair
- Offers a superior combination of positioning, manoeuvrability, and increased portability with the XLOCK® folding crossbrace option
- Shortest Possible Wheelbase
- This rotational system provides the shortest possible wheelbase without compromising stability to ensure a smooth and controlled rotation
- Comes with the True Fit Programme
- A complete growth package for free
- Cutting-edge rotation-in-space technology designed just for kids
- This tilting wheelchair uses patented rotation-in-space technology to rotate the seat frame around the user's centre of gravity



Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Product Weight
MWC427500	Zippie Iris Wheelchair (as per script)	254 - 457	254 - 457	320 - 500	17 kg

NDIS 05_122203821_0105_1_2

Ki Mobility Little Wave Flip XP Wheelchair

TIS



- Pediatric tilt-in-space solution designed to look more like a stroller with the technology of a wheelchair
- Over 2,000 colour combinations
- Simple “flip” and fold like a stroller, transforms into a small package that can fit in the trunk of the tiniest car



Code	Description	Seat Width	Seat Depth	Seat Height	Tilt Range
MWC427553	Little Wave Flip XP Wheelchair (as per script)	250 - 460	300 - 510	340 - 510	0 - 45°
PED491505	Little Wave Flip Folding Base (as per script)				

Ki Mobility Little Wave Arc Wheelchair

TIS



- Little Wave Arc is the lightest Pediatric Folding wheelchair
- Repositioning - Tilting the occupant has been made easier through an Arc shaped roller guide system helps translate the Centre of Gravity throughout the tilt process
- Attendant Wheel Lock - Easy access to the wheel lock for the attendant
- Drum Brakes
- Tilt options - Tilt the occupant up to 45°
- Folding recline back & Folding Adjustable Handle
- Folding Height Adjustable Rigidizer Bar



Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Product Weight
MWC427015	Arc Little Wave Wheelchair (as per script)	250 - 460	300 - 510	370 - 480	10kg
NDIS 05_122203821_0105_1_2					

Invacare Freedom NXT Wheelchair

TIS



- The NXT is the most rigid, folding, paediatric tilt-in-space manual wheelchair on the market
- It is strong and durable, yet compact and lightweight
- The NXT ensure both the child and carer remain comfortable by providing a smoother ride and a user friendly solution
- When the growth spurt hits and the frame won't fit, the NXT adjusts easily for a good fit today and tomorrow



Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Product Weight
MWC449180	Freedom NXT Wheelchair (as per script)	255 - 510	305 - 510	320 - 500	10kg
NDIS 05_122203821_0105_1_2					

Ottobock Kidevo Mini

TIS



- Perfect for early childhood
- Customisable seat width and depth starting from 16 cm
- Tailor-made to grow with the child in width, depth, and height
- No fixed sizes, maximum adaptability
- Adjustable centre of gravity for optimal driving behaviour
- Available with a tilting function (kidevo mini.t)



Code	Description	Seat Width	Seat Depth	Seat Tilt	Rear Wheel Size	Product Weight
MWC427572	Kidvevo Mini (As per Script)	160 - 320	160 - 350	-5° - 40°	508, 558.8 & 609.6	7.4kg
NDIS 05_122203821_0105_1_2						

Power-Assist Systems



Smoov One by Alber

- Unobtrusive, high quality design
- Robust, rustproof housing and rustproof fork
- Integrated lithium ion battery pack for a range up to 20km, and speed of up to 10km/h
- Ergonomic carry handle
- USB-C socket to charge the control unit or smartphone



Code	Description	Length	Width	Height	Product Weight
MWC449500	Smoov One by Alber	385	145	340	7.2kg
NDIS 05_122203821_0105_1_2					


Alber E-Fix E35 & E36

- Lightweight and compact add-on drive can be quickly removed from the wheelchair without the need for tools
- The modern, flat design of the control unit offers simple operation and perfect ergonomics
- Range of up to 16km with 6.0 Ah battery pack
- Weight of the person of up to 120kg
- Electronic automatic charge and switch off mechanism




Code	Description	SWL
MWC449520	Alber E-Fix E35	120kg
MWC449525	Alber E-Fix E36	160kg
NDIS 05_122409871_0105_1_2		

Benoit Systemes Light Drive2



- Instant installation
- Light and easily transportable
- Adapts to all wheelchairs
- Removable
- Discreet and unobtrusive design
- Can be transported by plane (12Ah Lithium Battery)
- Colour of your choice
- Tackles almost any terrain
- Indoors/outdoors




Code	Description	Charging Time	Range	Product Weight
WCA858507	Light Drive2 - Folding	5 - 6 hours	14 - 20km	4kg
WCA858512	Light Drive2- Rigid	5 - 6 hours	14 - 20km	4kg

NDIS 05_122442871_0105_1_2

Alber ViaGo Power Assist



- Provides extra power assistance for manual wheelchairs
- Compact and powerful design
- Easy to use and quick to install
- Ensures effortless manoeuvrability and confidence
- Helps with challenging terrains and inclines
- Supports independent mobility for various excursions



Code	Description	Suits
WCA858501	ViaGo Power Assist	Rigid and Folding Chairs

NDIS 05_122442871_0105_1_2

Permobil SmartDrive MX2+ Push Tracker



- Provides great freedom and power to manual chair users
- Ability to travel up hills and over thick carpet
- Lightweight and easy to hook on to manual wheelchairs
- Folding wheelchair compatibility, simply pull back on the spring-loaded flanges to quickly attach or remove your SmartDrive MX2
- PushTracker app available for iOS and Android devices
- PushTracker monitor with wristband allows easy tracking of day to day propulsion metrics - including push count, coast time, and more
- Includes battery and mounting hardware



Code	Description	Length	Width	Height	Product Weight
WCA858508	Smart Drive MX2 + Push Tracker	389	141	242	5.7kg

NDIS 05_122442871_0105_1_2

Quickie iXpress



- Japanese designed and made, highly reliable electronics
- Selectable speeds for easy transitions between environments
- Easy to use ergonomic joystick controller
- Full function attendant control option with innovative glide to steer function
- Programmable for more individual need
- Brush-less direct drive wheels add power to our manual wheelchairs
- Powered Wheels are removable and can be replaced with standard wheels



Code	Description	Charging Time	Range	Product Weight
WCA858506	iXpress Power Add On (as per script)	3 hours	16km/h	17.4kg

DVA AP25 NDIS 05_122442871_0105_1_2

Power Wheelchairs



Permobil Explorer Mini

- Regulatory-cleared powered mobility solution for children 12-36 months with mobility impairments
- Integrated table and backrest for support and freedom of movement
- Adjustable saddle seat supports hip development, functional sitting, and incremental standing
- Ergonomic design with multiple weight-bearing surfaces for safe, stable upright postures
- Promotes strength, endurance, and postural control
- Midline proportional joystick with bright yellow ball handle aids visual development, reaching, grasping, and releasing



Code	Description	Height	Max. User Height	Floor-to-seat Height	Ground Clearance	Base Width x Length	Max. Speed
PED491500	Explorer Mini (as per script)	730 - 930	1000	307 - 550	250	490 x 635	2.4 kph

NDIS 05_122203911_0105_1_2

Permobil K450 Power Wheelchair

- All-wheel suspension for a smooth ride indoors and outdoors
- Power seat-to-floor function for easy ground access, story time, transfers, and play
- Adjustable seat height from 3" to 26" for easy access to tables, counters, and low objects
- R-Net electronics with power functions: seat elevator, seat-to-floor, and optional power tilt
- Modular seating system with centre mount adjustable legrests and free growth kit
- Large flat tyres for improved traction and reduced maintenance



Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Base Width x Length	Max. Speed
PWC652554	K450 (as per script)	250 - 410	250 - 460	360 - 460	610 x 890	8 kph

NDIS 05_122306811_0105_1_2

Permobil Koala Power Wheelchair



- Designed to grow with children, ideal for toddlers
- Scaled for kids but built with rugged motor and front-wheel-drive technology
- Low seat-to-floor height for peer interaction and access to various surface heights
- Integrated elevation for easy access to tables and desks
- Tilt feature for improved stability and rest
- Optional features: front bumper and remote stop for added safety



Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Base Width x Length	Max. Speed
PWC652563	Koala (as per script)	250 - 360	250 - 360	375	585 x 810	6.5kph

NDIS 05_122306811_0105_1_2

Ottobock Skippi



- The Skippi is a fun, exciting, paediatric powerchair
- While our standard seat can be made as small as 24cm, it can also be paired with the much loved Squiggles and Mygo seating systems from LECKEY
- The Skippi Plus model allows for specialty controls, ensuring even the most complex clients' needs can be met



Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Base Width x Length	Max. Speed
PWC652265	Skippi (as per script)	300 - 380	280 - 380	430	570 x 850	6kph

Permobil K300 PS JNR



Zippie Q300 M Mini Kids



- Ultra-Compact Design, 520 mm width and 1120 mm turning circle for smooth navigation indoors and outdoors
- All-Wheel suspension, the patented suspension system for stability and safety on various terrains
- Powerful performance, equipped with super strong 4-pole motors
- Drive range of up to 25km on a single charge, ideal for extended exploration
- Customisable options numerous powered and styling options to match your child's personality
- Ideal for active children who need a compact, powerful, and customisable power chair for both indoor and outdoor adventures



Code	Description	Base Width	Length	Drive Wheels	Width	Depth	Max. Speed
PWC655620	Zippie Q300 M Mini Kids	520	1100	12.5"	10 - 16"	10 - 19"	10km/h

NDIS 05_122306811_0105_1_2



Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Base Width x Length	Max. Speed
PWC652500	K300 PS JNR (as per script)	280 - 410	250 - 460	420	610 x 1020	10kph

Zippie Z500 Mini



Sunrise Medical Q500 Mini Power Chair



- Available in Mid, Forward, Rear and Hybrid Wheel Drive
- Choice of rehab seating systems provide support and positioning to users from mild to medium clinical needs
- Robust outdoor chair with advanced electronics
- Compatible with Sedeo Pro and Sedeo Lite seats



Code	Description	Width	Length	Floor-to-seat Height	Max. Speed
PWC655330	Rear Wheel Drive (as per script)	640	1180	415 - 490	10kph
PWC655335	MID Wheel Drive Sedeo Lite	610 - 620	1280	415 - 490	10kph
PWC655340	MID Wheel Drive Sedeo Lite (as per script)	610 - 620	1280	415 - 490	10kph
PWC655350	Front Wheel Drive (as per script)	640	1340	415 - 490	10kph

- Ultra-compact design, 52cm wide, perfect for navigating tight indoor spaces such as narrow hallways and between desks
- Mid wheel drive base ensures intuitive handling and smooth turns
- Highly adjustable and customisable seating designed to meet the unique needs of active kids, promoting comfort and proper posture
- Powerful motors equipped with high-performance 10 kph motors, capable of handling tough terrain effortlessly
- Long-lasting batteries: provides extensive range on a single charge, ensuring your child can play and explore throughout the day without interruptions.
- Ideal for active children who need a reliable, agile, and customisable wheelchair for both indoor and outdoor use



Code	Description	Seat Width	Seat Depth	Floor-to-seat Height	Base Width x Length	Max. Speed
PWC655630	Zippie Z500 Mini	254 - 406	254 - 482	430 - 490	530 x 860	10km/h

NDIS 05_122306811_0105_1_2



Quantum Q6 Edge 3.0

QUANTUM



- Designed for an advanced power chair experience
- Upgraded SRS (Smooth Ride Suspension) for comfort and stability
- Optional iLevel® offers 7.2 km/h speed and up to 12" of power adjustable seat height
- Standard features: front and rear LED fender lights, USB charger, and drive wheel colour accents



Code	Description	Floor-to-seat Height	Max. Speed
PWC652223	Q500 Power Chair (as per script)	419 - 495	7.24kph
DVA AP16 NDIS 05_122306811_0105_1_2			

Quantum Edge 3 Stretto

QUANTUM



- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 3.5 mph)
- Fender lights and USB charger included
- Toolless removal of rear shroud with battery access door
- TRU-Balance® 3 seating configurations
- Available with Accu-Trac Advanced Tracking Technology



Code	Description	Floor-to-seat Height	Max. Speed
PWC652238	Edge 3 Stretto (as per script)	406 - 464	6.25kph

Quantum 4 Front

QUANTUM



- CASE (Castor Angle Sensor Encoder) technology
- Unmatched driving performance and ride quality
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Q-Logic 3 EX Advanced Drive Controls
- Advanced stability for smooth obstacle transitions
- Quiet ride technology



Code	Description	Floor-to-seat Height	Max. Speed
PWC652222	4 Front Power Chair (as per script)	387 - 500	9.6kph

Magic Mobility Frontier V6 Compact 40

magic
Mobility



- Narrow footprint for ultimate manoeuvrability in tight indoor spaces
- Nearly two decades of proven reliability within the Frontier family
- Easily handles kerbs and navigates potholes
- Enhanced freedom for use in and around streets
- The right wheelchair for everyday off-road recreation, urban outdoor living and working, roomy indoor spaces, and tight indoor spaces



Code	Description	Base Footprint	Floor-to-seat Height	Max. Speed
PWC655050	Magic Mobility Frontier V6 Compact 40 (as per script)	965 x 655	435	10kph

Permobil M300 PS JR Paediatric Power Wheelchair

permobil



- Tight turning radius allows users to maneuver through your home with ease, thanks to the M300 PS Jr's tight turning radius
- CATLink suspension system ensures a smoother ride over rough terrain by increasing traction and preventing high-centreing
- Expandable seat pan for a custom fit, allowing the chair to grow with the child
- Configurable legrests help maintain proper positioning as the user grows
- Versatile armrests with unique rotation and flip-back movements facilitate easier transfers and added comfort



Code	Description	Length	Width	Seat Width	Range	Max. Speed
PWC652540	M300 PS JR (as per script)	930	610	280 - 400	20km	10km/h
NDIS 05_122306811_0105_1_2						

Wheelchair Backrests, Headrests & Cushions

At Aidacare, we stock an extensive range of Wheelchair Backrests, Headrests and Cushions suitable for users with minimal to complex needs. For specialised needs or inquiries about our complex mobility and positioning options, please contact your Aidacare Complex Rehab consultant. To view our range of Wheelchair Backrests, Headrests and Cushions online, scan the QR codes below.



Wheelchair Backrests



Wheelchair Headrests



Cushions



Head Supports



Hensinger Head Support



- Hensinger features a soft foam wrap around the collar to support the head by using base support areas at the occiput and jaw line
- Hensinger head support must be used in conjunction with a chest support
- This head support promotes proper head positioning resulting in improved breathing, swallowing and simplified feeding
- Hensinger enhances more social interaction and awareness during daily activities
- Revolutionary in the world of cervical collars, this soft foam wraparound collar is designed to brace the head by the jawline
- The Hensinger Head support can help maintain mid-line head positioning. Proper head positioning facilitates improved breathing, swallowing, and simplified feeding
- It also aids in social interaction and awareness during classroom activities
- The Hensinger is built with a strong yet pliable reinforcement material around the outside of the support
- A simple quick-release buckle allows easy on/off
- The Hensinger Head Support can be used with a chest support or mounted to a seating system
- Soft covers are included

Code	Description	Inner Circumference
FBS279900	Hensinger Head Support - X-Small	250
FBS279901	Hensinger Head Support - Small	300
FBS279902	Hensinger Head Support - Medium	350
FBS279903	Hensinger Head Support - Large	400

NDIS 05_060315111_0135_1_2

Night Time Positioning



Medifab Jenx Dreama Sleep System

- Dreama provides 24/7 postural management with a pressure-reducing mattress and customisable support pads for complex postures, complementing conventional positioning
- It supports muscle tone management, reduces tonal spasms, and prevents musculoskeletal deformities
- Ideal for all ages, Dreama enhances skeletal development and quality of life through consistent high-quality positioning day and night



Code	Description	Footprint	Height	Cushion Size	No. of Mattress Units
PED491360	Cot Size (9 months - 4 years)	560 x 1200	170	Single Mattress: 560 x 95 x 150 Double Mattress: 560 x 195 x 150	Single: 8 Double: 2
PED491362	Single Size (9 months - Adult)	900 x 1900	170	Single Mattress: 1200 x 95 x 150 Double Mattress: 1200 x 95 x 150	Single: 15 Double: 2

NDIS 05_181224811_0103_1_2

Medifab Simple Stuff Works Packages

- The Simple Stuff Works way provides high quality education and therapeutic positioning equipment with an aim to eradicate the preventable problems associated with body shape distortion




Code	Description
PTA534400	Cushion - Positioning - Simple Stuff Work (as per script)

NDIS 05_122506111_0105_1_2

HospEquip Symmetriket Range



The 24-hour Postural Care Programme integrates learning to lie straight, crucial for preventing or correcting body shape distortion. This programme includes gradual progression with adjustable elements suitable for all ages and bed positions, accommodating various mattress types like airflow, dynamic, turning, and hybrid. It supports multiple positions (supine, side lying, semi-side lying, and prone) and features Airmantle Technology for pressure relief, alongside a CoolOver Active Temperature Management system. Assessment and setup are efficiently completed in one visit, thanks to its lightweight and portable design, ensuring ease of handling and use.



- A 1 x Single Hook & Loop Flat Receptor Sheet (with guide marks)
- B 1 x Single Airmantle with CoolOver Cover and PU (wipedown) backing
- C 1 x Pair of Small Brackets
- D 1 x Pair of Medium Brackets
- E 1 x Pair of Medium/Large Brackets
- F 1 x Medium Roll Cushion including CoolOver Cover
- G 1 x Large Horseshoe Cushion including CoolOver Cover
- H CoolOver Bracket Pads - Small (x2), Medium (x2), Medium/Large (x2)
- I 1 x Large Symmetrisleep Travel Bag (Blue)

Code	Description
Available upon request	Adult (as per script)
NDIS 05_180909811_0103_1_2 Items could be purchased individually	



- A 1 x Single Hook & Loop Flat Receptor Sheet (with guide marks)
- B 1 x Single Airmantle with CoolOver Cover and PU (wipedown) backing
- C 1 x Pair of Small Brackets
- D 1 x Pair of Medium Brackets
- E 1 x Medium Roll Cushion including CoolOver Cover
- F 1 x Medium Horseshoe Cushion including CoolOver Cover
- G CoolOver Bracket Pads - Small (x2), Medium (x2)
- H 1 x Large Symmetrisleep Travel Bag (Blue)

Code	Description
Available upon request	Junior (as per script)
NDIS 05_180909811_0103_1_2 Items could be purchased individually	

Acheeva Learning Station



- Enables learning to lie straight, aiding in preventing or correcting body shape distortion
- Part of a 24-hour Postural Care Programme for maintaining body symmetry and joint protection
- Gradual progression with adjustable elements for all ages and bed positions
- Suitable for all mattress types: airflow, dynamic, turning, hybrid
- Supports supine, side lying, semi-side lying, and prone positioning
- Features Airmantle Technology for pressure relief
- CoolOver Active Temperature Management system
- Assessment and setup completed in one visit

Code	Description	Length	Width	Mattress Length	Mattress Width	Power Height Range	Power Tilt-in-space	Back Rest Raise/Lower	Product Weight	SWL
PED491364	Student	1650	700	1550	600	760 - 1160	+/-20°	+50°/-10°	73kg	100kg
PED491366	Graduate	2180	700	2080	600	760 - 1160	+/-20°	+50°/-10°	75kg	120kg

NDIS 03_120000911_0105_1_1

Medifab Care Wave Cushions



Cylindrical	Semi-Fowler	Decubitus	Heel	Delta			
<ul style="list-style-type: none">• Body positioning system for lying and seating support for infants and adults• Various cushion shapes for different positioning needs• Unique micro-particle padding for flexibility, molding effect, and comfort• Outer cover: Waterproof, breathable Polyamid bi-elastic fabric with antibacterial treatment• Filling: 100% polystyrene micro-particles					Code	Description	One Size
					PTA535075	Cylindrical (as per script)	700 x 170
					PTA535095	Semi-Fowler (as per script)	700 x 500 x 140
					PTA535080	Decubitus (as per script)	300 x 1800
					PTA535090	Heel (as per script)	900 x 180
					PTA535085	Delta (as per script)	550 x 300 x 300
					NDIS 05_122506111_0105_1_2		

Medifab Poz 'in' Form Lying Positioning System



- Proper positioning prevents pressure ulcers; Poz'in'Form cushions excel in bedsores prevention
- Micropore fabric combined with bi-elastic, breathable Pharmatex or Lenzing FR cover for maximum breathability
- Thin microfiber filling fits body shape, optimizes pressure distribution, and reduces friction and shear
- Designed for support and comfort of long-term bedridden individuals
- Microfibers prevent creep effects and early deterioration
- Available in various forms, sizes, and models for diverse positioning needs



Code	Description
PED470078	Poz'in'Form' (as per script)
NDIS 05_044839811_0103_1_2	

DO Ability Nino Series Bed



- The Nino Series allows children to perceive the world around them while at the same time providing protection and security
- The transparent Perspex merges with the room allowing occupants to have an unrestricted view of their immediate surroundings yet still protects against falls and injuries
- Also available in different powdercoat colours, or with higher door elements or doors on both sides
 - Nino series - single width
 - Electric hi-low
 - Electric backrest
 - Electric kneebreak
 - Standard locking castors
 - Gates one side with Perspex infills
 - Fixed panel one side with Perspex infills
 - 750mm high sides (from mattress base)
 - Head and foot infills in natural oak matt finish
 - Battery backup



Code	Description	Length	Width	Sleeping Surface	Height Range	Back Angle	Knee Angle
PED491830	Single	2125	1025	2000 x 900	390 - 830	77°	16°
PED491700	King Single	2125	1195	2000 x 1070	390 - 830	77°	16°

NDIS 05_181210811_0103_1_2

Bathroom



Firefly GottaGo Portable Toilet Seat



- GottaGo is an innovative, compact, portable toilet seat which provides you with the consistency needed for toileting success
- Complete with a neat backpack for stress-free travel and a toilet training guide, making toilet time easier for everyone
- Retractable legs for easy folding, stowing and different height settings
- Can be used over any toilet or potty
- Snug hug-like fit and warm touch materials for extra comfort
- Accessories include the footplate, straps, seat strap, potty and backpack
- Lightweight and sturdy

2 YEARS WARRANTY

SWL 30KG

Code	Description	Seat Width	Seat Depth*	Back Height
PED491155	Size 1 with Seat - Frame & all Accessories	225	220	240
PED491160	Size 2 with Seat - Frame & all Accessories	270	290	360
PED491165	Toilet Seat Accessory - GottaGo Size 1			
PED491170	Toilet Seat Accessory - GottaGo Size 2			

NDIS 05_091200111_0103_1_2

*Minimum Seat Depth

Firefly Splashy Portable Bath Seat

- Flexible use in various environments
- 26 recline options: 106° to 140°, compatible with GoTo Seat
- Adaptable Splashy harness: 5-point or 3-point pelvic support
- Warm-touch, quick-dry foam covering, soft wet or dry
- Sculpted shape for comfortable, supportive sitting
- Adjustable harness and bumper supports with toggles
- 4 detachable bumper supports for headrest,ommel, or anti-slide shelf
- Suitable for ages 1-8 years



Optional Firefly Backpack



SWL 30KG

2 YEARS WARRANTY

GREEN

PINK

BLUE

Code	Description	Seat Back Rest Height x Width	Seat Base Depth x Width	Floorsitter Height x Width
PED491000*	Blue Base with Green Harness	610 x 385	300 x 350	630 x 420
PED491005	Blue Base with Blue Harness	610 x 385	300 x 350	630 x 420
PED491010	Blue Base with Pink Harness	610 x 385	300 x 350	630 x 420
PED491015	Green Base with Pink Harness	610 x 385	300 x 350	630 x 420
PED491020	Green Base with Blue Harness	610 x 385	300 x 350	630 x 420
PED491025	Green Base with Green Harness	610 x 385	300 x 350	630 x 420
PED491030	Pink Base with Green Harness	610 x 385	300 x 350	630 x 420
PED491035	Pink Base with Pink Harness	610 x 385	300 x 350	630 x 420
PED491050*	Pink Base with Blue Harness	610 x 385	300 x 350	630 x 420
PED491250	Firefly Backpack for use with GoTo & Splashy			

NDIS 05_122500821_0105_1_2 *05_093305121_0103_1_2

Firefly Splashy Big Portable Bath Seat

- Felxible use in an environment
- Choose from 26 recline options, from upright at 106° to a fully reclined at 140°
- The Splashy Big harness adapts to your child's postural support needs, with the option of a 5-point or the 3-point pelvic support
- 4 x detachable 'bumper' supports for use anywhere on the seat where support is required (i.e. head rest,ommel or anti-slide shelf)
- The new Splashy Big stabiliser feature provides additional stability when used outside of the bath
- Age range: 8-16 years



Optional Firefly Backpack



SWL 55KG

2 YEARS WARRANTY

BLUE

LIME

CORAL

Code	Description	Seat Back Rest Height x Width	Seat Base Depth x Width	Floorsitter Height x Width
PED491100*	Blue Base with Coral Harness	720 x 450	400 x 450	795 x 420
PED491150*	Blue Base with Lime Harness	720 x 450	400 x 450	795 x 420
PED491250	Firefly Backpack for use with GoTo & Splashy			


NDIS 05_122500821_0105_1_2 *05_093305121_0103_1_2

Leckey Advanced Bath Chair



SWL 72KG
3 YEARS WARRANTY


- Supports and holds the child to assist with bathing
- Easy-to-use adjustable legs can be set up prior to use
- Simple lever can be used to adjust the backrest recline position in the bath
- Folds flat for easy storage and transportation
- Washable, mesh fabric is breathable and soft to touch



Optional Head Support



Optional Hip Belt



Optional Advanced Chest Belt



Optional Abduction Straps

Code	Description	User Height	Seat Depth	Back Height
PED478480	Size 1 1 - 5 years (as per script)	750 - 1050	225	350
PED478500	Size 2 4 - 9 years (as per script)	950 - 1250	275	400
PED478510	Size 3 8 - 14 years (as per script)	1150 - 1450	335	460
PED478530	Size 4 12 - 18 years (as per script)	1350 - 1650	400	550

NDIS 05_093305121_0103_1_2

Code	Description
PED478572	Head Support
PED478574	Hip Belt
PED478576	Advanced Bath Chest Belt
PED478580	Abduction Belt Size 1 & 2
PED478578	Abduction Belt Size 3 & 4

NDIS 05_093305121_0103_1_2

Medifab Sharky Bathroom Chair



SWL 60KG
2 YEARS WARRANTY

- Versatile, practical, and functional, these aspects have been the primary focus throughout the development
- Backrest recline and 45° tilt-in-space
- Hygiene opening - removable section in rear area to assist with hygiene
- Height adjustable
- Compact folding - folds down to a height of 28cm for easy transport and storage
- Removable padding for easy care and cleaning

Code	Description	Seat Width	Seat Depth	Back Height	Seat Angle	Backrest Angle	Knee Angle	Footplate Angle	Wheel
BTS105970	Size 1	240, 270 & 300	270 - 360	420	-10° - 45°	72° - 120°	78° - 138°	80° - 100°	750 & 1000
GTS105972	Size 2	290, 320 & 350	270 - 400	420	-10° - 45°	72° - 120°	78° - 138°	80° - 100°	750 & 1000

NDIS 05_093305121_0103_1_2

Otter Bathing System



2 YEARS WARRANTY

- Durable, lightweight plastic frame
 - Easy-to-clean surfaces and machine washable, removable fabric
 - Allows bathing a child at the bottom of the bath with back support
 - Adjustable, slip-resistant legs and height
 - Leg straps for abduction and adduction
- Two positioning straps for enhanced safety
 - Adjustable head cushions for comfort and support
 - Folds flat for storage and transport
 - Maintains upper body balance
 - Available in three sizes with optional accessories

Code	Description	Back Height	Seat Width	Seat Depth	SWL
PED480829	Small	450	370	330	34
PED480830	Medium	640	370	330	75
PED470072	Otter Tub Stand				91

NDIS 05_093303078_0103_1_2

Apex Mobility Rifton HTS



1 YEAR WARRANTY

- The unique design of the optional seat pad opening - open to the rear for larger clients, facilitates clear access for cleaning and hygiene
- The portability kit provides clients with greater independence - and their families greater freedom to travel
- Optional molded foam pads for the back and seat provide comfort for clients and easy disinfecting for caregivers
- Optional gas-assisted tilt-in-space enables the HTS to tilt 5° forward for easier transfer and better toileting position, as well as 25° back for showering and hair washing
- On, Over, Off, no tools required, no matter what your toileting situation or special need, the HTS can meet it with its unmatched versatility

Code	Description	Seat Depth	Seat Width	Backrest Height	Backrest Angle*	User Height	SWL
BTS105976	Small	203 - 280	305	419	6° 6° & 18°	300 - 460	34kg
BTS105978	Medium	292 - 268	356	495	6° 6° & 18°	420 - 560	68kg
BTS105980	Large	368 - 445	432	622	6° 6° & 18°	500 - 740	113kg

NDIS 05_093305121_0103_1_2

*Forward | Backward

Apex Mobility Rifton The Wave



1 YEAR WARRANTY

- Extra wide seat with open sides
- Large size accommodates adults up to 6'2" and 250 lb.
- Completely configurable; all sizes of the Rifton Wave fit on any of the 3 base options
- Calf rest adjusts independently of seat angle
- Tub transfer base enables clients to be placed in the bath chair using a mechanical lift and then rolled back and rotated over the tub for showering
- Conversion kit allows user to move up from small to medium size by just purchasing a new back, seat and cover
- Least costly alternative to remodeling a bathroom or assisted living, without resorting to bed baths

Code	Description	Seat Depth	Seat Width	Back Height	Back & Seat Angle	User Height	SWL
PED491800	Small	267	381	584	0° - 75°	760 - 117	34kg
PED491810	Medium	330	381	673	0° - 75°	1020 - 1420	68kg
PED491820	Large	432	445	813	0° - 75°	1270 - 1880	113kg

NDIS 05_044839811_0103_1_2

RAZ - CAT Mobile Shower Commode

- Designed for pediatrics, adolescents, or smaller adults
- Compact design for a better fit for smaller clients
- 18" width base and seat frames for toilet positioning
- Reduced seat frame depth with depth-adjustable back (12" to 15" or up to 17" with overhang)
- Custom seats available
- Standard seat with elongated aperture for 12" to 15" depths
- Shortened base frame by 3" for better manoeuvrability
- Depth-adjustable back for growth capability



SWL 115KG

Code	Description	Width	Length	Seat Depth	Back Height
BTC091800	One Size (as per script)	550 - 460	810	280 - 400	380 / 460

NDIS 05_091203821_0103_1_2

1 YEAR WARRANTY Components

2 YEARS WARRANTY Frame

