





#### **Powerful Suspension**

Provides a smoother ride, even over rough terrain, to help you handle the adventures ahead.



### **Highly Configurable**

Allows you to maintain proper positioning as you grow.



#### **Multi-Motion Armrest**

Perfectly positioned armrest with unique rotation and flip back movements that allow for easier transfers.

#### **Colour Options**



Volcano Red

Purple

Rain



Pop Star Pink

Phantom

Gray



Cobalt Blue



Onyx Black

## **Additional Hub Colour Options**



Chili Pepper Red



Titan Blue



Extreme Yellow



Vintage Wine



Crystal Blue

# Specifications

Maximum User Weight	75 kg
Maximum Speed	8 km/h or 10 km/h
Range	25 km
Base Width	24"
Base Length	40"
Minimum Turning Radius	25"
Ground Clearance	3"
Recommended Battery Type	Group 34
Weight Including Batteries	123.3 kg
Drive Electronics	90A R-Net
Maximum Slope	10° up, 6° down†
Seat-to-Floor Height	16.5" (elevator)
	16.5". 17.25". 18" (seat tube)

Optional Seat Elevator	8"
Tilt Options	0° – 45° (power) 0° – 45° (manual)
Manual Recline	85° – 120°
Legrest Elevation	70° (swing-away) 90° – 150° (center mount elevating)
Backrest Height	16.5" – 22.5"
Seat Depth	10" – 18"
Seat Width	11" – 16"
Armrest Pad Lengths	10", 13"

<sup>\*</sup>Actual driving range will vary based on driving and battery conditions.

<sup>&</sup>lt;sup>†</sup>Refer to Owner Manual for safe driving practices.