





Miniflex Power Seating

Unique power seat elevation and tilt designed to provide valuable positioning and exceptional low seat to floor height.



Parental Approved

Optional features such as front bumper and remote stop give any parent that extra piece of mind for their little explorer.

Colour Options



Volcano Red



Cobalt Blue



Popstar Pink



Purple Rain



Extreme Yellow



Platinum Silver



Green Apple



Onyx Black

Additional Hub Colour Options



Chili Pepper Red



Vintage Wine



Titan Blue



Crystal Blue



Phantom Gray



Sublime Green

Specifications

Maximum User Weight	35 kg
Maximum Speed	6.5 km/h
Range	24 km*
Base Width	23"
Base Length	32"
Minimum Turning Radius	21"
Ground Clearance	2.5"
Recommended Battery Type	M40
Weight Including Batteries	102.5 kg
Drive Electronics	R-Net
Maximum Slope	10° up, 6° down†
Seat-to-Floor Height	15"

Optional Seat Elevator	8"
Power Tilt	0° – 25°
Manual Recline	85° – 120°
Manual Legrest Elevation	90° – 150°
Backrest Height	14", 16", 20"
Seat Depth	10", 12", 14"
Seat Width	10", 12", 14"
Armrest Pad Lengths	8", 10"

^{*}Actual driving range will vary based on driving and battery conditions.

†Refer to Owner Manual for safe driving practices.