





Catlink Suspension

Provides a smoother ride, even over rough terrain, to help you handle the adventures ahead.



Highly Configurable

Allows you to maintain proper positioning as you grow.



Multi-motion Armrest

Perfectly positioned armrest with unique rotation and flip back movements that allow for easier transfers.

Shroud Colour Options





Pop Star Pink



Sublime Green



Cobalt Blue

Chili Pepper Red

Additional Hub Colour Options

Titan Blue



Extreme Yellow















Onyx Black

Vintage Wine



Crystal Blue

Specifications

Maximum User Weight	75 kg
Maximum Speed	10 km/h
Range	20 km*
Base Width	24"
Base Length	36.5"
Minimum Turning Radius	20"
Ground Clearance	3"
Recommended Battery Type	Group 24, Group 34
Weight Including Batteries	142 kg G34 / 151 kg G24
Drive Electronics	120A R-Net
Maximum Slope	9° up, 9° down†
Seat-to-Floor Height	17" (elevator)
	17". 17.75". 18.5" (seat tube)

•	0° – 45° (power) 0° – 45° (manual)
Manual Recline	85° – 120°
•	70° (swing-away) 90° – 150° (centre mount elevating)
Backrest Height	16.5" – 22.5"
Seat Depth	10" – 18"
Seat Width	11" – 16"
Armrest Pad Lengths	10", 13"

^{*}Actual driving range will vary based on driving and battery conditions. †Refer to Owner Manual for safe driving practices.