



TiFit. Personal Ergonomics. +

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. A TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do. TiFit...not standard.

Custom Configuration. +

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimise your riding experience. This includes choices of wheels, handrims, tyres, wheel locks, colour, tattoos, and more. TiLite offers more options for the individual.

Technical Specifications ⊢

Frame Features	Aluminium Dual-Tube Rigid TiFit Full Adjustability
Total Chair Width*	23"
Seat Width	12" – 20"
Seat Depth	12" – 20"
Front Seat Height	15" – 21"
Rear Seat Height	13" – 21"
Front Frame Angle	70°, 75°, 80°, 85°, or 90°
Seat Back Type	Folding Seat Back Adj. Height 8.5" – 21" Adj. Angle 80° – 101°
Camber	0°, 2°, 4°, 6°, 8°, or 12°

Centre of Gravity Adjustment	5"
Front Seat Width (Seat Taper)	0", 1", 2", 3"
Transit Options	Transport Bracket
Transport Weight*	5.4 kg.
Weight Limit	120.2 kg Standard 158.7 kg Heavy Duty
HCPCS Code	K0005
Warranty	Five Years on Frame One Year on TiLite Components

^{*} Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.27 kg.

Colour Options

High Gloss Paint











Super Wet Black



Titanium







Azure



Pearlescent Paint

White

Midnight Blue Pearl Pearl

Metallic Paint



Salsa Red Tangerine Acid Green Ocean Blue Metallic Metallic



Metallic





Pepper













Metallic



Anodised Colours

Orange

Sapphire Ultraviolet Fuchsia

permobil