## permobil



# Multifaceted simplicity.

The Aero X features an innovative modular design with clean modern lines and configurability. Featuring multiple configurations that include heavyduty, hemi, and amputee options. Add superior rolling dynamics and full adjustability, and you get an Aero X that excels in both performance and price.



#### **Technical Specifications**

Frame Features	Aluminium Box Frame Folding Modular Design Full Adjustability
Total Width Folded*	10.5"
Total Chair Width*	22.5"
Seat Width	12" – 22"
Seat Depth	12" – 20"
Front Seat Height	13.5" – 21.5"
Rear Seat Height	13" – 20"
Front Frame Angle	Fixed Front - 70°, 80°, or 85° Swing Away - 70°, 80°, 85°, or 90°
Seat Back Type	Folding Seat Back Adj. Height 11" – 20.5" Adj. Angle 80° – 101°
	<i>Fixed Seat Back</i> Adj. Height 11" – 20"

Camber	0°, 2°, 4°, or 6°
Centre of Gravity Adjustment	6"
Transit Options	Transit Tie-Down Transport Bracket
Transport Weight*	Fixed Front - 7 kg Swing Away - 8.2 kg
Weight Limit	120.2 kg Standard 158.7 kg Heavy Duty
HCPCS Code	K0005
Warranty	Five Years on Frame One Year on TiLite Components

\* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.27 kg.

#### Modular Design

Optimise your ride with custom configuration. This chair's modular design allows you to add options including after market seating systems, switching fixed and swing away front ends, and growing chair width and / or depth to meet changing user needs.

#### Custom Configuration +

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimise your riding experience. This includes choices of wheels, handrims, tyres, wheel locks, colour, tattoos, and more. TiLite offers more options for the individual.

#### **Colour Options**



#### Pearlescent Pai



### permobil

Ph: 1300 845 483 | orders.au@permobil.com | www.permobil.com/en-au