



TiFit. Personal Ergonomics. -

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. A TiFit chair brings unparallelled performance and function. You will feel the difference in everything you do. TiFit...not standard.

Custom Configuration. -

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimise your riding experience. This includes choices of wheels, handrims, tyres, wheel locks, colour, tattoos, and more. TiLite offers more options for the individual.

Technical Specifications

Frame Features	Aluminium Mono-Tube Rigid TiFit Full Adjustability
Total Chair Width*	23"
Seat Width	10" – 20"
Seat Depth	10" – 20"
Front Seat Height	16" – 21"
Rear Seat Height	12.5" – 20"
Front Frame Angle	70°, 75°, 80°, 85°, or 90°
Seat Back Type	Folding Seat Back Adj. Height 8.5" – 21" Adj. Angle 80° – 101°

Camber	0°, 2°, 4°, 6°, 8° or 12°
Centre of Gravity Adjustment	6"
Front Seat Width (Seat Taper)	0", 1.5", 2.5"
Transit Options	Transport Bracket Transit Tie-Down
Transport Weight*	5.6 kg
Weight Limit	120.2 kg
HCPCS Code	K0005
Warranty	Five Years on Frame One Year on TiLite Components

^{*} Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.27 kg.

Colour Options

High Gloss Paint

















Blue Ice





Lava







Blue



Plum

Pearlescent Paint

Blue Pearl

Midnight White

Pearl

Metallic Paint



Metallic





Pink









Grey













Metallic

Red Pepper

Rurnt

Orange

Sapphire Ultraviolet Fuchsia

Black

