




Designed to move.

The clean frame design features a seat pan and side panels that provides easy access for adjustments, making changes simple and intuitive. Unlike on other chairs, all of these adjustments can be done with



Frame designed for maximum stability



Specifications

Frame features	Carbon frame with 1" lower aluminium tubing
Seat width	8" - 14" (20 - 35 cm)
Seat depth	9" – 15" (23 - 38 cm)
Front seat height	12" - 16.5" (30 - 42 cm)
Rear seat height	12" - 16" (30 - 41 cm)
Front frame angle	Unlimited
Seat back height	11" to 13" (8 - 33 cm) with standard back

Centre of gravity adjustment	0.5" – 3.5" (1 - 9 cm)
Transit options	Transit Tie-Down, Transport
Transport weight	4.9 kg*
Weight limit	43 kg transit, 54 kg standard
Seat to footrest	3" – 14.5" (8 - 37 cm) without high mount parts

Frame Colours



Acid Green Metallic



Neon Pink



Plum



Blue Metallic

Anodised Colours









Red Pepper

Ultraviolet



[†]Specifications may vary based upon configuration.

* Weight varies by configuration and does not include rear wheels.