



TiFit. Personal Ergonomics.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. A TiFit chair brings unparallelled performance and function. You will feel the difference in everything you do. TiFit... not standard.

Custom Configuration. +

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimise your riding experience. This includes choices of wheels, handrims, tyres, wheel locks, colour, tattoos, and more. TiLite offers more options for the individual.

Technical Specifications

Frame Features	Titanium Dual-Tube Rigid TiFit Limited Adjustability
Total Chair Width*	22"
Seat Width	12" - 20"
Seat Depth	12" – 20"
Front Seat Height	15" – 21"
Rear Seat Height	12.5" – 21"
Front Frame Angle	60°, 65°, 70°, 75°, 80°, 85°, or 90°
Seat Back Type	Folding Seat Back Adj. Height 8.5" – 21" Adj. Angle 80° – 101°

Fixed Seat Back Fixed Height 5" - 20" Fixed Angle 80° - 101°

Camber	0°, 2°, 4°, 6°, 8° or 12°
Centre of Gravity Adjustment	5"
Front Seat Width (Seat Taper)	0", 1", 2", 3"
Ergonomic Seat	4"-10"
Transit Options	Transport Bracket
Transport Weight*	4.2 kg
Weight Limit	120.2 kg
Warranty	Five Years on Frame One Year on TiLite Components

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.27 kg.

Colour Options

High Gloss Paint









Neon

Pink









Black









Blue



Plum

Natural

Satin Titanium

Metallic Paint



Metallic





Acid Green Metallic



Ocean Blue Metallic











Pearlescent Paint



Midnight Blue Pearl

White Pearl

Anodised Colours





Blue Ice



Orange







