



Specifications

Frame features	Aluminium mono-tube rigid Growable Full adjustability
Total chair width*	19"
Seat width	8" – 15" (2" of built-in Seat Width Adj.)
Seat depth	8" – 16" (3" of built-in Seat Depth Adj.)
Front seat height	13" – 20" (up to 2.25" of built-in front seat height adj.)
Rear seat height	12" – 19.5" (up to 3.5" of built-in rear seat height adj.)
Front frame angle	70°, 80°, or 90°
Seat back type	Folding seat back Adj. height 8.5" – 21" Adj. angle 80° – 101°

Camber	0°, 2°, 4°, 6°, 8°, or 12°
Centrer of gravity adjustment	6"
Transport options	Transit tie-down
Transport weight*	5.4 kg
Weight limit	74.8 kg
Warranty	Five years on frame One year on TiLite components

 $^{^{*}}$ Weights and widths based on 12" x 12" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weights. For reference, a pair of standard 22" rear wheels weighs 3.17 kg.

Frame Colours

High Gloss Paint











Blue

Electric Plum

Super Wet Black

Metallic Paint







Matte Textured Paint



Flame Orange Green



Hyper

Pearlescent Paint



Midnight White Blue Pearl Pearl

Matte Paint



Grey Lava

Component Colours

Anodised



Red

Pepper



Orange



Sapphire









Black