permobil



Unrivaled design. Undeniable beauty.

The ZR's minimalist 1.25" monotube frame combined with TiFit craftsmanship, results in a chair that can only be a TiLite. Titanium's advanced vibration dampening and our unique frame design provide the perfect balance between strength and comfort.

Ti



TiFit. Personal Ergonomics. +

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. A TiFit chair brings unparallelled performance and function. You will feel the difference in everything you do. TiFit...not standard.

Custom Configuration. +

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimise your riding experience. This includes choices of wheels, handrims, tyres, wheel locks, colour, tattoos, and more. TiLite offers more options for the individual.

Technical Specifications ⊢

Frame Features	Titanium Mono-Tube Rigid TiFit Limited Adjustability
Total Chair Width*	23"
Seat Width	10" – 20"
Seat Depth	10" – 20"
Front Seat Height	16" – 21"
Rear Seat Height	12.5" – 20"
Front Frame Angle	70°, 75°, 80°, 85°, or 90°
Seat Back Type	Folding Seat Back Adj. Height 8.5" – 21" Adj. Angle 80° – 101°
	Fixed Seat Back

Fixed Height 8" – 20" Fixed Angle 80° – 101°

Camber	0°, 2°, 4°, 6°, 8° or 12°
Centre of Gravity Adjustment	2.25"
Front Seat Width (Seat Taper)	0", 1.5", 2.5"
Ergonomic Seat	4"–10"
Transit Options	Transport Bracket
Transport Weight*	4.2 kg
Weight Limit	120.2 kg Standard
	158.7 kg Heavy Duty
Warranty	Five Years on Frame
	One Year on TiLite Components

* Weights and widths based on 16" x 16" chair with O^o camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.27 kg.



perm_obil