



POWER CHAIRS	ADULT						
	Mid Wheel Drive		Rear Wheel Drive		Front Wheel Drive		Front Wheel Drive / Standing
	M300 Corpus 3G	M300 Corpus HD	C350 Corpus PS	C350 Corpus 3G	F5 Corpus	F3 Corpus	F5 VS
Seat Width (mm)	430 - 580	470 - 620	410 - 560	430 - 580	430 - 580	430 - 580	430 - 580
Seat Depth (mm)	360 - 560	470 - 620	410 - 560	360 - 560	360 - 560	360 - 560	360 - 560
Floor to Seat Ht (mm)	445	460 - 510	465	445	445 - 495	445 - 495	445 - 495
Max User Wt (kg)	136	205	136	136	136	136	136
Based W x L (mm)	610 x 925	660 x 965	610 x 915	610 x 915	650 x 1070	610 x 1020	650 x 1090
Max Speed (kph)	10, 12	8	10	10	10	10	10
Features	Turning radius of 51 cm, PG 120 amp, R-Net electronics, Climbing and Traction Link system, Provides quick access to components	High torque motors; strengthened suspension; Corpus HD seating system; strong actuators; adjustable armrest and legrests	Compact base; adjustable suspension; R-net electronics; height adjustable and flip back armrests; removable leg rests	Compact base; adjustable suspension; R-net electronics; corpus 3G seating system; memory positioning fully adjustable armrests and leg rests	Adjustable suspension, 14" seat elevate, articulated power footplate	Adjustable suspension, 12" seat elevate, power elevating legrests	Adjustable suspension, 14" seat elevate, automatic anti-tippers, articulated power footplate

POWER CHAIRS	PAEDIATRIC			
	Front Wheel Drive	Rear Wheel Drive	Mid Wheel Drive	Front Wheel Drive
	K300 PS Jnr	K450	M300 PS JR	Koala
Seat Width (mm)	280 - 410	250 - 410	280 - 410	250 - 360
Seat Depth (mm)	150 - 460	250 - 460	150 - 460	250 - 360
Floor to Seat Ht (mm)	420	360 - 460	430	375
Max User Wt (kg)	75	60	75	35
Based W x L (mm)	610 x 1020	620 x 840	610 x 925	585 x 810
Max Speed (kph)	10	8	10	6.5
Features	Compact base; adjustable suspension; R-net electronics; fully growable seating system; height and angle adjustable, flip back armrests; removable leg rests	Unique 'terra' function allows user access to the ground, Fully growable seating system, Advanced suspension system	Turning radius of 51 cm, PG 120 amp. R-Net electronics, Climbing and Traction Link system, Provides quick access to components	True paediatric size power base and seating; fully growable seating systems; tilt and elevate standard

MANUAL CHAIRS	ADULT	
	Standing	
	Helium LSA	Helium LS/LSE
Seat Width (mm)	360 - 460	360 - 480
Seat Depth (mm)	430 - 540	410 - 550
Floor to Seat Ht (mm)	510	510
Weight (kg)	16	25 - 32
Max User Wt (kg)	100	120
Features	Lightweight, dual component knee pads, chest belt & armrests	LS offers manual standing, LSE offers a power function. Includes dual component knee pads, chest belt & armrests