PHOENIX"



The Phoenix[™] is the most supportive and functional of the Seating Matters chairs and is available in both adult and children's sizes.

Designed to give maximum support to the upper body through lateral and head supports, the Phoenix[™] will also adapt to the changing needs of a person being rehabilitated or meet the long term needs of a progressive condition. The unique adjustability allows the Phoenix[™] to be used with individuals of different dimensions and with varying needs which is ideal for use in acute wards or changing environments.

The adjustable lateral, shoulder and head supports can conform to the person's individual posture and position with fast and easy adjustments. With tilt in space and adjustable components, the Phoenix™ will meet a range of individual support needs, sizes, and levels of mobility making it ideal for multiple users. Used in rehabilitation settings, it can assist in the development and maintenance of independence or for those with progressive conditions; it can provide additional support as necessary.

The Phoenix[™] is effective for adults with the following disabilities and conditions:

Cerebral palsy Spina bifida Spinal cord injury Parkinson's disease Motor neurone disease Stroke Multiple sclerosis Alzheimer's disease Advanced dementia Osteoporosis

The Phoenix[™] is effective in seating clients with the following postural and medical conditions:

Pressure ulcers Reduced upper body tone Kyphosis Scoliosis Muscle contractures Oedema Adjustable head support

Adjustable shoulder support

Back height adjustment

Pressure management materials

Arm angle adjustment

Removable seat cushion

Seat width adjustment –

Forward tilt

Elevating leg rest

Footplate height adjustment

Swivel castors



Push handle

Back angle recline

> *Lateral supports

Arm height adjustment

Tilt in space

Seat depth adjustment

Removable arms

Negative angle leg rest

Braking castors

Patient lift/ Hoist accessible

> Directional lock castor

Angle adjustable footplate