PRISM MEDICAL UK

Exclusively Distributed by





Portable P200 Ceiling Hoist

Portable Hoist that can be used on the go

PRISM MEDICAL UK

Exclusive Distribution with High Quality Service & Support

Prism Medical UK are the global leader in innovation and safety focused medical products. This trusted brand has a proven history in Australia and is the preferred choice of ceiling hoist solutions for many Aged Care and Hospital groups. Key safety features and stringent product testing ensure Prism Ceiling Hoists are the top performers in their category.

Aidacare's partnership with Prism Medical UK ensures that extremely high quality and internationally tested products are readily available in Australia and backed with excellent sales, installation and preventative maintenance support nation wide.

P200 Portable Ceiling Hoist

The P200 with its unique innovative design is one of the lightest portable hoists on the market, yet it is capable of lifting up to 200 kg. The P200 has numerous safety features to ensure safer lifting and handling for both the client and caregiver. Suitable for use with fixed ceiling tracks or mobile gantry systems. The design and features of the P200 makes it perfect for use in community or nursing environments.

Quick Release System (QRS)

The Freeway QRS® allows the P200 to be safely and easily attached or removed from the gantry or track system by hand or by means of a new and improved reacher.

Item Code: LIC354170 - QRS Reacher Bar (sold as an optional accessory)



To Attach

Attatch the QRS carabiner into the reacher and pull on the lift tape to open the red latch. Insert the carabiner into the trolley set and remove the reacher – which closes and locks the carabiner.





To Remove

Holding the lift tape, slide the reacher up the tape and insert into the QRS carabiner. Push to reacher upwards to fully open the QRS latch and remove from the trolley set.

Exclusively Distributed by



Key Features



Simple On Board Controls

- Complements the pneumatic hand set
- Quiet operation with `soft start/ stop' for client comfort

On Board Display

- Audio and visual battery indicators warn when the battery needs recharging, ensuring the hoist is always ready to use and avoids the battery running out at a critical moment
- Integrated long life battery offers the convenience of on board charging. Eliminates the risk of batteries and charging units getting lost or damaged (keeping running costs down). Battery fully charged in 4hrs
- 'Auto sleep' facility prolongs battery life



Emergency Lowering Facility

- Allows the safe lowering of the hoist in the unlikely event of the failure of other controls
- Other safety features: Over speed governor, built in safety limits & slack tape sensor
- The slack tape sensor will automatically stop the unit if accidentally lowered onto a resident



Integrated Full Width Carry Bar with Unique Contoured Shape

• Creates optimum space for the client and the caregiver



'Room to Room' Transfer Facility

• Allows the client to be transferred from 'room to room' without the need to be disconnected from the hoist or for costly building modifications



Carry Bar

Folds for ease of transporting and storage



ON/OFF Switch (on Rear of Hoist)

Ergonomically Designed HandsetSimple 2 button operation

Safe to use in wet conditions

• Allows the hoist to be switched off quickly in the case of an emergency



Carry Bar Latch Sections

- Client safety is assured with positive locking latches to secure the sling
- Anodised coating ensures longevity of operation and ease of cleaning

Pneumatically controlled, requiring minimum pressure to operate



- Clip on Rear of Handset
 - Allows the handset to be conveniently positioned on the hoist carry bar when not needed

PRISM MEDICAL UK

Exclusively Distributed by



Hoisting Range



P200 Ceiling Hoist Range

Code	Description	S.W.L
LIC350110	P200 Portable Ceiling Hoist - Manual Traverse & Quick Release System	200kg
LIC350120	P200 Portable Ceiling Hoist - Manual Traverse	200kg
LIC354170	QRS Reacher Bar (sold as an optional accessory)	-

P200 Specifications

Specification			
Lift Motor	24 VDC		
Charger Input	100-240 VAC, 0.35 Amps		
Charger Output	29.4 VDC, 0.5 Amps		
Batteries	24 VDC (2x12 VDC) 2.3 AH		
Hoist Case	Flame Retardant ABS		
Hand Control	Pneumatic		
Hoisting Range	Up to 2070mm		
Hoist Weight	4.8kg		
Maximum Load	200kg		
Internal Fuses	15 Amp		
Rated Performance	15 - 20 lifts at 200kg from fully charged		

Optional Accessories





Clax Transport Cart WAT758020

Free Standing Gantry LIC352100

Quality & Compliance

- All systems are tested and installed in compliance with Australian Standard AS/NZS ISO 10535-2011: Hoists for the transfer of disabled persons
- Aidacare is an ISO9001 compliant company
- Aidacare uses an accredited network of trained installers, engineers and load testers.

