



<15 km
Range per Charge



120 kg
Weight Capacity



145 cm
Turning Radius



GO GO *My-Lite*

Ultra lightweight, portable.

Experience ultimate portability with the Go Go Mg Lite mobility scooter, designed for travel, convenience, and independence. Reaching speeds of up to 6.4 km/h, it features an ultra-light magnesium foldable frame weighing just 17.9 kg.





Easy disassemble



Battery level display



Digital display with contactless lock

RANGE PER CHARGE (MAX) ¹	Up to 15 km
MAXIMUM WEIGHT CAPACITY	120 kg
TURNING RADIUS	145 cm
MAXIMUM SPEED ¹	6.4 km/h
LENGTH & WIDTH	L: 99 cm W: 62.5 cm
SEAT TYPE WIDTH & DEPTH (STANDARD)	Compact foldable W: 42 cm (16.5") D: 34 cm (13")
SEAT-TO-DECK HEIGHT	44.1 cm
SEAT-TO-GROUND HEIGHT	56 cm
GROUND CLEARANCE ³	7.5 cm at motor

OVERALL WEIGHT (NOT INC BATTERY) ²	16.3 kg
WEIGHT OF HEAVIEST PIECE	12 kg (frame section)
SUSPENSION	N/A
3W/4W OPTION	4W
FRONT TYRES	7.5" Solid
REAR TYRES	7.5" Solid
BATTERY REQUIREMENTS	12V - 11Ah Lithium (264Wh)
BATTERY CHARGER	Off-Board, 2A
BATTERY WEIGHT	1.6 kg each
PRODUCT CODE	■ Black - SCT670617

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Includes standard seat weight.
3. Does not account for anti-tip clearance

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.

Exclusively distributed by

aidacareTM
HEALTHCARE EQUIPMENT

1300 133 120 | aidacare.com.au

PRIDE[®]
MOBILITY

