

**PRIDE**<sup>®</sup>  
MOBILITY



**<15 km**  
Range per Charge



**136 kg**  
Weight Capacity



**85 cm**  
Turning Radius



**Jazzy**  
CARBON MAX

## Exceptional strength lightweight frame

The Jazzy Carbon Max, a revolutionary powerchair that combines exceptional strength with ultra-lightweight design. Engineered for performance and comfort, it supports users up to 136 kg while maintaining an incredibly light frame of just 11.5 kg.





Easy folding mechanism



Folding footplate



Lightweight strong carbon fibre frame

RANGE PER CHARGE (MAX) <sup>1</sup>	Up to 15 km
MAXIMUM WEIGHT CAPACITY <sup>4</sup>	136 kg
TURNING RADIUS	85 cm
MAXIMUM SPEED <sup>1</sup>	6 km/h
LENGTH & WIDTH	L: 89 cm W: 59 cm
WEIGHT OF HEAVIEST PIECE	11.5 kg
GROUND CLEARANCE	2.5 cm at motor
SEAT SIZE WIDTH & DEPTH	W: 44 cm x D: 42 cm (17" x 16")
HIGH BACK SEAT-TO-FLOOR RANGE	48.5 cm fixed

STANDARD ELECTRONICS	35 A
DRIVE TRAIN	Two motor, rear-wheel drive
BRAKING SYSTEM	Regenerative and electro-mechanical
SUSPENSION TYPE	Limited front only
FRONT CASTORS	6.89" solid
DRIVE WHEELS	8.25" solid
BATTERY REQUIREMENTS	(1) 12Ah, 24V (lithium-ion, 288Wh)
BATTERY CHARGER	Off-board. 2A lithium-ion
BATTERY WEIGHT	2.1 kg each
COLOUR AND CODE	■ Blue - PWS654160

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Width of base (not including seat).
3. Deep cycle (AGM or gel cell type recommended)
4. Does not account for anti-tip clearance.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.

Exclusively distributed by

**aidacare**<sup>TM</sup>  
HEALTHCARE EQUIPMENT

1300 133 120 | aidacare.com.au

**PRIDE**<sup>®</sup>  
MOBILITY

Certified System  
Quality  
ISO 9001  
SAI GLOBAL