

**PRIDE**<sup>®</sup>  
MOBILITY



**<20 km**  
Range per Charge



**182 kg**  
Weight Capacity



**92.5 cm**  
Turning Radius



**Jazzy**  
CARBON MAX HD

## Heavyduty Powerchair lightweight frame

The Jazzy Carbon Ultra Max HD, a revolutionary heavy-duty powerchair with additional seat width of 20". With all the unique features of the Jazzy Carbon HD. Exceptional strength with ultra-lightweight design. Engineered for performance and comfort, it supports users up to 182 kg while maintaining an incredibly light frame of just 18.1kg





Full front suspension



Quick disconnect joystick with LCD display



Puddle light

RANGE PER CHARGE (MAX) <sup>1</sup>	Up to 20 km
MAXIMUM WEIGHT CAPACITY <sup>4</sup>	182 kg
TURNING RADIUS	92.5 cm
MAXIMUM SPEED <sup>1</sup>	6 km/h
LENGTH & WIDTH	L: 93.5 cm W: 67 cm
WEIGHT OF HEAVIEST PIECE	14.6 kg
GROUND CLEARANCE	2.5 cm (at motor)
SEAT SIZE WIDTH & DEPTH	W: 50 cm x D: 50 cm (20" x 20")
HIGH BACK SEAT-TO-FLOOR RANGE	48cm fixed

STANDARD ELECTRONICS	35 A
DRIVE TRAIN	Two motor, rear-wheel drive
BRAKING SYSTEM	Regenerative and electro-mechanical
SUSPENSION TYPE	front
FRONT CASTORS	6.5" solid castors
DRIVE WHEELS	10.6" solid castors
BATTERY REQUIREMENTS	(1) 20Ah, 24V (lithium-ion, 480Wh)
BATTERY CHARGER	Off-board, 2A lithium-ion
BATTERY WEIGHT	3.4 kg each
COLOUR AND CODE	■ Blue - PWS654165

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Width of base (not including seat).
3. Deep cycle (AGM or gel cell type recommended)
4. Maximum weight capacity is reduced to 181 kg on 18" wide seat configuration.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.

Exclusively distributed by

**aidacare**<sup>TM</sup>  
HEALTHCARE EQUIPMENT

1300 133 120 | aidacare.com.au

**PRIDE**<sup>®</sup>  
MOBILITY

