



15 km Range per Charge*



Weight Capacity

120 cm Turning Radius



Effortless, Lightweight and Foldable Travel Companion.

With carbon fibre technology, the Litestream® Carbon Power Wheelchair is a great travel companion. In just a few quick steps, the Litestream® Carbon folds down, making it a super lightweight effortless option to take when travelling in planes, trains and cars. It's a convenient mobility solution with the ability to fold and go.

Visit aidacare.com.au for more information.







Conveniently folds down for transportation.



Lightweight carbon fibre frame with under seat storage.



Joystick controller.

Code	PWS654185		
Weight Capacity	136 kg (300 lbs.)	Ground Clearance ^{2,7,8}	5.5 cm (at the motor)
Maximum Speed ⁵	6 km/h (3.7 mp/h)	Seat-to-ground Height	47.5cm (18.7") (without cushion)
Turning Radius ²	120 cm (47.3″)	Seat Type (Standard)	Compact, foldable, mesh L 43 cm x W 38 cm (17"x 15")
Range per Charge ^{1,3}	15 km (9.3 mp)	Weight (excl. battery) ²	16 kg (35.3 lbs.)
Overall Length ²	93 cm (36.6")	Battery Weight ⁹	1.8 kg (4 lbs.)
Overall Width ^{2,6}	57.7 cm (22.7")	Battery Charger	24v 2amp lithium charger
Front Tyres	17.8 cm (7")	Drive System	24V/250W with brushless motor (aluminium alloy material)
Back Tyres	21.5 cm (8.5")	Battery Requirements	24V12Ah with lithium battery (288wh)
Max. Climbing Ability	6°	Colour/s	Carbon

- Varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of (+/- 10%).
- 2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of (+ or 3%).
- Tested in accordance with ANSI/RESNA, WC Vol 2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battern specifications and drive system performance. Test conducted at maximum weight capacity.
- 4. Battery weight may vary based on manufacturer.
- Varies with user weight, terrain type, battery charge, battery condition and tyre condition. This specification can be subject to a variance of (+ or -) 10%.
 Width of base (not including seat).
- Width of base (not including seat).
 Does not account for anti-tip clearance.
- Boes not account for anti-tip cl
 At midframe.
- 9. Battery weight may vary based on manufacturer.



Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Litestream Carbon® and its lithium ion battery are permitted.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range is normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.



