



<25 km
Range per Charge



181 kg
Weight Capacity



113 cm
Turning Radius



The Platinum Difference

Victory® Platinum 4-wheel scooter offers a premium mobility experience. Patented EZ Turn technology helps you navigate tight spaces with 3-wheel precision. Exclusive Charge360 technology keeps you powered up with magnetic charging. And with full front and rear suspension, take on uneven terrain with confidence.





Interchangeable colour shrouds available in Red or Blue



360 magnetic charger



EZ turn technology

RANGE PER CHARGE (MAX) ¹	Up to 40 km
MAXIMUM WEIGHT CAPACITY	181 kg
TURNING RADIUS	113 cm
MAXIMUM SPEED ¹	10 km/h
LENGTH & WIDTH ³	L: 118 cm W: 55.9 cm
SEAT TYPE WIDTH & DEPTH (STANDARD)	Pillow-top, high-back seat W: 46 cm (18") D: 50.8 cm (20")
SEAT-TO-DECK HEIGHT	44.7cm - 48.7cm
SEAT-TO-GROUND HEIGHT	64.7 cm - 69.7.8 cm
GROUND CLEARANCE ⁵	7 cm at motor

OVERALL WEIGHT (NOT INC BATTERY) ²	65 kg
WEIGHT OF HEAVIEST PIECE	30.4 kg (front section)
SUSPENSION	front & rear
3W/4W OPTION	4W
FRONT TYRES	10.4" Solid
REAR TYRES	10.75" Solid
BATTERY REQUIREMENTS ⁴	(2) 12V - 50Ah
BATTERY CHARGER	Off-Board, 5A (magnetic)
BATTERY WEIGHT	11.6 kg each
PRODUCT CODE	SCT670645

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Includes standard seat weight.
3. Width of base (not including seat).
4. Deep cycle (AGM or gel cell type recommended)
5. Does not account for anti-tip clearance

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.

Exclusively distributed by

aidacareTM
HEALTHCARE EQUIPMENT

1300 133 120 | aidacare.com.au

PRIDE[®]
MOBILITY

