## Enjoy Comfort with Advanced Suspension





Full advanced suspension

Be free to go where you want to go! The **GT** offers full advanced suspension for a smoother, more comfortable ride over varied terrain. Front and rear LED lighting illuminates your path and the adjustable delta tiller allows for easy, one-handed operation. Experience the GT feeling for yourself!



Visit pridemobility.com.au to learn more



## Specifications

SCT670670

Weight Capacity	159 kg (350.5 lbs.)
Turning Radius	120 cm (47.2")
Maximum Speed	Up to 10 km/h (6.2 mph)
Overall Length	120 cm (47.2")
Overall Width	60.5 cm (23.8")
Overall Height	90 cm (35.4")
Total Weight	82 kg (108.7 lbs.)
Maximum Range	Up to 30 km (18.6 miles)
Max. Degree of Climb	8°
Motor	24V, 400W
Battery Capacity	2 x 35 AH
Battery Weight	21.3 kg (46.9 lbs.) (battery pack)
Charger	3.5A, 24V
Wheel Size	Front: 22.8 cm (9") Rear: 25.4 (10") (pneumatics)
Ground Clearance	8 cm (3.1")
Colour	

## Features

- Full advanced suspension for optimum handling
- Full LED lighting, including directional signals and headlight
- Easy-drive tiller with wraparound handles
- Convenient front basket
- Easy-access battery compartment



International tiller shown

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.



