OUTSTANDING OUTDOOR PERFORMANCE WITH LUXURIOUS FEATURES

PATHRIDER® 130 XL PURSUIT

The **Pathrider**[®] **130 XL Pursuit** features a powerful drivetrain and large 33 cm (13") pneumatic wheels to handle rugged outdoor terrain. Luxury touches like a Victory tiller, low-profile non-marking tyres, seat post suspension for a smoother ride and an LED kerb light blend seamlessly giving the Pathrider 130 XL Pursuit a sleek, sporty style.



Pathrider [®] 13	O XI Pursuit	Features	
Weight capacity:	158.75 kg (350 lbs.)	Large 33 cm (13") pneumatic tyres for excellent outdoor performance	
Maximum speed ¹ :	Up to 10 km/h (6.2 mph)	 Full lighting package including directional signals and angle- adjustable, dual bulb lower headlight 	
Ground clearance ² :	11.4 cm (4.5")		
Turning radius ² :	167.6 cm (66")		
Overall length ² :	135.9 cm (53.5")		
Overall width ² :	63.5 cm (25")	 Seat post suspension for a smoother ride Victory tiller Front basket standard Easily accessible tie down points (for transport of unoccupied scooter) Rearview mirror standard Options 	
Front tyres:	10 cm x 33 cm (4" x 13") pneumatic		
-			
Rear tyres:	10 cm x 33 cm (4" x 13") pneumatic		
Seat-to-deck:	46 cm (18")		
Seat-to-ground:	63 cm - 67 cm (25" - 26")		
Suspension:	Front & rear		
Range per charge ^{1,4} :	Up to 24.99 km (15.53 miles)		
Maximum obstacle climbing ability:	12 cm (4.72")		
Total weight without batteries ^{2,5} :	101.1 kg (223 lbs.)		
Heaviest piece when disassembled:	42.9 kg (94.5 lbs.) (front section)		
Standard seating:	Type: CRS Deluxe	Rear basket	
	Weight: 22.5 kg (49.5 lbs.)	Safety flag	
	Material: Black Vinyl		
	Dimensions: width - 45 cm (18") (overall width) depth - 45 cm (18") (effective depth) height - 47.6 cm (18.75")	Walker holder	
Drive system:	Rear-wheel drive, sealed transaxle, 24-volt DC motor		
Dual braking system:	Regenerative and electromechanical		
Battery requirements ³ :	Two 12-volt, deep cycle		
battery requirements.	Size: 55 AH (standard)		
	Weight ⁶ : 55 AH: 16.6 kg (36.5 lbs.) each		
Battery charger:	5-amp, off-board		
Warranty:	2-year limited		
Colours:	¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition	and tyre condition. This specification can be	
	subject to a variance of (+ or -) 10%. ² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).		
	³ AGM or gel cell type required.		
Candy Re	Apple theoretical calculations based on battery specifications and drive system performance test co	⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum wight canacity.	
	⁵ Includes standard seat.		
	⁶ Battery weight may vary based on manufacturer.		
	NOTE: This product conforms to all applicable ANSI/RESNA testing requirements and ISO 7176 s	eries EN12184 standards. All specifications subject to change without not	
	The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.		
	Theoretical range is calculated under test conditions in accordance with E in normal use will depend on many factors, including the condition of the of the driver, correct tyre pressures, ambient temperature, and the gradien	vehicle and its batteries, the weight	

of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.





Distributed by:



1300 133 120 www.aidacare.com.au

AFY118_180618