

Q6 edge[®] HD



MAXIMUM STABILITY

The Q6 Edge[®] HD is a heavy-duty mid-wheel drive power chair built for stability, strength, and control. With a high weight capacity and powerful 4-pole motors, it delivers confident handling in both indoor and outdoor settings. Its wide base offers excellent manoeuvrability, while compatibility with advanced seating, electronics, and optional iLevel[®] seat elevation up to 25.4 cm (10") makes it a versatile choice for complex needs.



| quantumrehab.com.au

Specifications

DRIVE WHEELS (Mid-Wheel Drive)	14" Pneumatic or 14x4"wide Flat Free Tyres
CASTER WHEELS	Front: 15.2 cm (6") Solid Rear: 15.2 cm (6") Solid
ANTI -TRIP WHEELS	N/A
MAXIMUM SPEED ¹	Up to 7.6 km/h (5.6 km/h w/ iLevel)
GROUND CLEARANCE	7.6 cm (3") (battery tray)
TURNING RADIUS	50 cm (19.7")
OVERALL LENGTH ²	90.4 cm (35.6")
BASE WIDTH	66.5 cm (26.2")
SEATING SIZES	TB2 HD - TB3 HD
TRU-Balance [®]	N/A
TRU-Balance [®] 3 Static Seat	TB2 - W: 40.5cm - 71cm (16-28") - D: 5.8cm - 71cm (20-28") TB3 - W: 50.8cm - 71cm (20-28") - D: 40.6cm - 71cm (16-28")
TRU-Balance [®] 3 Power Tilt	TB2 - W: 40.5cm - 71cm (16-28") - D: 50.8cm - 71cm (20-28") TB3 - W: 40.5cm - 71cm (16-28") - D: 50.8cm - 63.5cm (20-25")
TRU-Balance [®] 3 Power Lift & Tilt	TB2 - N/A TB3 - W: 50.8cm - 61cm (20-24") - D: 40.6cm - 63.5cm (16-22")
SEAT-TO-FLOOR HEIGHTS (TRU-BALANCE [®] 3)	
TRU-Balance [®] 4	TB2 - N/A TB3 - N/A
TRU-Balance [®] 3 Static Seat	TB2 - 42cm - 47cm (16.5-18.5") TB3 - 44.5cm - 49.5cm (17.5-19.5")
TRU-Balance [®] 3 Power Tilt ^{5/7}	TB2 - 46cm - 53cm (18.2-20.9") TB3 - 46cm - 52cm (18.5-20.5")
TRU-Balance [®] 3 Tilt & iLevel [®] ^{5/7}	TB2 N/A TB3 44.5cm - 50cm (17.7-19.7")
BATTERY SIZE	MK 70Ah(2)
BATTERY WEIGHT ^{5/4}	23.3 kg each
AVAILABLE ELECTRONICS	120A Q-Logic 3
RANGE PER CHARGE ³	Up to 25km with 14" solid and pneumatic tyres
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	4-poles motors
WEIGHT CAPACITY ⁴	204.1 kg
BASE WEIGHT	216 kg
CODE	PWC652227

- Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- Without front riggings
- Battery capacity may change when adding power positioning systems.
- Weight capacity is dictated by the selection of seating.
- All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 136 kg. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- Battery weight may vary +/- 1.5kg. based on manufacturer.
- All Quantum Rehab [®] units can be ordered with specialty controls.
- Call for availability
- Seat-to-floor heights are measured with the unit unloaded.

Standard Features

- Mid-Wheel 6@ allows six wheels on the ground for maximum stability
- Low impact OMNI-Casters (nylon,spherical-shaped casters)
- ATX Suspension (active-Trac[®] 3 with extra stability) incorporates front OMNI-Caster beam for enhanced performance over more varied terrain
- Q-Logic 3 Drive Controls
- ISO 7176-19 transit crash tested approved with all configurations

Options

- iLevel[®] power adjustable seat height ⁷
- TRU-Balance[®] 3 Power Positioning: 10" of power seat elevation while at walking speed (5.6 km/h) Power Tilt, Recline, and Elevating Legrests
- 10cm (14") wide flat-free or pneumatic knobby extended traction tyres^{9,4}
- TRU-Balance 2 bariatric power tilt
- USB mobile device charger
- Swing-away joystick mount
- LED fender lights
- LED lighting Package

Colours



Deep Purple

Midnight Blue

Electric Blue

Back in Black

Lime Light

Red Corvette



QUANTUM [®]
#1 FOR REHAB POWER

 **PROUDLY BUILT IN AUSTRALIA**

Distributed by
aidacare [™]
HEALTHCARE EQUIPMENT

1300 133 120 | aidacare.com.au