

Q6 EDGE® 2.0 shown in *Lime Light* with 30.5 cm (12") iLevel® technology and Q-Logic 3 Advanced Drive Controls with Bluetooth®

With features like standard 9.6 km/h (6 mph) motors and advanced responsiveness, the Q6 Edge® 2.0 is engineered to meet the performance needs of the most active user. The Q6 Edge 2.0 accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel® technology.



Specifications

I	
DRIVE WHEELS	Standard: 35.56 cm (14") Solid Optional: 35.56 cm (14") Pneumatic
CASTER WHEELS	Front: 15.24 cm (6") Solid Rear: 15.24 cm (6") Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Up to 9.6 km/h (6 mph)
GROUND CLEARANCE	6.9 cm (2.7") (battery tray)
TURNING RADIUS ²	52.1 cm (20.5") with 15.24 cm (6") casters
OVERALL LENGTH ²	32.23" with 15.24 cm (6") casters 40.48" with 20 cm (8") casters
BASE WIDTH	62.2 cm (24.4")
SEATING SIZES	
Synergy®/TRU-Balance® 3 Static Seat	Width: 30 - 56 cm (12-22") Depth: 30 - 56 cm (12-22")
TRU-Balance [®] 3 Power Tilt	Width: 30.5 - 56 cm (12-22") Depth: 30.5 - 56 cm (12-22")
iLevel®	Width: 30.5 - 56 cm (12-22") Depth: 30.5 - 56 cm (12-22")
SEAT-TO-FLOOR HEIGHTS	
Synergy®/TRU-Balance® 3 Static Seat ⁵	42 - 47 cm (16.5-18.5")
TRU-Balance [®] 3 Power Tilt ⁵	44 - 49.5 cm (17.5-19.5")
TRU-Balance® 3 Tilt & iLevel ^{® 5}	44 - 49.5 cm (17.5-19.5")
MANUAL RECLINE 6	Yes
BATTERY SIZE ³	NF-22, 55 AH
BATTERY WEIGHT 7	17.24 kg (38 lbs.)
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole Mid-Wheel 6®
WEIGHT CAPACITY ⁴	136 kg (300 lbs.)
BASE WEIGHT	63 kg (138.6 lbs.) with 15.24 cm (6") casters

1) Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure

- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/-.3175 cm (1/8") with 46 x 46 cm (18" x 18") seat, fixed towers, 136 kg (300 lbs.) weight capacity and solid tyres (add .95 cm (3/8") with pneumatic tyres). Widths may increase
- seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified. 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 1.4 kg (3 lbs.) based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability.



Features

- Standard 9.6 km/h (6 mph) motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Bold aesthetics .
- iLevel® option available with up to 30.5 cm (12") of seat elevation at walking speed (up to 5.6 km/h (3.5 mph))
- iLevel is retrofittable on the Q6 Edge[®] 2.0

Options

- iLevel[®] option allows 30.5 cm (12") of seat elevation while at walking speed (5.6 km/h (3.5 mph))
- Q-Logic 3 EX Drive Controls TRU-Balance[®] 3 Seating

TRU-Comfort 2 Seating

 Trilogy vent tray⁹ Compact vent tray⁹

LED fender lights

Accu-Trac Advanced

Tracking Technology

Colours

Grey

Street

Back

in Black



