



**Narrowest** complex rehab base on the market









quantumrehab.com

## Specifications

DRIVE WHEELS	35.6 cm or 31.8 cm pneumatic or flat free tires
CASTER WHEELS	Front: 12.7 cm flat free Rear: 15.2 cm flat free
MAXIMUM SPEED <sup>1</sup>	Up to 10 km/h Up to 5.6 km/h w/iLevel activated
GROUND CLEARANCE	7.4 cm
TURNING RADIUS <sup>2</sup>	53.3 cm
OVERALL LENGTH <sup>2</sup>	82.5 cm without front riggings
BASE WIDTH	55.25 cm with 35.6 cm drive wheels 52.71 cm with 31.8 cm drive wheels
SEATING SIZES	
TRU-Balance® 3 Static Seat	Width: 30.5 cm-51 cm Depth: 30.5 cm-51 cm
TRU-Balance® 3 Power Tilt	Width: 30.5 cm-51 cm Depth: 30.5 cm-51 cm
TRU-Balance® 3 Tilt and iLevel®	Width: 30.5 cm-51 cm Depth: 30.5 cm-51 cm
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Static Seat <sup>5</sup>	40.6 cm - 45.7 cm
TRU-Balance® 3 Power Tilt <sup>5</sup>	44.5 cm - 49.5 cm
TRU-Balance® 3 Tilt & iLevel® 5	44.5 cm- 49.5 cm
MANUAL RECLINE 6	Yes (TRU-Balance® Flex only)
BATTERY SIZE <sup>3</sup>	40 AH (2) or 55 AH (2)
BATTERY WEIGHT <sup>7</sup>	40 AH: 14 kg each 55 AH: 18.1 kg each
AVAILABLE ELECTRONICS 8	120A Q-Logic 3 EX 75A Q-Logic NE 75A Q-Logic NE+
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel 6®
WEIGHT CAPACITY <sup>4</sup>	136 kg
BASE WEIGHT WITHOUT SEATING SYSTEM	57.3 kg

## **Features**

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 30.5 cm of seat elevation at walking speed (up to 5.6 km/h)
- Fender lights and USB charger included
- Made in America of foreign and domestic parts

## Options

- iLevel® option allows 30.5 cm of seat elevation while at
- Trilogy vent tray<sup>9</sup>
- walking speed (5.6 km/h)
- Compact vent tray<sup>9</sup>
- Q-Logic 3 EX Drive Controls
- 30.5 cm x 30.5 cm seat sizes with iLevel (SOF only)
- TRU-Balance® 3 Seating

## Colours\*



- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 0.3 cm with 45.7 cm x 45.7 cm seat, fixed towers, 136 kg weight capacity and solid tires (add 0.95 cm with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 1.4 kg based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability.





