

SEDEO PRO SEATING SYSTEM

GOOD SEATING STARTS WITH THE PERSONAL NEEDS OF THE CLIENT.

ERGONOMIC SEATING IS THE BASIS OF AN OPTIMAL MOBILITY SOLUTION.

It is all about creating a seating posture which makes it possible to achieve your goals and aspirations.

All of this is fully reflected in our new SEDEO PRO seating system. For moderate to complex clinical needs, the seating system combines a highly adjustable seat frame with an excellent choice of postural support options, including a wide range of:

- Seat & Backrests cushions
- Armrests, Side Guards & Thoracic Lateral Supports
- Legrests & Knee support pads
- Headrests

For the first time both JAYTM and WHITMYERTM are incorporated to offer advanced positioning capabilities. Truly modular in design, SEDEO PRO provides a new level of functionality, comfort and design. It can be easily customised to suit individual requirements and quickly adapted to changing needs of the user.














SEDEO PRO QUICK AND EASY TO ADJUST



CREATING THE RIGHT SEATING POSTURE FOR YOU — EASY TO DO FOR THE HEALTHCARE PROFESSIONAL

- Total flexibility / fully adjustable
- Intuitive adjustments. One tool adjustment (5 mm Allen key) / one simple working routine / self-explaining measuring scales / no extra parts are required
- Parts can be steplessly adjusted, repositioned, widened, turned to create the best support
- Adjustments can be done with the user in the chair, easily reachable from the side.



2 = Additional armrest width: + 2 cm to - 2 cm each side

3 = Seat depth: 40 - 56 cm

4 = Back height: 44 - 62 cm

SEDEO PROSEATING FRAME

THE WIDE RANGE OF SEAT WIDTH AND SEAT DEPTH ADJUSTMENT MAKES THAT SEDEO PRO IS IDEAL. ALWAYS THE RIGHT FIT, YOU CAN'T ORDER WRONG.









HIGHLY CONFIGURABLE

SEATING FRAME

FOUR BACK OPTIONS



Tool Angle Adjustable Back, 90°-120° Stepless adjustment, easy to adjust with the user in the chair.



User Adjustable
Manual Recline Back, 90°-130°
Shifting body weight made easy with
gas spring. Operating lever for user
adjustment under the armrest.



Powered Recline, 90°-130° Comfortable and easy reclining with the press of a button.



Powered Recline with 100 mm Anti-Shear, 90°-150° lts patented unique mechanism enhances skin protection and allows the positioning elements and even chin control to stay in place during the full recline operation.

FOUR ARMREST TYPES WITH SEVERAL ARM PAD AND SIDE GUARD OPTIONS



Armrest, Single Post, Flip Back



Armrest, Reclining, Flip Back




Armrest, Single Post, Removable



CONFIGURE THE SEDEO PRO TO PERFECTLY
FIT THE USER BY SELECTING COMPONENTS
FROM THE WIDE VARIETY OF OPTIONS. WITH
NO COMPROMISES – ALWAYS A PERFECT FIT.

TWO CENTER MOUNT LEGRESTS WITH CALF PLATE OPTIONS







Tool Adjustable Centre Mount Legrest

THREE SWING AWAY LEGRESTS WITH CALF PLATE OPTIONS



SEDEO PRO Swing Away Legrest, 70°



SEDEO PRO Manual Elevating, Swing Away Legrest



SEDEO PRO Power Elevating, Swing Away Legrest

SEATING OPTIONS:

JAY CUSHIONS & BACKS











your needs.

A wide range of JAY cushions are available and compatible

with the SEDEO PRO frame to meet higher clinical needs

Designed for pelvic stability, skin protection and comfort, JAY cushions and backs are low maintenance, highly customisable and are available in a large variety to meet

and ensure we offer the best seating possible.



Fusion Fluid & Air

J2







Superior Clinical Seating.

SEDEO Back Cushion

J3 PDC

J3 PDL

J3 PD

JAY J3 BACKS

J3 PL

SEDEOHEADRESTS



SEDEO HEADREST PAD, SMALL

240 x 150 mm, Cover optional



SEDEO HEADREST PAD, FLEXIBLE 400 x 110 mm



=WHITMYER[®]

A WIDE RANGE OF WHITMYER HEADRESTS ARE
AVAILABLE WITH THE SEDEO PRO SEATING SYSTEMS
TO MEET EVERY USER'S CLINICAL NEEDS AND ENSURE
WE OFFER THE BEST SEATING POSSIBLE.



LINX





AXYS



AXYS Offset




Cobra XTRA Flip Back





SEDEO PRO

POSITIONING SUPPORT ITEMS




LATERAL THORACIC SUPPORTS

- Recommended with Profiled back or Profiled adjustable back
- Depth, width, height and angle adjustable
- Integral slim-line swing away mechanism

Pad sizes:

Short pad: 90 x 120 mm Long pad: 150 x 130 mm




POMMEL

- Depth and height adjustable
- Swing away mechanism


Pad (wxdxh): 70 x 190 x 100

mm



POSITIONING BELT

- Carseat Bucke
- Short: 50 x 1700 mm
- Long: 50 x 2700 mm



KNEE & THIGH SUPPORTS

- Depth, width and height adjustable
- Swing away mechanism

Pad sizes:

Small pad: 120 x 120 mm Large pad: 180 x 120 mm



TRANSFER BAR

THE POWER OF QUICKIE

A complete power range from basic mobility to complex rehab. Whatever you need - without compromise.





Q300

SEATING SEDEO PRO SEDEO LITE

BASES MWD RWD



Q400

SEATING SEDEO LITE **BASES** MWD FWD RWD



Q500

SEATING SEDEO PRO SEDEO LITE BASES MWD FWD RWD



Q700

SEDEO ERGO SEDEO PRO

BASES MWD **FWD** RWD



Q700-UP

SEATING SEDEO ERGO UP MWD

FWD

BASE PERFORMANCE

